

Product datasheet for SC317719

PNMA8A (NM_018215) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PNMA8A (NM_018215) Human Untagged Clone
Tag:	Tag Free
Symbol:	PNMA8A
Synonyms:	PNMAL1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC317719 representing NM_018215. Blue=Insert sequence Red=Cloning site Green=Tag(s)

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
 GATCCGGTACCGAGGAGATCTGCCGCC**CGATCGCC**
 ATGTCCAAGACCATGGCGATGAACCTTCTGGAGGATTGGTGCAGGGGAATGGAAGTGACATCCACAGG
 TCCTTGTGGTACAGGCATCCAGAGGACTGTGGGCAGGCAGAAATTGAGGAGACCTTGAATGGGGTC
 CTCTCCCCACTGGGCCGTACCGCGTGCTCAACAAGATTTTGTGAGGGAAGAGAATGTTAAAGCTGCC
 CTCATTGAGGTTGGTGAAGGTGTGAATCTGAGCACCATCCCCGTGAATCCCAGGAAGGGTGGTGTG
 TGGAGAGTGGTCTGTAGAGACCTACCCAGGATGCCGAGTTTTTAAAAATCTGAATGAATTCCTGGAT
 GCCGAGGGGCGCACCTGGGAGGATGTGGTCCGCCTGCTCCAGCTCAACCACCCACCTGTCCAGAAC
 CAGCATCAGCCCCAGAGAACTGGGCAGAAGCTCTGGGGTGCTTCTGGGAGCAGTGGTGCAGATCATC
 TTCTGCATGGATGCCGAGATCCGCAGCCGGGAGGAAGCCAGGGCCAGGAGGCCGCTGAATTCGAGGAG
 ATGGCAGCCTGGGCTTTAGCAGCAGGGAGGAAGGTGAAGAAAGAACCGGGGCTTGACGACGAGGTGGGT
 TCTGCCTTAAAGGCAGAGACCCCAACAACCTGGAATGCCACGGAAGACCAGCATGAGCCTACCAACCT
 TTGGTTCGAGGGCTGGAGCTAAGTCTCGCTCCAGGAGAAAGAAGCAGAAAGAAGAACTCCAGGCAGGAA
 GCAGTGCCTGGAAGAAACCCAAAGGCATCAATTCCAACAGCACAGCTAATTGGAGGATCCTGAGGTG
 GGTGATGCTGAAAGCATGGCGATCTCAGAGCCGATCAAGGGCAGCAGAAAGCCCTGTGTGAATAAGGAG
 GAGTTGGCTTTGAAGAAGCCCATGGCGAAATGTGCCTGGAAGGGTCCCAGAGAGCCACCTCAGGATGCC
 CGGGCAGAAGCCGAGAGCCAGGAGGCGCCTCTGAGTCAGACCAAGATGGTGGCCATGAAAGCCCA
 AAGAAGAAGGCCGTGGCCTGGGTGTCTGCCAAGAACCCGCTCCCATGAGGAAGAAGAAGGTGAGC
 TTGGGCCCTGTCTCCTACGTCTTGGTTGACTCAGAAGATGGCAGGAAGAAGCCGGTGATGCCAAAGAAA
 GGGCCAGGCTCAAGAAGGGAGGCATCAGATCAGAAGGCCCTCGGGGCCAGCAGCCTGCCGAGGCAACA
 GCCTCAACCTCTAGGGTCCGAAGGCCAAGCCAGAAGGCTCTCTCGGCGTGCCACCAATGAATCCAGA
 AAGTT**TGA**
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
 TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC


[View online »](#)

Restriction Sites:	Sgfl-Mlul
ACCN:	NM_018215
Insert Size:	1320 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_018215.3</u>
RefSeq Size:	3708 bp
RefSeq ORF:	1320 bp
Locus ID:	55228
UniProt ID:	<u>Q86V59</u>
Cytogenetics:	19q13.32
MW:	48.2 kDa