

Product datasheet for SC317728

CHST9 (NM_031422) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CHST9 (NM_031422) Human Untagged Clone
Tag:	Tag Free
Symbol:	CHST9
Synonyms:	GALNAC4ST-2; GalNAc4ST2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC317728 representing NM_031422. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTGTAGTGAACCGTCAGAATTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGCAGCCATCTGAAATGGTTCATGAACCCAAACAAGTCTTCTCTCTGTGCTGATATTTGGAGTAGCT
GGGCTACTCCTCTTCATGATTTTGAAGTCTGGATTGAAGAACAACATACAGGGAGAGTGGAGAAGAGA
AGAGAACAAAAGTAACCTCAGGATGGGGACCAGTGAAGTACTTGCAGCCTGTACCCAGAATCATGAGT
ACAGAAAAATCCAGGAACATATCACCAACCAGAACCCCAAGTTTCACATGCCTGAGGATGTACGAGAA
AAAAAGGAAAATCTTCTACTCAATCTGAGAGATCTACTAGGCTCTTAACAAAGACCAGTCATTCACAA
GGAGGGGATCAAGCTTTAAGTAAGTCCACAGGGTCACCAACAGAGAAGTTGATTGAAAAAGTCAAGGA
GCTAAGACTGTTTTAAACAAGTTCAGCAACATGAATTGGCCAGTGGACATTCACCCTTTAAACAAAAGT
TTAGTCAAAGATAATAAATGGAAGAAAAGTCTGAGGAGACCCAAAGAGAAAACGAAGGCTTTCCCTTCAGGAG
TTTTGCAAGAAATACGGTGGGGTGAATCATCATCAGTCACATCTTTTTATACAGTATCCAGAATCTAT
GTAGAAGATAACACAAAATCTTATATTGTGAGGTACCTAAGGCTGGCTGTTCCAATTGGAAAAGAATT
CTGATGGTACTAAATGGATTGGCTTCTCTGCATACAACATCTCCCAATGCTGTCCACTACGGGAAG
CATTTGAAGAAGCTAGATAGCTTTGACCTAAAAGGGATATATACCCGCTTAATACTTACACCAAAGCT
GTGTTTGTTCGTGATCCCATGGAAGATTAGTATCAGCCTTTAGGGACAAATTTGAACACCCCAATAGT
TATTACCATCCAGTATTCGAAAAGGCAATTATCAAGAAATATCGACCAAAATGCCTGTGAAGAAGCATT
ATTAATGGATCTGGAGTCAAGTTCAAAGAGTTTATCCACTACTTGTGGATTCCCACCGTCCAGTAGGA
ATGGACATTCAGTGGAAAAGGTCAGCAAACTCTGCTATCCGTGTTTGTCAACTATGATTTTGTAGGG
AAATTTGAGACTTTGGAAGAAGATGCCAATTACTTTTTACAGATGATCGGTGCTCCAAAGGAGCTGAAA
TTTCCCAACTTTAAGGATAGGCACTCTCCGATGAAAGAACCAATGCTCAAGTCGTGAGACAGTATTTA
AAGGATCTGACTAGAAGTGAAGACAATTAATCTATGACTTTTACTTGGACTATTTAATGTTTAAT
TATACAACCTCCATTTTGTAG
ACGCGTACGCGGCCGCTCGAGCAGAAAACATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```



[View online »](#)

Restriction Sites: Sgfl-MluI

Plasmid Map:



ACCN: NM_031422

Insert Size: 1332 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_031422.5](#)

RefSeq Size: 2280 bp

RefSeq ORF: 1332 bp

Locus ID: 83539

UniProt ID: [Q7L1S5](#)

Cytogenetics: 18q11.2

Protein Families: Transmembrane

MW: 52.1 kDa

Gene Summary: The protein encoded by this gene belongs to the sulfotransferase 2 family. It is localized to the golgi membrane, and catalyzes the transfer of sulfate to position 4 of non-reducing N-acetylgalactosamine (GalNAc) residues in both N-glycans and O-glycans. Sulfate groups on carbohydrates confer highly specific functions to glycoproteins, glycolipids, and proteoglycans, and are critical for cell-cell interaction, signal transduction, and embryonic development. Alternatively spliced transcript variants have been described for this gene. [provided by RefSeq, Aug 2011]

Transcript Variant: This variant (1) represents the longer transcript and encodes the longer isoform (1).