

## Product datasheet for SC317748

### KRT35 (NM\_002280) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	KRT35 (NM_002280) Human Untagged Clone
Tag:	Tag Free
Symbol:	KRT35
Synonyms:	Ha-5; HA5; hHa5; K35; KRTHA5
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC317748 representing NM_002280. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTGTAGTGAACCGTCAGAATTTTGTAAACGACTACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCTTCCAATGCCTCAAGGCCGGTCTCTTCTGGGTCTCTCAAGAGCCAGGAGGGGCCAGTGGG
GGCTCCACTCGTGTGTCCGCAATGTAAGTCCAGCAGCTCTTGAAGCTTCCAAGTCTCTCCCTGTGGCC
AGAAGTTTCTCTGCCTGCTCAGTGGTCTGGGCAGAAGCAGCTACAGGGCCACCAGCTGCCTCCCTGCT
CTCTGCCCTCCCTGCTGGAGGCTTCGCTACCAGCTACAGTGGGGTGGGGGCTGGTTGGGGAGGGCAGC
CTCACTGGCAATGAGAAGGAGACCATGCAATCCCTGAACGACCGCTGGCCGGCTACCTGGAGAAGGTG
CGTCAGCTGGAGCAGGAGAACGCCAGCCTGGAGAGCCGATCCGTGAGTGGTGTGAGCAGCAGGTCCCC
TACATGTGCCCTGACTACCAGTCTACTTCCGGACCATCGAGGAGCTCCAGAAGAAGACTCTATGCAGC
AAGGCTGAGAATGCCAGGCTGGTGGTGGAGATTGACAATGCCAAATTGGCTGCAGATGACTTCAGGACC
AAGTATGAGACGGAGGTGTCCCTGCGGCAGCTGGTGGAGTCCAGACATCAACGGCCTGCGCAGGATCCTG
GATGACCTGACCTGTGCAAGTCTGACCTGGAGGCCAGGTGGAGTCCCTGAAGGAGGAGCTGCTCTGC
CTGAAGAAGAACCATGAGGAGGAAGTGAAGTCACTGCGCTGCCAACTTGGTGACCGCCTCAATGTTGAG
GTGGAGCTGCCCCACCTGTTGACCTGAACCGAGTCTGGAGGAGATGAGGTGCCAGTATGAAACCTG
GTGGAGAATAACCGCCGGATGCTGAAGACTGGTTGGACACCCAGAGTGAAGGAGTGAACCCAGGAGTG
GTGCCAGCTCAGAGCAGTTGCAGTCTGCCAGGCAGAGATCATCGAGCTGAGACGCACGGTCAACGCC
CTGGAGATTGAGCTGCAGGCTCAGCACAGCATGAGAGATGCTTTGGAATCCACCCTGGCAGAGACGGAG
GCCCGCTATAGCTCCCAGCTGGCCAGATGCAGTGCATGATCACCAACGTGGAGGCCAGCTGGCCGAG
ATCCGGGCTGACCTGGAGCGGCAGAACCAGGAGTACCAGGTGCTGCTGGACGTCGGGCCCGGCTGGAG
TGTGAGATCAACACGTACCGGGGCTGCTGGAGAGTGAAGACAGCAAGCTCCCCTGTAACCCATGTGCA
CCTGACTACTCACCTCCAAGTCATGCCTTCCCTGTCTTCTGCGGCTCCTGCGGTCTAGTGCAGCC
CGCACAACTGCAGCCCCGCCCATTTGTGTGCCCTGCCAGGGGGTTCGGTTCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTAAACGGCCGCGC
```



[View online »](#)

<b>Restriction Sites:</b>	Sgfl-Mlul
<b>ACCN:</b>	NM_002280
<b>Insert Size:</b>	1368 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<u><a href="#">NM_002280.4</a></u>
<b>RefSeq Size:</b>	1688 bp
<b>RefSeq ORF:</b>	1368 bp
<b>Locus ID:</b>	3886
<b>UniProt ID:</b>	<u><a href="#">Q92764</a></u>
<b>Cytogenetics:</b>	17q21.2
<b>MW:</b>	50.4 kDa
<b>Gene Summary:</b>	The protein encoded by this gene is a member of the keratin gene family. This type I hair keratin is an acidic protein which heterodimerizes with type II keratins to form hair and nails. The type I hair keratins are clustered in a region of chromosome 17q12-q21 and have the same direction of transcription. [provided by RefSeq, Jul 2008]