

Product datasheet for SC317749

C14ORF140 (ZC2HC1C) (NM_024643) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	C14ORF140 (ZC2HC1C) (NM_024643) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZC2HC1C
Synonyms:	C14orf140; FAM164C
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC317749 representing NM_024643. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTGTAGTGAACCGTCAGAATTTTGTAAACGACTACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCTGGTCTCCAGCGTTGGCGTCACATCTGCCTGTGGCGTTATGCTCCACATAATACAACGGAA
GCTCCAGGGCCCCACTCAGCCAAGCAAGACTCTTACGAACAAGGTGACTCTCCAGCAGTCCTGAAG
GGGCACCTGAGGAACAATTTCCAGAAGCAGCTTTTGGCAACAAGAGTTGATTCTGGATAAAGTCTAT
ACTCACCCAAAATGGAACACCCAAAACAAAAGCCGGAGCTACTCTATCCCACTGTACTGGAATCAGC
CAGCAAGATCCAGAAAGTATTCCAGGGCCAAGGAAATGGTTTGTACTCGTCAGGCCCTCAATCC
TGGTATCCAAAGCCAATAACCAGGACTTTATCCCTTTACAAAGAAACGAGTTGGAGTGGACCGGGCG
TTCCATTGAAACCCATGGTCCACAGGAAGTCGTGCAAGTACAGGTGAGGCTGGCACTGATGGGGACCAT
AATGTCTACCCAAGGCCCTGAGCCGAGAGATTTTCATCTAGGAACTTTGGTGTGAGGAACAGGGC
AACTTTTCTGTGGTTGGTACTGTTCTTGTGCCACGCAGGCGGAGAAGGCCGTGGCAAATTTGACAGG
ACGGAGTGGGTGCAGATCCGAAGACTAGAAGCTGCAGGGGAGAGCTTAGAGGAGGAAATCCGAAGAAAG
CAGATTCTCTGAGGGGAAAGCTGAAGAAGACAGAGGAGGAACTCAGAAGGATCCAGACGCAAAAGGAA
CAGGCCAAGGAAAATGAAAACGGAGAGCTACAGAAAATATACTCCCAAGGAGCAGAGTTAAAGGTAAT
AAAAGCAACACCATGTACAAACCTATCTTCTCCCAAGAAATTTGAGTTTGGGAAGAATTTAGTAGAGAC
AGGAGAGAGGATGAAACTTGGGACGGTCTCAACAAAATTCAGGTCCATTCCAGTCTCTGATTATAGA
ATCCAGAGGCTCAAAAGGGAAAGGCTGGTAGCAAGCAATAATAAAATTCGAGACCCAGTCTCAGAGCCA
TCGGTGGAGAAATTTCTCCCGCCTTCAGAAACACCAAGTCCGGTGTCTTGCAGGGATCCGCCAGAAATTC
AGCCTGTCCATGGCACCAGACTCCTCAGGTTCCAGCGGCTCCATTGAAGAGCCACAGCTGGGTGAGTGC
AGCCACTGTGGGCGCAAATTCCTCTCGTTCAGGCTGGAGAGACTCCAACATCTGCAGCAGGATGCGG
GGTTCCAAGAGGAAAGTGTGACTCCTCCAGGGCCCGGCTAAGGGCACAGAACTAGAGCAGTACTTG
AACTGGAAGGGCCAGCTTCAGCCAAGGCTGAACCTCCTCAGAAGAGCAACTGGAGATAG
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```



View online »

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_024643

Insert Size: 1371 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_024643.3](#)

RefSeq Size: 3986 bp

RefSeq ORF: 1371 bp

Locus ID: 79696

UniProt ID: [Q53FD0](#)

Cytogenetics: 14q24.3

MW: 51.7 kDa