

Product datasheet for **SC317775**

ZNF385B (NM_152520) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF385B (NM_152520) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF385B
Synonyms:	ZNF533
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC317775 representing NM_152520.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGAATATGGCAAATTTCTACGGGGCTTTGAAGAAAAGGGGATAAAGAACGACAGGCCTGAGGACCAG
TTGAGCAAAGAGAAAAAGAAAATTTCTTCTCCTTCTGTGAGGTGTGAACATCCAGCTGAACCTCTGCA
GCCCAGGCTCAGGTGCATTCCAACGGCAAATCCCACCGCAAACGAGTGAAGCAGCTGAGTGATGGGCAG
CCCCACCAACCCGCCAGGCCAGCCCCAGCAGCAACAGCAGCACAGGCAGTACATGCCACACTACTACC
CTTCTGCTCTTGTGCGCACACCTACCCTGATGATGCAGCCTTCACTGGATATCAAACCTTTATGTCT
TTTCCAGTGGACAGTAGTTCTGCTGTTGGGCTCTTTCAAATTTTAACACAATGGATCCTGTGCAAAAA
GCGGTTATTAACCATACATTTGGAGTATCCATTCCTCCCAAGAAAGAAACAAGTTATTTCTGTAAATGTC
TGTCAGCTTCGCTTTAACTCAGATAGCCAGGCCGAGGCCCACTACAAAGGAAGTAAACATGCCAAGAAG
GTCAAAGCACTAGACGCAACGAAAAATAAACCCAAAATGGTTCCTTCAAGGACAGCGCAAAGGCTAAT
CCCAGCTGCTCCATCACTCCAATCACAGGCAACAACCTCTGACAAATCAGAAGATAAAGGGAAGTTAAAA
GCCAGCAGTTCAGTCAGCCATCAAGCTCTGAAAGTGGCTCATTTCTCCTCAAATCTGGCACAACACCC
CTGCCACCTGGAGCAGCCACTTCTCCCTCCAAGAGCACAAATGGAGCTCCCGGTACTGTTGTTGAATCA
GAAGAAGAAAAAGCCAAAAAATTACTTTATTGTTCACTATGCAAAGTGGCTGTGAATCCCTGTCACAG
CTAGAGGCACACAACACAGGATCTAAACACAAGACCATGGTTGAAGCTCGTAATGGGGCTGGTCCAATT
AAATCCTATCCTAGACCTGGATCAAGATTAAAGATGCAGAATGGCAGTAAGGGGTCAGGACTACAGAAC
AAGACATTTTATTGTGAAATCTGTGATGTTTCAATGTTAATTCAGAAATCAACTCAAACAGCACATTTCT
AGCCGAAGGCATAAAGATCGAGTTGCAGGGAACCACTGAAGCCAAAATACAGCCCTTACAACAACTC
CAGCGGAGCCCGAGTATTCTAGCAGCAAACTTGCAATTCAGAAAGATATGATGAAGCCTTTGGCCCCA
GCCTTCCTGTCTCCTCACCTCTCGCAGCGGCGGAGCCGTGCTCCTCAGCGCTGTCACTCCACCCCGGCC
TCTGCCTCGCTCTTCCAGGCTCCAGCCATTCTCCAGCTCTTCTGAGGCCTGGGCATGGGCCCATCCGC
GCCACTCCTGCCTCCATCCTCTTTGCTCCGTACTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_152520

Insert Size: 1416 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_152520.4</u>
RefSeq Size:	3320 bp
RefSeq ORF:	1416 bp
Locus ID:	151126
UniProt ID:	<u>Q569K4</u>
Cytogenetics:	2q31.2-q31.3
Domains:	ZnF_U1, zf-C2H2
MW:	50.4 kDa
Gene Summary:	May play a role in p53/TP53-mediated apoptosis.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (1) represents the longest transcript and encodes the longest isoform (1).