

Product datasheet for **SC317791**

SLC25A24 (NM_013386) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SLC25A24 (NM_013386) Human Untagged Clone
Tag:	Tag Free
Symbol:	SLC25A24
Synonyms:	APC1; SCAMC-1; SCAMC1
Vector:	<u>pCMV6 series</u>
Fully Sequenced ORF:	>NCBI ORF sequence for NM_013386, the custom clone sequence may differ by one or more nucleotides <pre> ATGTTGCGCTGGCTGCGGGACTTCGTGCTGCCACCGCGGCTGCCAGGACGCGGAGCAG CCGACGCGCTACGAGACCCTCTCCAGGCACTGGACCGCAATGGGGACGGAGTGGTGGAC ATCGGCGAGCTGCAGGAGGGGCTCAGGAACCTGGGCATCCCTCTGGGCCAGGACGCCGAG GAGAAAATTTTACTACTGGAGATGTCAACAAAGATGGGAAGCTGGATTTTGAAGAATTT ATGAAGTACCTTAAAGACCATGAGAAGAAAATGAAATTGGCATTAAAGAGTTTAGACAAA AATAATGATGGAAAAATTGAGGCTTCAGAAATTGTCCAGTCTCTCCAGACACTGGGTCTG ACTATTTCTGAACAACAAGCAGAGTTGATTCTTCAAAGCATTGATGTTGATGGGACAATG ACAGTGGACTGGAATGAATGGAGAGACTACTTCTATTTAATCCTGTTACAGACATTGAG GAAATTATCCGTTTCTGGAACATTCTACAGGAATTGACATAGGGGATAGCTTAACATT CCAGATGAATTCACGGAAGACGAAAAAAATCCGGACAATGGTGGAGGCAGCTTTTGGCA GGAGGCATTGCTGGTGTCTCTCGAACAAGCACTGCCCTTTGGACCGTCTGAAAATC ATGATGCAGGTTACGGTTCAAATCAGACAAAATGAACATATTTGGTGGCTTTTCGACAG ATGGTAAAAGAAGGAGGTATCCGCTCGCTTTGGAGGGGAAATGGTACAAACGTCATCAA ATTGCTCCTGAGACAGCTGTTAAATTCTGGGCATATGAACAGTACAAGAAGTTACTTACT GAAGAAGGACAAAAATAGGAACATTTGAGAGATTTATTTCTGGTTCCATGGCTGGAGCA ACTGCACAGACTTTTATATATCCAATGGAGGTTATGAAAACCAAGGCTGGCTGTAGGCAAA ACTGGGCAGTACTCTGGAATATATGATTGTGCCAAGAAGATTTTGAACATGAAGGCTTG GGAGCTTTTACAAAGGCTATGTTCCCAATTTATTAGGTATCATACCTTATGCAGGCATA GATCTTGCTGTGTATGAGCTCTTGAAGTCCTATTGGCTGGATAATTTGCAAAAGATTCT GTAAACCCTGGAGTCATGGTGTGCTGGGATGCGGTGCCTTATCCAGCACCTGTGGTCAG CTGGCCAGCTACCCATTGGCTTTGGTGAGAACTCGCATGCAGGCTCAAGCCATGTTAGAA GGTTCCCCACAGCTGAATATGGTTGGCCTCTTTCGACGAATTATTTCAAAGAAGGAATA CCAGGACTTTACAGAGGCATCACCCAACTTCATGAAGGTGCTCCCTGCTGTAGGCATC AGTTATGTGGTTTATGAAAATATGAAGCAAACCTTAGGAGTAACCCAGAAA </pre>
Restriction Sites:	Please inquire
ACCN:	NM_013386
Insert Size:	3600 bp


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_013386.3</u> , <u>NP_037518.3</u>
RefSeq Size:	3600 bp
RefSeq ORF:	3600 bp
Locus ID:	29957
UniProt ID:	<u>Q6NUK1</u>
Cytogenetics:	1p13.3
Protein Families:	Druggable Genome
Gene Summary:	This gene encodes a carrier protein that transports ATP-Mg exchanging it for phosphate. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, May 2012] Transcript Variant: This variant (1) represents the longer transcript and encodes the longer isoform (1).