

Product datasheet for **SC317795**

ZNF187 (ZSCAN26) (NM_001023560) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF187 (ZSCAN26) (NM_001023560) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZSCAN26
Synonyms:	SRE-ZBP; SREZBP; ZNF187
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC317795 representing NM_001023560.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCAACAGCATTGGTGAGTGCCCATCCCTGGCTCCCCTGAATCTGAAGAAGGAGGGCTTCGGGTA
GTGAGGGAGGATCACTACTCTACTTTGGGAACAGGGATTCAAGCTGCAAGGAAACAGTAAAGGCCTTGGA
CAGGAGCCATTGTGCAAACAATTCAGGCAGTTGCGTTATGAAGAGACCACAGGACCTCGAGAAGCACTA
AGTCGGCTCCGGGAGCTCTGTCAACAGTGGCTACAGCCGAGACCCATACCAAGGAGCAGATCCTGGAG
CTGCTGGTCTGGAGCAGTTTCTGATCATCCTGCCTAAGGAGCTCCAGGCCCGGGTGACAGGAGCATCAC
CCAGAGAGCAGGGAGGACGTGGTTGTTGTTCTGGAGGATTTGCAGCTGGATCTTGAGAAACAGGACAA
CAGGTGGACCCAGACCAGCCAAAGAAACAAAAATACTTGTGGAGGAGATGGCCCTCTGAAAGGAGTA
CAGGAACAGCAGGTTCCGCATGAGTGTGAAGTTACAAAGCCTGAGAAAGAGAAGGTGAGGAGACAAGG
ATTGAGAAATGGGAAGCTTATTGTAGTAACAGACTCTTGTGGAAGAGTAGAGTCATCTGGGAAAATATCT
GAACCCATGGAGGCTCATAATGAGGGCTCTAACTTGGAAAGGCATCAGGCCAAGCCCAAAGAGAAGATT
GAGTATAAATGCTCAGAACGTGAGCAGAGATTCATCCAGCACTTGACCTGATTGAACATGCGAGTACA
CACACGGGAAAGAACTCTGCGAGTCTGATGTGTGTCAGAGTTCAGTCTTACAGGACATAGAAAGTC
CTCTCTAGAGAGAAAGGTCATCAGTGTGATGAGTGTGGGAAAGCCTTTCAGAGGAGTTCACACCTCGTC
AGACATCAGAAAAATCCATCTTGGTGAGAAGCCTTATCAGTGCAATGAGTGTGGCAAAGTCTTTAGCCAG
AATGCAGGCCTTTTGGAACATCTCAGAATTCATACTGGAGAGAAACCTTATCTATGTATCCATTGTGGA
AAAAATTTAGGGCAGCTCTCACCTTAATCGACATCAGAGAATTCACAGTCAGGAGGAGCCCTGTGAG
TGCAAGGAGTGTGAAAAACCTTTAGTCAGGCCCTACTCCTCACCCACCATCAGAGAATCCATAGTCAC
TCCAAAAGCCATCAATGTAACGAGTGTGAAAAGCTTTTCAGTTTGACCTCAGACCTTATTCGACACCAC
AGAATTCATACTGGAGAAAAACCTTTCAAGTGTAAACATATGCCAGAAAGCCTTCGGACTAAACTCACAC
CTTGCTCAGCATGTAAGAAATCCACAATGAAGAAAAACCTATCAGTGTAGTGAATGTGGAGAAGCCTTC
AGGCAAAGGTCAGGTCTTTTTCAACATCAGAGATATCACCAAAAGACAAACTGGCTTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites:

SgfI-MluI

Plasmid Map:



ACCN: NM_001023560

Insert Size: 1440 bp

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001023560.3</u>
RefSeq Size:	2850 bp
RefSeq ORF:	1440 bp
Locus ID:	7741
UniProt ID:	<u>Q16670</u>
Cytogenetics:	6p22.1
Protein Families:	Transcription Factors
MW:	55.4 kDa
Gene Summary:	<p>May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (a) represents the longest transcript and encodes the longest isoform (a).</p>