

Product datasheet for **SC317803**

SYNC (NM_030786) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SYNC (NM_030786) Human Untagged Clone
Tag:	Tag Free
Symbol:	SYNC
Synonyms:	SYNC1; SYNCOILIN
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)


Fully Sequenced ORF: >SC317803 representing NM_030786.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGCCAGCCCGGAGCCCCGGCGCGGGGACGGCCCGCCAGGCCGCGAGGAAAAACAAGAGTAGAG
GCCAATTCTCCTTTCCAAAGAAGCTGGATCCCTAAATGAGGCAGAAGCCTTGAACCCAGAAGTTACT
CTATCTTCAGAGGGTCCCTAAACCTCGAAGACATTCTCTACCTGGAGGACACAGGTGACCTTGATGAG
ACACTCTATGTGCAAGAGACTGAGAAGGCAGAGGAGGCCCTGTATATTGAGGAGGCCATGCAGCCAGAT
GAGGCTCTGCATGTGGAGGAGCCTGGGAATCCAGAGGAGACAGTGTGTGTGGAGGAAACCACGGAGCCA
GATCGGATACAGTTTGTGGAGGGGCCCGTGGAGCCAGGAAAGCCACAAGCCAGAGCAGCTTGTAT
GAGGGAGAGACAGTCACAAGGGCGGAGAAATCTAACCTGAGGAGAGCCTCAGAGCCGAGCAGAGCCCC
AGCGTGGAGGAGAACCTGAGCATAGAGGACCTGGAATTGCTAGAGGGCGCTTCCAGCAGTGTGTCCAA
GCTGTGGCCAGCTGGAAGAGGAGAGGGATCAGCTCATCCATGAGCTTGTATTGCTCCGGGAACAGCC
CTGAGGAGGTACAGCAAGTCCATCAAGACATCCTGGCTGCCTACAAGCTCCATGCCAAGCAGAGCTG
GAGAGAGATGGCCTAAGGAGGAGATCCGGCTGGTCAAGCAGAAGCTTTTCAAAGTGACAAGGAATGT
GTGGCCTACCAATACCAGCTGGAGTGCCGCCAGCAGGACGTGGCTCAGTTTGCCGATTTCCGGGAAGTG
CTGACTACAAGGGCAACCCAGCTCTCAGAGGAAGTGGCCAGCTCCGGGATGCCTATCAGAAGCAGAAG
GAGCAGCTGCGGCAACAAGTGAAGCCCTCAAGCCAGAGGGATGGGCACCTTCTCCAGGAAAGCCGG
CGACTCTCTGCCAGTTTGAAAATCTCATGGCAGAGAGCCGCCAGGACCTGGAGGAGGATATGAGCCT
CAGTTCCTGCGGCTCCTAGAGAGGAAAGAAGCTGGGACCAAGCTCTGCAGAGAACCAGGCTGAGATC
CAGGAAATGAAGGAGGCTCTGAGACCCCTGCAAGCAGAGGCCCGGCAGCTCCGCTGCAAAAACAGGAAC
CTGGAGGACCAGATCGCACTTGTGAGGCCAAAACGAGATGAAGAGGTGCAGCAGTACAGGGAACAGCTG
GAGGAAAATGGAAGAACGCCAGAGGAGTAAAGAAATGGGGTGCAACTCCAGCAACAGAAGAACAAAGAG
ATGGAACAGCTAAGGCTCAGTCTTGTGAAGAGCTCTCTACTTATAAGGCTATGCTACTACCCAAGAGC
CTGGAACAGGCTGATGCTCCCACTTCTCAGGCAGGTGGAATGGAGACACAGTCTCAAGGGCTGTTTAG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_030786

Insert Size:	1449 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_030786.2
RefSeq Size:	3493 bp
RefSeq ORF:	1449 bp
Locus ID:	81493
UniProt ID:	Q9H7C4
Cytogenetics:	1p35.1
MW:	55.3 kDa
Gene Summary:	<p>This gene encodes a member of the intermediate filament family which contains an N-terminal head domain, followed by a central coiled-coil region and a short C-terminal tail. The protein is highly expressed in skeletal and cardiac muscle. The protein links the dystrophin associated protein complex (DAPC) to desmin filaments in muscle and may have a structural role in striated muscle. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jun 2009]</p> <p>Transcript Variant: This variant (1) encodes the longer isoform (1).</p>