

## Product datasheet for **SC317807**

### DDX6 (NM\_004397) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	DDX6 (NM_004397) Human Untagged Clone
Tag:	Tag Free
Symbol:	DDX6
Synonyms:	HLR2; IDDILF; P54; RCK; Rck/p54
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**Fully Sequenced ORF:**

>SC317807 representing NM\_004397.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGAGCACGGCCAGAACAGAGAACCCTGTTATAATGGTCTGTCCAGTCAAATGGTCAGCTGAGAGGC
CCTGTGAAACCCACTGGTGGCCCTGGAGGAGGGGCACACAGACACAGCAACAGATGAACCAGCTGAAA
AACACCAACACAATCAATAATGGCACTCAGCAGCAAGCACAGAGTATGACCACCCTATTAACCTGGT
GATGACTGGAAAAAGACTTTAAACTCCCTCCAAAGGATCTAAGAATCAAACCTCGGATGTGACCTCC
ACAAAAGGAAATGAGTTTGAAGATTACTGTTTGAACGGGAGTTACTGATGGAAATTTTGAATGGGC
TGGGAAAAGCCATCTCCTATTCAGGAGGAGAGCATTCCCATTTGCTTTATCTGGTAGGGATATCTTAGCT
AGAGCAAAAAATGGAACAGGCAAGAGCGGTGCCTACCTCATTCCCTTACTGAACGGCTAGACCTGAAG
AAGGACAATATAACAAGCAATGGTATTGTTCCCACTAGAGAACTGCTCTACAGGTCAGTCAAATTTGC
ATCCAGGTCAGCAAACACATGGGAGGGCCAAAGTGATGGCAACCACAGGAGGAACCAATTTACGAGAT
GACATAATGAGGCTTGATGATACAGTGCACGTGGTATTGCTACCCCTGGGAGAATCCTGGATCTTATT
AAGAAAGGAGTAGCAAAGTTGATCATGTCCAGATGATAGTATTGGATGAGGCAGATAAGTTGCTGTCA
CAGGATTTTGTGCAGATAATGGAGGATATTATTCTCACGCTACCTAAAAACAGGCAGATTTTACTATAT
TCCGCTACTTTCCCTCTTAGTGTACAGAAGTTCATGAATCCCATTTGCAGAAACCCTATGAGATTAAC
CTGATGGAGGAACTAECTCTGAAGGGAGTAACCCAGTACTACGCATATGTAECTGAGCGCAAAAAAGTA
CACTGCCTCAACACACTTTTCTCCAGGCTTCAGATAAACCCAGTCGATCATTTTCTGTAECTCTCAG
CGAGTTGAATTGCTAGCCAAGAAGATTTCTCACTGGGTATTCTTGCTTCTATATTCATGCTAAAATG
AGGCAGGAACATCGAAATCGTGTATTTTCATGATTTCCGAAATGGCTTATGCCGCAATCTTGTTGCACT
GATCTGTTTACCCGAGGTATTGATATAACAAGCTGTAATGTGGTAATAAACTTTGATTTCCCAAAGCTG
GCAGAGACCTATCTCCATCGTATTGGAAGATCAGGTCGCTTTGGTCATCTTGCTTAGCCATCAACTTG
ATCACATATGATGATCGTTCAACCTGAAAAGTATTGAGGAGCAGCTGGGAACAGAAATTAACCTATT
CCGAGCAACATTGATAAGAGCCTGTATGTGGCAGAATACCACAGCGAGCCTGTAGAAGATGAGAAACCT
TAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
    
```

**Restriction Sites:**

SgfI-MluI

**Plasmid Map:**



**ACCN:**

NM\_004397

<b>Insert Size:</b>	1452 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_004397.5</a>
<b>RefSeq Size:</b>	6246 bp
<b>RefSeq ORF:</b>	1452 bp
<b>Locus ID:</b>	1656
<b>UniProt ID:</b>	<a href="#">P26196</a>
<b>Cytogenetics:</b>	11q23.3
<b>Domains:</b>	DEAD, helicase_C
<b>Protein Pathways:</b>	RNA degradation
<b>MW:</b>	54.4 kDa
<b>Gene Summary:</b>	This gene encodes a member of the DEAD box protein family. The protein is an RNA helicase found in P-bodies and stress granules, and functions in translation suppression and mRNA degradation. It is required for microRNA-induced gene silencing. Multiple alternatively spliced variants, encoding the same protein, have been identified. [provided by RefSeq, Mar 2012] Transcript Variant: This variant (1) represents the longer transcript. Variants 1 and 2 encode the same protein.