

## Product datasheet for **SC317809**

### THNSL2 (NM\_018271) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	THNSL2 (NM_018271) Human Untagged Clone
Tag:	Tag Free
Symbol:	THNSL2
Synonyms:	SOFAT; THS2; TSH2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**Fully Sequenced ORF:** >SC317809 representing NM\_018271.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGTGGTATGTCAGCACCAGGGGCGTAGCCCCACGGGTCAACTTTGAGGGGGCCCTTCTCTGCGCTAT
GCACCTGACGGGGCCCTTTTATGCCTGAAGAGCTCCACAGTTGGACAGAGGGACCCTGTGCCAGTGG
AGCACACTCTCCTATCCTGGCCTGGTGAAGGAGCTGTGTGCCCTTTCATTGGCTCTGAGCTCCTCCA
AAAGATGAATTAATGATCTGATCGACCAGCCTTCAGCAGATCCGTCACAGAGAAGTGGTCCATCTG
TCCAGGTTGAGGAATGGGCTGAACGTGTTGGAGCTGTGGCATGGCGTCACATATGCATTTAAGGACCTG
TCCCTGTCCTGCACAACACAGTTCCTGCAGTACTTCTGGAGAAGAGGGAGAAGCACGTCAGTGTGGTT
GTAGGAACATCTGGGACACAGGAAGTGTGCCATTGAGAGTGTCAAGGGGCAAAGAACATGGACATT
ATCGTTCTGCTGCCCAAAGGTCAGTGCACAAGATTCCAGAGCTCCAGATGACAACGGTGTGAAGCAG
AACGTACATGTGTTGGAGTGGAGGAAACAGCGATGAGCTCGATGAGCCGATCAAGACTGTGTTTGCC
GATGTGGCTTTTGTCAAGAAGCACAATCTGATGAGCCTGAATTCGATCAACTGGTCCCGGGTCTGGTG
CAGATGGCCCATCACTTCTTTGCTACTTCCAGGTACGCCATCCTTGGACACACATCCCCTACCCCTG
GTGGAGGTGGTTGTGCCAACAGGGGCTGCCGTAACCTTGCAGCTGGGTACATTGCTCAAAAAGTAGGC
CTGCCATCCGTCTGGTCTGGCAGTGAACCGCAATGACATCATCCACAGGACTGTCCAGCAGGGAGAC
TTCTCTCTCTGAGGCTGTTAAATCAACCTTGGCATCAGCTATGGACATTCCAGTGCCTACAACATG
GAGAGGGTGTCTGGCTGCTCTGCTGCTGACAGCCAGGTGACAAGAGCCCTCATGGAGCAGTTTGAA
AGGACCCAAAGTGAATCTGCCAAGGAAGTGCACAGCAAGCTTTCAGAGGCAGTGACATCCGTGTCA
GTGTCGGATGAAGCCATCACCCAGACCATGGGCCGCTGCTGGGATGAGAACCAGTACTTGTGTGCCCC
CACTCAGCGGTGGCCGTGAAGTACCATTACCAGCAGATAGACAGGCAGCAGCCAGCAGCCTCCCGGTGC
TGCTCGCCCTGCCTCTGCAGCAAGTTCCCGAAGCTGTCTGGCTGCTGGCCTGACCCCTGAGACT
CCCGCGGAGATCGTAGCCCTGGAGCACAAGGAGACACGCTGCACCCTGATGCGGAGAGGTGACAACCTGG
ATGCTGATGCTTCGGGACACCATTGAGGACCTTAGCCGACAGTGGAGGAGTCATGCCCTCAACACCTCC
CAGTAG
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

**Restriction Sites:**

SgfI-MluI

**Plasmid Map:**



**ACCN:**

NM\_018271

<b>Insert Size:</b>	1455 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<u><a href="#">NM_018271.4</a></u>
<b>RefSeq Size:</b>	1985 bp
<b>RefSeq ORF:</b>	1455 bp
<b>Locus ID:</b>	55258
<b>UniProt ID:</b>	<u><a href="#">Q86YJ6</a></u>
<b>Cytogenetics:</b>	2p11.2
<b>MW:</b>	54.1 kDa
<b>Gene Summary:</b>	<p>This gene encodes a threonine synthase-like protein. A similar enzyme in mouse can catalyze the degradation of O-phospho-homoserine to a-ketobutyrate, phosphate, and ammonia. This protein also has phospho-lyase activity on both gamma and beta phosphorylated substrates. In mouse an alternatively spliced form of this protein has been shown to act as a cytokine and can induce the production of the inflammatory cytokine IL6 in osteoblasts. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Sep 2011]</p> <p>Transcript Variant: This variant (1) encodes the longest isoform (1). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>