

## Product datasheet for **SC317811**

### Aspartyl Aminopeptidase (DNPEP) (NM\_012100) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Aspartyl Aminopeptidase (DNPEP) (NM_012100) Human Untagged Clone
Tag:	Tag Free
Symbol:	DNPEP
Synonyms:	ASPEP; DAP
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)


**Fully Sequenced ORF:** >SC317811 representing NM\_012100.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGAGCGGACACAGCCCCACGCGGGCCATGCAGGTGGCCATGAACGGTAAGGCCCGCAAAGAGGCG
GTGCAGACTGCGGCTAAGGAACCTCAAGTTCGTGAACCGGAGTCCCTCTCCTTTCCATGCTGTGGCT
GAATGCCGCAACCGCCTTCTCCAGGCTGGCTTCACTGAACCTCAAGGAGACTGAGAAATGGAATATTAAG
CCCAGAGCAAGTACTTCATGACCAGGAACTCCTCCACCATCATAGCTTTTGCTGTAGGGGGCCAGTAC
GTTCTGGCAATGGCTTCAAGCTCATCGGGGCCACACGGACAGCCCCTGCCTCCGGGTGAAACGTCGG
TCTCGCCGAGCCAGGTGGGCTTCCAGCAAGTCGGTGTGGAGACCTATGGTGGTGGGATCTGGAGCACC
TGGTTTGACCGTGACCTGACTCTGGCTGGACGCGTCATTGTCAAGTGCCTACCTCAGGTGCGCTGGAG
CAGCAGCTGGTGCACGTGGAGCGGCCATTCTTCGCATCCACACCTGGCCATCCATCTGCAGCGAAAT
ATCAACGAGAACCTTTGGGCCAACACAGAGATGCATCTAGTCCCATTCTTGCCACAGCCATCCAGGAG
GAGCTGGAGAAGGGGACTCCTGAGCCAGGGCCTCTCAATGCTGTGGATGAGCGGCACCATTCCGGTCTC
ATGTCCCTGCTCTGTGCCATCTGGGGCTGAGCCCAAGGACATAGTGGAGATGGAGCTCTGCCTTGCA
GACACCCAGCCTGCGGTCTTGGTGGTGCCTATGATGAGTTCATCTTTGCTCCTCGGCTGGACAATCTG
CACAGTGCTTCTGTGCCCTGCAGGCCTTGATAGATTCTGTGCAGGCCCTGGTCCCTGGCCACAGAG
CCTCACGTGCGCATGGTCACACTCTATGACAACGAAGAGGTGGGGTCTGAGAGTGACAGGGAGCACAG
TCACTGCTGACAGAGCTGGTGTGCGGGCGATCTCAGCCTCGTCCAGCACCCGACAGCCTTCGAGGAA
GCCATACCAAGTCTTCATGATCAGCGCAGACATGGCCATGCTGTGCATCCCAACTACCTGGACAAG
CATGAGGAGAACCACCGCCTTTATTCCACAAGGGCCCGTGATCAAGGTGAACAGCAACGCAACGCTAT
GCTTCAAACGCGGTGTCAGAGGCCCTGATCCGAGAGGTGGCCAACAAAGTCAAGTCCCCCTGCAGGAT
CTCATGGTCCGGAATGACACCCCTGTGGAACCACTTGGACCTATCTGGCTTCTCGGCTGGGGCTG
CGGTGCTGGATTTAGGAGCCCAACTGGCCATGCACTCTATCCGGGAGATGGCTGCACCACAGGA
GTCCTCCAGACCCTCACCTCTCAAGGCTTCTTTGAGCTGTTCCCTTCTAAGCCATAATCTCTTA
GTGGATTGA
AGCGGACCGACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGAT
ATCCTGGATTACAAGGATGACGACGATAAGGTTTAA
  
```

**Restriction Sites:** SgfI-RsrII

**Plasmid Map:**



**ACCN:** NM\_012100

<b>Insert Size:</b>	1458 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<u><a href="#">NM_012100.3</a></u>
<b>RefSeq Size:</b>	3884 bp
<b>RefSeq ORF:</b>	1458 bp
<b>Locus ID:</b>	23549
<b>UniProt ID:</b>	<u><a href="#">Q9ULA0</a></u>
<b>Cytogenetics:</b>	2q35
<b>Domains:</b>	Peptidase_M18
<b>Protein Families:</b>	Druggable Genome, Protease
<b>MW:</b>	53.4 kDa
<b>Gene Summary:</b>	<p>The protein encoded by this gene is an aminopeptidase which prefers acidic amino acids, and specifically favors aspartic acid over glutamic acid. It is thought to be a cytosolic protein involved in general metabolism of intracellular proteins. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jan 2016]</p> <p>Transcript Variant: This variant (3) differs in the 5' UTR and coding sequence compared to variant 1. The resulting isoform (c) has a shorter and distinct N-terminus compared to isoform a. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>