

Product datasheet for **SC317816**

LMBR1L (NM_018113) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	LMBR1L (NM_018113) Human Untagged Clone
Tag:	Tag Free
Symbol:	LMBR1L
Synonyms:	LIMR
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC317816 representing NM_018113.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGAAGCACCTGACTACGAAGTGCTATCCGTGCGAGAACAGCTATTCCACGAGAGGATCCGCGAGTGT
ATTATATCAACACTTCTGTTTGCAACACTGTACATCCTCTGCCACATCTTCTGACCCGCTTCAAGAAG
CCTGCTGAGTTCACCACAGTGGATGATGAAGATGCCACCGTCAACAAGATTGCGCTCGAGCTGTGCACC
TTTACCTTGGCAATTGCCCTGGGTGCTGTCTCTGCTCTTCTCCATCATCAGCAATGAGGTGCTG
CTCTCCCTGCCTCGGAACCTACTACATCCAGTGGCTCAACGGCTCCCTCATCCATGGCCTCTGGAACCTT
GTTTTCTCTTCTCCAACTGTCCCTCATCTTCTCATGCCCTTTCATATTTCTTCACTGAGTCTGAG
GGCTTTGCTGGCTCCAGAAAGGTGTCTGGGCCGGGTCTATGAGACAGTGGTGTGTTGATGCTCCTC
ACTCTGCTGGTGCTAGGTATGGTGTGGGTGGCATCAGCCATTGTGGACAAGAACAGGCACAGAGAG
TCACTCTATGACTTTTGGGAGTACTATCTCCCTACCTCTACTCATGCATCTCCTTCTTGGGGTTCTG
CTGCTCCTGGTGTGTACTCCACTGGGTCTCGCCCGCATGTTCTCCGTCACTGGGAAGCTGCTAGTCAAG
CCCCGGCTGCTGGAAGACCTGGAGGAGCAGCTGTAAGTCTCAGCCTTTGAGGAGGCAGCCCTGACCCGC
AGGATCTGTAATCCTACTTCTGCTGGCTGCCTTTAGACATGGAGCTGCTACACAGACAGGTCTGGCT
CTGCAGACACAGAGGTCTGCTGGAGAAGAGGCGGAAGGCTTCAGCCTGGCAACGGAACCTGGGCTAC
CCCCGGCTATGCTGTGCTTGTGGTGTGCTGACGGGCTGTCTGTGCTCATTGTGGCCATCCACATCTG
GAGCTGCTCATCGATGAGGCTGCCATGCCCGAGGCATGCAGGTACCTCCTTAGGCCAGGTCTCCTTC
TCCAAGCTGGGCTCCTTTGGTGCCGTATTACAGTTGTACTCATCTTTACCTAATGGTGTCTCAGTT
GTGGGCTTCTATAGCTCTCCACTCTTCCGGAGCCTGCGGCCAGATGGCACGACACTGCCATGACGCAG
ATAATTGGGAACGTGTGCTGTCTCTGGTCTTAAGCTCAGCACTTCTGTCTTCTCTCGAACCTGGGG
CTCACTCGCTTTGACCTGCTGGGTGACTTTGGACGCTTCAACTGGCTGGGAATTTCTACATTGTGTTT
CTCTACAACGCAGCCTTTGCAGGCCTCACCACACTCTGTCTGGTGAAGACCTTCACTGCAGCTGTGCGG
GCAGAGCTGATCCGGGCCTTTGGGCTGGACAGACTGCCGCTGCCGCTCTCCGGTTTCCCCAGGCATCT
AGGAAGCCACGACCCAGTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC

Restriction Sites: SgfI-MluI

ACCN: NM 018113

Insert Size: 1470 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_018113.3</u>
RefSeq Size:	2361 bp
RefSeq ORF:	1470 bp
Locus ID:	55716
UniProt ID:	<u>Q6UX01</u>
Cytogenetics:	12q13.12
Domains:	LMBR1
Protein Families:	Druggable Genome, Transmembrane
MW:	55.2 kDa
Gene Summary:	Probable LCN1 receptor. May mediate LCN1 endocytosis.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (1) encodes the longest isoform (1).