

Product datasheet for **SC317819**

CYP2W1 (NM_017781) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CYP2W1 (NM_017781) Human Untagged Clone
Tag:	Tag Free
Symbol:	CYP2W1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)


Fully Sequenced ORF: >SC317819 representing NM_017781.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGCCCTGCTGCTTTGCTGTTCTGGCCCTCTGGGGCTCTGGGGCTGCTCTGCGCCTGCGCCCAA
GACCCCTCCCAGCTGCCGGTGGCCCCGGGGCTCGCCCGCTGCCGCTCGTGGGAACCTGCACCTTG
CTGCGTCTGTCGCAACAGGACCGGTCCCTGATGGAGCTCTCAGAACGCTACGGGCCGGTGTTCACCGTG
CACCTGGGGCGCCAGAAGACGGTGGTGTGACGGGTTTCAGGCGGTCAAAGAGGCGCTGGCGGGCCCC
GGGCAGGAGCTGGCCGACCGGCCTCCCATCGCCATCTTCCAGCTCATCCAGCGAGGTGGAGGCATCTTC
TTCTCATCTGGGCGCGCTGGAGGGTGGCCGCCAGTTCACGGTGCCTGCCCTGCACAGCCTGGCGGTG
GGCCGGGAGCCGGTGGCTGACAAGATTCTGCAGGAGCTGAAATGCCTCTCTGGCAGCTGGATGGCTAC
AGAGGCCGGCCCTTCCCCTGGCCCTACTGGCTGGGCTCCCTCAATATCACCTTCGCGCTCCTCTTC
GGCCGCCGATTTGACTACCGGGACCCCGTGTGGTGTCCCTGCTGGGTCTCATCGATGAGGTGATGGTC
CTTTGGGGTCCCCTGGCCTGCAGCTGTCAACGTCTACCATGGCTCGGGGCCCTGCTCCAGCTGCAC
CGCCCGTCTGCGCAAGATCGAGGAGTCCGTGCCATTCTGAGGACCCTCTGGAGGCGCGGAGGCC
CACGTGTGCCCGGGGACCCCGTGTGCAGCTATGTGGACGCCCTGATCCAGCAGGACAGGGGGATGAC
CCCAGGGCCTGTTTGTGAGGCCAACGCGGTGGCCTGCACCCTGGACATGGTATGGCCGGGACGGAG
ACGACCTCGGCCACGCTGCAGTGGGCCGCACTTCTGATGGGCCGGCACCCGGACGTGCAGGGCCGGGTG
CAGGAGGAGCTAGACCGGTGCTGGGCCCTGGGCGGACTCCCCGGCTGGAGGACCAGCAGGCTGCCCC
TACACAAGCGCGTGTCCACGAGGTGCAGCGTTCATCACGCTCCTGCCGACGTGCCCCGCTGCACC
GGCGCCGACACACAGCTGGCGGCTTCTGTCCCCAAGGGACGCCCCGTGATCCCTGCTGACCTCG
GTGCTCCTGGATGAGACACAGTGGCAGACCCAGCCAGTTC AACCCCGGCCATTTCTGGAGCGAAT
GGGCACTTTGTGAAGCGGGAGGCCCTTCTGCCTTCTCTGCAGGCCCGCGCTGTGTGGGGAGCGC
CTGGCCAGGACCGAGCTTCTCTGCTGTTTGGCGCCTCTGCAGAGGTACCGCCTGCTGCCCCCGCT
GGCGTCAGTCCGGCTCCCTGGACACCACGCCCGCCGGCTTTTACCATGAGGCCGAGGGCCAGGCC
CTGTGTGCGGTGCCAGGCCTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_017781

Insert Size: 1473 bp

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_017781.2
RefSeq Size:	2354 bp
RefSeq ORF:	1473 bp
Locus ID:	54905
UniProt ID:	Q8TAV3
Cytogenetics:	7p22.3
Protein Families:	Druggable Genome, P450
MW:	53.8 kDa
Gene Summary:	This gene encodes a member of the cytochrome P450 superfamily of enzymes. The cytochrome P450 proteins are monooxygenases which catalyze many reactions involved in drug metabolism and synthesis of cholesterol, steroids and other lipids. [provided by RefSeq, Jul 2008]