

Product datasheet for **SC317821**

PLD3 (NM_012268) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PLD3 (NM_012268) Human Untagged Clone
Tag:	Tag Free
Symbol:	PLD3
Synonyms:	AD19; HU-K4; HUK4; SCA46
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)


Fully Sequenced ORF: >SC317821 representing NM_012268.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGAAGCCTAAACTGATGTACCAGGAGCTGAAGGTGCCTGCAGAGGAGCCCGCCAATGAGCTGCCCATG
AATGAGATTGAGCGTGGAAAGGCTGCGGAAAAGAAAGCCCGCTGGGTCTGCTGGTCTCATTCTGGCG
GTTGTGGGCTTCGGAGCCCTGATGACTCAGCTGTTTCTATGGGAATACGGCGACTTGCATCTCTTGGG
CCCAACCAGCGCCAGCCCTGCTATGACCCTTGCGAAGCAGTCTGGTGGAAAGCATTCTGAGGGC
CTGGACTTCCCAATGCCTCCACGGGAACCTTCCACCAGCCAGGCTGGCTGGGCTGCTCGCCGGT
GCGCACAGCAGCCTGGACATCGCTCTTCTACTGGACCCTACCAACAATGACACCCACACGCAGGAG
CCCTCTGCCAGCAGGGTGAAGGCTCCTCCGGCAGCTGCAGACCCTGGCACAAAGGGCGTGAACGTC
CGCATCGTGTGAGCAAGCCAGCGGGCCAGCCACAGCGGACCTGCAGGCTCTGCTGCAGAGCGGT
GCCAGGTCCGCATGGTGGACATGCAGAAGCTGACCCATGGCGTCTGCATACCAAGTTCTGGTGGTG
GACCAGCCACTTCTACCTGGGAGTGCACAATGGACTGGCGTTCCTGACCCAGGTCAAGGAGCTG
GGCGTGGTCATGTACAACCTGCAGCTGCCTGGCTCGAGACCTGACCAAGATCTTTGAGGCCTACTGGTTC
CTGGGCCAGGCAGGCAGCTCCATCCCATCAACTTGGCCCCGGTCTATGACACCCGCTACAACCAAGAG
ACACCAATGGAGATCTGCCTCAATGGAACCCCTGCTCTGGCTACCTGGCGAGTGCGCCCCACCCCTG
TGTCCAAGTGGCCGCACTCCAGACCTGAAGGCTCTACTCAACGTGGTGGACAATGCCCGAGTTTCATC
TACGTCGCTGTATGAATACCTGCCACTCTGGAGTTCTCCACCCTCACAGGTTCTGGCTGCCATT
GACGATGGGCTGCGGGCGGCCACCTACGAGCGTGGCGTCAAGTGGCCTGCTCATCAGCTGCTGGGA
CACTCGGAGCCATCCATGCGGGCCTTCTGCTCTCTGGCTGCCCTGCGTGACAACCATACCCACTCT
GACATCCAGTGAAACTCTTTGTGGTCCCGCGGATGAGGCCAGGCTCGAATCCCATATGCCCGTGC
AACCAACAAGTACATGGTGACTGAACGCGCCACCTACATCGGAACCTCCAACCTGGTCTGGCAACTAC
TTCACGGAGACGGCGGCACCTCGCTGCTGGTGCAGCAGAATGGGAGGGCGGCCTGCGGAGCCAGCTG
GAGGCCATTTCTGAGGGACTGGGACTCCCTTACAGCCATGACCTTGACACCTCAGCTGACAGCGTG
GGCAACGCCTGCCGCTGCTTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_012268

Insert Size:	1473 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_012268.3
RefSeq Size:	2431 bp
RefSeq ORF:	1473 bp
Locus ID:	23646
UniProt ID:	Q8IV08
Cytogenetics:	19q13.2
Domains:	PLDc
Protein Families:	Transmembrane
MW:	54.7 kDa
Gene Summary:	<p>This gene encodes a member of the phospholipase D (PLD) family of enzymes that catalyze the hydrolysis of membrane phospholipids. The encoded protein is a single-pass type II membrane protein and contains two PLD phosphodiesterase domains. This protein influences processing of amyloid-beta precursor protein. Mutations in this gene are associated with Alzheimer disease risk. Alternatively spliced transcript variants encoding the same protein have been found for this gene. [provided by RefSeq, Apr 2014]</p> <p>Transcript Variant: This variant (2) represents the longer transcript. Variants 1, 2 and 3 encode the same protein.</p>