

Product datasheet for **SC317822**

SLC43A3 (NM_014096) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SLC43A3 (NM_014096) Human Untagged Clone
Tag:	Tag Free
Symbol:	SLC43A3
Synonyms:	EEG1; ENBT1; FOAP-13; PRO1659; SEEEG-1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC317822 representing NM_014096.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCGGGCCAGGGCCTGCCCTGCACGTGGCCACACTGCTGACTGGGCTGCTGGAATGCCTGGGCTTT
GCTGGCGTCCTCTTTGGCTGGCCTTCACTAGTGTCTTCAAGAATGAAGATTACTTTAAGGATCTG
TGTGGACCAGATGCTGGGCCGATTGGCAATGCCACAGGGCAGGCTGACTGCAAAGCCAGGATGAGAGG
TTCTCACTCATCTTCACCTGGGGTCTTCATGAACAACCTTCATGACATTCCTCCACTGGCTACATCTTT
GACCGGTTCAAGACCACCGTGGCAGCCTCATAGCCATATTTTTCTACACCACCGCCACACTCATCATA
GCCTTACCTCTGCAGGCTCAGCCGTGCTGCTCTTCTGGCCATGCCAATGCTCACCATTGGGGGAATC
CTGTTTCTCATCACCAACCTGCAGATTGGGAACCTATTTGGCCAACACCGTTCGACCATCATCACTCTG
TACAATGGAGCATTTGACTCTTCTCGGCAGTCTTCTTATTATTAAGCTTCTTTATGAAAAGGCATC
AGCCTCAGGGCCTCTTCATCTTCTGTCTGCAGTACCTGGCATGTAGCACGCACTTCTCTCTG
ATGCCCCGGGGCACATCCCATACCCACTGCCCCCAACTACAGCTATGGCCTGTGCCCTGGGAATGGC
ACCACAAAGGAAGAGAAGGAAACAGCTGAGCATGAAAACAGGGAGCTACAGTCAAAGGAGTTCCTTTCA
GCGAAGGAAGAGACCCAGGGGCAGGGCAGAAGCAGGAACCTCCGCTCCTTCTGGAGCTACGCTTTCTCT
CGGCGCTTTGCTGGCACCTGGTGTGGCTGTCTGTGATACAGTTGTGGCACTACCTTTCATTGGCACT
CTCAACTCCTTGCTGACCAACATGGCCGGTGGGGACATGGCAGAGTCAACACCTACACAAATGCCTTT
GCCTTCACTCAGTTCGGAGTGTGTGTGCCCTTGGAAATGGCCTGCTCATGGACCGGCTTAAACAGAAG
TACCAGAAGGAAGCAAGAAAGACAGGTTCTCCACTTTGGCGTGGCCCTCTGCTCGACGGTGCCTTCG
CTGGCCCTGACATCCCTGCTGTGCTGGGCTTCGCCCTCTGTGCCTCAGTCCCCATCTCCCTCTCCAG
TACCTCACCTTTCATCCTGCAAGTATCAGCCGCTCCTTCTCTATGGGAGCAACGGCCCTTCTCACC
CTTGCTTTCCCTCAGAGCACTTTGGCAAGCTCTTTGGGCTGGTATGGCCTTGTGGCTGTGGTGTCT
CTGCTCCAGTTCCCATCTTACCTCATCAAAGGCTCCCTTCAAGATGACCCATTTTACGTAATGTG
ATGTTTATGCTTGCCATTCTTCTGACATTTCTCCACCCCTTTCTGGTATATCGGGAATGCCGTAATTGG
AAAGAAAGTCCCTCTGCAATTGCATAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_014096

Insert Size: 1476 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_014096.3](#)

RefSeq Size: 2904 bp

RefSeq ORF: 1476 bp

Locus ID: 29015

UniProt ID: [Q8NBI5](#)

Cytogenetics: 11q12.1

Protein Families: Transmembrane

MW: 54.5 kDa

Gene Summary: Putative transporter.[UniProtKB/Swiss-Prot Function]

Transcript Variant: This variant (1) and variants 2, 3 and 4 encode the same isoform (1).