

Product datasheet for **SC317838**

TRIM22 (NM_006074) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TRIM22 (NM_006074) Human Untagged Clone
Tag:	Tag Free
Symbol:	TRIM22
Synonyms:	GPSTAF50; RNF94; STAF50
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)


Fully Sequenced ORF: >SC317838 representing NM_006074.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGATTTCTCAGTAAAGGTAGACATAGAGAAGGAGGTGACCTGCCCATCTGCCTGGAGCTCTGACA
GAACCTTGAGCCTAGATTGTGGCCACAGCTTCTGCCAAGCCTGCATCACTGCAAAGATCAAGGAGTCA
GTGATCATCTCAAGAGGGAAAGCAGCTGTCTGTGTGTCAGACCAGATTCCAGCCTGGGAACCTCCGA
CCTAATCGGCATCTGGCCAACATAGTTGAGAGAGTCAAAGAGGTCAAGATGAGCCACAGGAGGGGCAG
AAGAGAGATGTCTGTGAGCACCATGAAAAAACTCCAGATCTTCTGTAAGGAGGATGAAAAAGTCATT
TGCTGGGTTTGTAACTGTCTCAGGAACACCAAGGTCACCAAACATTCCGCATAAACGAGGTGGTCAAG
GAATGTCAGGAAAAGCTGCAGGTAGCCCTGCAGAGGCTGATAAAGGAGGATCAAGAGGCTGAGAAGCTG
GAAGATGACATCAGACAAGAGAGAACCCTGGAAGAATTATATCCAGATCGAGAGACAGAAGATTCTG
AAAGGTTCAATGAAATGAGAGTCATCTTGACAATGAGGAGCAGAGAGAGCTGCAAAGCTGGAGGAA
GGTGAGGTGAATGTGCTGGATAACCTGGCAGCAGCTACAGACCAGCTGGTCCAGCAGAGGCAGGATGCC
AGCAGCTCATCTCAGATCTCCAGCGGAGTTGAGGGGATCGTCACTAGAGATGCTGCAGGATGTGATT
GACGTCATGAAAAGGAGTGAAAGCTGGACATTGAAGAAGCAAAAATCTGTTTCCAAGAACTAAAGAGT
GTATTCGAGTACCAGATCTGAGTGGGATGCTGCAAGTCTTAAAGAGCTGACAGATGTCCAGTACTAC
TGGGTGGACGTGATGCTGAATCCAGGCAGTGCCACTTCGAATGTTGCTATTTCTGTGGATCAGAGACAA
GTGAAAAGTGTACGCACCTGCACATTTAAGAATCAAATCCATGTGATTTTTCTGCTTTTGGTGTCTTC
GGCTGCCAATATTTCTCTCGGGGAAATATTACTGGGAAGTAGATGTGTCTGGAAAGATTGCCTGGATC
CTGGGCGTACACAGTAAAATAAGTAGTCTGAATAAAAAGGAAGAGCTCTGGGTTTGCTTTTGATCCAAGT
GTAAATTTATCAAAGTTTACTCCAGATATAGACCTCAATATGGCTACTGGGTTATAGGATTACAGAA
ACATGTGAATATAATGCTTTTGAGGACTCTCCTCTTCTGATCCCAAGGTTTTGACTCTCTTTATGGCT
GTGCTCCTGTGCTATTGGGTTTTCTAGACTATGAGGCAGGATTGTCTATTTTTCAATGTACA
AACCACGGAGCACTCATCTACAAGTTCTCTGGATGTCGCTTTTCTCGACCTGCTTATCCGATTTTCAAT
CCTTGGAACTGCCTAGTCCCATGACTGTGTGCCACCGAGCTCCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_006074

Insert Size:	1497 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_006074.4
RefSeq Size:	3111 bp
RefSeq ORF:	1497 bp
Locus ID:	10346
UniProt ID:	Q8IYM9
Cytogenetics:	11p15.4
Domains:	zf-B_box, RING
Protein Families:	Druggable Genome, Transcription Factors
MW:	56.9 kDa
Gene Summary:	<p>The protein encoded by this gene is a member of the tripartite motif (TRIM) family. The TRIM motif includes three zinc-binding domains, a RING, a B-box type 1 and a B-box type 2, and a coiled-coil region. This protein localizes to the cytoplasm and its expression is induced by interferon. The protein down-regulates transcription from the HIV-1 LTR promoter region, suggesting that function of this protein may be to mediate interferon's antiviral effects. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Dec 2010]</p> <p>Transcript Variant: This variant (1) represents the longer transcript and encodes the longer isoform (1).</p>