

## Product datasheet for **SC317839**

### **SIGLEC8 (NM\_014442) Human Untagged Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	SIGLEC8 (NM_014442) Human Untagged Clone
Tag:	Tag Free
Symbol:	SIGLEC8
Synonyms:	SAF2; SIGLEC-8; SIGLEC8L
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)


**Fully Sequenced ORF:** >SC317839 representing NM\_014442.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGCTGCTGCTGCTGCTGCTGCTGCCCTGCTCTGGGGGACAAAGGGGATGGAGGGAGACAGACAATAT
GGGGATGGTTACTTGTGCAAGTGCAGGAGCTGGTGACGGTGCAGGAGGGCCTGTGTGCCATGTGCC
TGCTCCTTCTCCTACCCCAAGGATGGCTGGACTGACTCTGACCCAGTTTCATGGCTACTGGTTCCGGGCA
GGAGACAGACCATAACCAAGACGCTCCAGTGGCCACAAACAACCCAGACAGAGAAGTGCAGGCAGAGACC
CAGGGCCGATTCCAACCTCCTTGGGGACATTTGGAGCAACGACTGCTCCCTGAGCATCAGAGACGCCAGG
AAGAGGGATAAAGGGGTCAATTTCTTTTCGGCTAGAGAGAGGAAGCATGAAATGGAGTTACAAATCACAG
TTGAATTACAAAATAAGCAGCTGTCTGTGTTTGTGACAGCCCTGACCCATAGGCCTGACATCCTCATC
CTAGGGACCCTAGAGTCTGGCCACTCCAGAACCTGACCTGCTCTGTGCCCTGGGCTGTAAGCAGGGG
ACACCCCCATGATCTCCTGGATTGGGGCTCCGTGTCCTCCCCGGGCCCACTACTGCCCGCTCCTCA
GTGCTCACCTTACCCCAAAGCCCCAGGACCACGGCACCAGCCTCACCTGTCAGGTGACCTTGCCCTGGG
ACAGGTGTGACCAGCAGCAGTACCGTCCGCCTCGATGTGTCTACCTCCTTGGAACTTGACCATGACT
GTCTTCCAAGGAGATGCCACAGCATCCACAGCCCTGGGAAATGGCTCATCTCTTTAGTCCTTGAGGGC
CAGTCTCTGCGCCTGGTCTGTGCTGTCAACAGCAATCCCCTGCCAGGCTGAGCTGGACCCGGGGGAGC
CTGACCTGTGCCCTCACGGTCTCAAACCTGGGCTGCTGGAGCTGCCTCGAGTGCACGTGAGGGAT
GAAGGGGAATTCACCTGCCGAGCTCAGAACGCTCAGGGCTCCCAGCACATTTCCCTGAGCCTCTCCCTG
CAGAATGAGGGCAGGCACCTCAAGACCTGTATCACAAGTGACACTGGCAGCAGTCGGGGGAGCTGGA
GCCACAGCCCTGGCCTTCTGTCTTCTGCATCATCTTCATCATAGTGAGGTCTGCAGGAAGAAATCG
GCAAGGCCAGCAGCGGGCTGGGGGATACAGGCATGGAAGATGCAAAGGCCATCAGGGGCTCGGCCTCT
CAGGGACCCCTGACTGAATCCTGGAAAGATGGCAACCCCTGAAGAAGCCTCCCCAGCTGTTGCCCC
TCGTCAAGGGAGGAAGGAGAGCTCCATTATGCAACCTCAGTTCATATAAGTGAAGCCTCAGGACCCG
CAGGGACAGGAGGCCACTGACAGTGAATACTCGGAGATCAAGATCCACAAGCGAGAAACTGCAGAGACT
CAGGCCTGTTGAGGAATCACAACCCCTCCAGCAAAGAAGTCAGAGGCTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

**Restriction Sites:** SgfI-MluI

**Plasmid Map:**



**ACCN:** NM\_014442

<b>Insert Size:</b>	1500 bp
<b>OTI Disclaimer:</b>	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at <a href="mailto:custsupport@origene.com">custsupport@origene.com</a> or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a></p>
<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<a href="#">NM_014442.2</a>
<b>RefSeq Size:</b>	2967 bp
<b>RefSeq ORF:</b>	1500 bp
<b>Locus ID:</b>	27181
<b>UniProt ID:</b>	<a href="#">Q9NYZ4</a>
<b>Cytogenetics:</b>	19q13.41
<b>Protein Families:</b>	Druggable Genome, Transmembrane
<b>MW:</b>	54 kDa

**Gene Summary:**

Sialic acid-binding immunoglobulin (Ig)-like lectins, or SIGLECs (e.g., CD33 (MIM 159590)), are a family of type 1 transmembrane proteins each having a unique expression pattern, mostly in hemopoietic cells. SIGLEC8 is a member of the CD33-like subgroup of SIGLECs, which are localized to 19q13.3-q13.4 and have 2 conserved cytoplasmic tyrosine-based motifs: an immunoreceptor tyrosine-based inhibitory motif, or ITIM (see MIM 604964), and a motif homologous to one identified in signaling lymphocyte activation molecule (SLAM; MIM 603492) that mediates an association with SLAM-associated protein (SAP; MIM 300490) (summarized by Foussias et al., 2000 [PubMed 11095983]).[supplied by OMIM, May 2010]