

Product datasheet for **SC317840**

ZNF253 (NM_021047) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF253 (NM_021047) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF253
Synonyms:	BMZF-1; BMZF1; ZNF411
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC317840 representing NM_021047.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGGACCATTTGCAATTTAGAGATGTGGCCATAGAATTCTCTCTGGAGGAGTGGCATTGCCTGGACACT
GCACAGCGGAATTTATATAGGGATGTGATGTTAGAGAACTACAGAACTTGGTCTTCCTTGGTATTGTT
GTCTCTAAGCCAGACCTGGTTACCTGTCTGGAGCAAGGAAAAAACCTTTAACTATGGAAGACATGAG
ATGATTGCCAAACCCCGAGTTATGAGTTCTCATTGTTGCCCAAGACCTTTGGCCAGAGAACATACAAAAT
TCTTTCCAAATAGGGATGCTGAGAAGATATGAAGAATGCAGACATGACAATTTACAGTTAAAAAAGGC
TGTAAGCGTGGGTGAGCATAAGGTGCACAAAGGAGTTATAATGGACTTAACCAATGTTTGACAACT
ACCCAGAAAGAAATATTTCAATGTGATAAATATGGAAAAGTCTTCATAAGTTTTCAAATCAAACACA
TATAAGACAAGACATACTGGAATAAATCTTTCAAATGTATAATATGTGGCAAAGCTTTAAACGGTCC
TCAACCCTTACTACACATAAGAAAATTCATACTGGAGAGAAACCTTACAGATGTGAAGAATGTGGCAAA
GCTTTTAAACCAATCTGCAACCTTACTACACATAAGAGAATTCATACCGAGAGAAACCTTACAGATGT
GAAGAATGTGGCAAAGCCTTTAAGCAGTCTCAAACCTTACTACACATAAGAAAATTCATACTGGAGAG
AAACCCCTACAAATGTGAAGAATGTGGCAAAGCCTTCAACCGATCCACAGACCTTACTACACATAAGATA
GTTCACTACTGGAGAGAAACCTACAAATGTGAAGAATGTGGCAAAGCCTTTAAGCACCCTCACACGTT
ACCACACATAAGAAAATTCATACTAGAGGGAAACCTACAACCTGTGAAGAATGTGGCAAATCCTTTAAG
CACTGCTCAACCTTACTATACATAAGAGAATTCATACAGGAGAGAAACCTACAATGTGAAGAATGT
GGCAAAGCCTTTACCTATCCTCACACCTTACTACACATAAGATACTTCATACTGGAGAGAAACCTAC
AGATGTAGAGAATGTGGCAAAGCCTTTAACCATTCCACAACCTTTTTCACATGAGAAAATTCATACT
GGAGAGAAACCTACAAATGTGATGAATGTGGCAAACCTTTACCTGGCCCTCAATCCTCTCCAAACAT
AAAAGAACTCATACTGGAGAGAAACCTACAAATGTGAAGAATGTGGCAAATCCTTTACTGCATCCTCA
ACTCTAACTACACATAAGAGAATTCATACTGGAGAGAAACCTTACAAATGTGAAGAATGTGGCAAAGCT
TTAACTGGTCTCAGACCTTAATAAACATAAGAAAATTCATATTGAACGAAACCTACATAGTGAAG
AATGTGACAGATCTTTTAAATGTTCTCCACTTTTAATTAGCATAAGATA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_021047

Insert Size: 1500 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_021047.2</u>
RefSeq Size:	2351 bp
RefSeq ORF:	1500 bp
Locus ID:	56242
UniProt ID:	<u>O75346</u>
Cytogenetics:	19p13.11
Domains:	KRAB, zf-C2H2
Protein Families:	Transcription Factors
MW:	57.6 kDa
Gene Summary:	May function as a transcription factor. Seem to have a transcriptional repression activity. [UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (1) represents the longest transcript and encodes the longest isoform (1).