

Product datasheet for **SC317842**

RNF36 (TRIM69) (NM_182985) Human Untagged Clone

Product data:

| | |
|---------------------------|---|
| Product Type: | Expression Plasmids |
| Product Name: | RNF36 (TRIM69) (NM_182985) Human Untagged Clone |
| Tag: | Tag Free |
| Symbol: | TRIM69 |
| Synonyms: | HSD-34; HSD34; RNF36; Trif |
| Mammalian Cell Selection: | Neomycin |
| Vector: | pCMV6-Entry (PS100001) |
| E. coli Selection: | Kanamycin (25 ug/mL) |



[View online »](#)

Fully Sequenced ORF: >SC317842 representing NM_182985.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGGAGGTATCCACCAACCCCTCTCCAACATCGATCCAGGCGACTATGTTGAAATGAATGATTCAATC
ACCCACTACCTCTAAAGTGGTGATACAAGATATTACTATGGAGCTACACTGCCCTCTGTGCAATGAT
TGGTCCGAGACCACTGATGCTAAGCTGTGGCCACAACCTCTGTGAAGCCTGTATCCAAGACTTTTGG
AGGCTGCAAGCAAAGGAAACATTCTGTCCTGAGTGAAGATGCTATGTCAGTATAACAACGTACATTC
AACCTGTACTGGACAAGTTGGTAGAGAAGATTAAGAAGTTACCCTTACTCAAGGGCCATCCACAGTGC
CCAGAGCATGGAGAGAACCTGAACTGTTCAAGTAAACCAGATGGGAAACTGATCTGCTTTCAATGCAAG
GATGCTCGGTTGTCTGTGGGCGAGTCTAAGGAGTTCCTGCAAATCTGTATGCTGTCCATTTCTTACAG
GAGGAGCTTGCCATCCAACAGGGTCAACTGGAGACAACCTGAAGGAGCTTCAGACCCTGAGGAACATG
CAGAAGGAAGCTATTGCTGCTACAAGGAAAACAAGCTACATCTGCAGCAACATGTGTCCATGGAGTTT
CTAAAGCTGCATCAGTTCCTGCACAGCAAAGAAAAGGACATTTTAACTGAGCTCCGGGAAGAGGGGAAA
GCCTTGAATGAGGAGATGGAGTTGAATCTGAGCCAGCTTCAGGAGCAATGTCTCTTAGCCAAGGATATG
TTGGTGAGCATTAGGCAAAGACGGAACAACAGAACTCCTTCGACTTTCTCAAAGACATCAAACTCTC
TTACATAGCTTGGAGCAAGGAATGAAGGTGCTGGCAACCAGAGAGCTTATTTCCAGAAAGCTGAACCTG
GGCCAGTACAAAGGTCTATCCAGTACATGGTATGGAGGGAAATGCAGGACACTCTGCCAGGCGCTG
TCTCCACTAACTGGACCCTAAAACAGCTCACCCAAATCTGGTCTCTCCAAAAGCCAAACCAGCGTC
TGGCATGGTGACATTAAGAAGATAATGCCTGATGATCCTGAGAGGTTTACTCAAGTGTGGCTGTACTG
GGCTCAAGAGGCTTCACTCTGGAAAGTGGTACTGGGAAGTAGAAGTAGCAAAGAAGACAAAATGGACA
GTTGGAGTTGTGAGAGAATCCATCATTCCGAAGGCGAGCTGCTCTAACTCTGAGCAAGGATTCTGG
CTTTTAAGACTAAGGAACCAAACCTGATCTAAAGGCTCTGGATTTGCCTTCTTTTCAGTCTGACACTGACT
AACAACTCGACAAGGTGGGCATATACCTGGATTATGAAGGAGGACAGTTGCTCTTACATGCTAAA
ACCATGACTCACATTTACACCTCAGTAACACTTTTATGGAGAACTTTTATCCCTACTTCTGCCCTGC
CTAATGATGGTGGAGAGAATAAAGAACCATTGCACATTTACATCCACAGTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites:

SgfI-MluI

Plasmid Map:



ACCN:

NM_182985

| | |
|-------------------------------|---|
| Insert Size: | 1503 bp |
| OTI Disclaimer: | Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP). |
| OTI Annotation: | This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA. |
| Components: | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water). |
| Reconstitution Method: | <ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C. |
| RefSeq: | NM_182985.4 |
| RefSeq Size: | 1914 bp |
| RefSeq ORF: | 1503 bp |
| Locus ID: | 140691 |
| UniProt ID: | Q86WT6 |
| Cytogenetics: | 15q21.1 |
| Protein Families: | Druggable Genome |
| MW: | 57.4 kDa |
| Gene Summary: | <p>This gene encodes a member of the RING-B-box-coiled-coil (RBCC) family and encodes a protein with an N-terminal RING finger motif, a PRY domain and a C-terminal SPRY domain. The mouse ortholog of this gene is specifically expressed in germ cells at the round spermatid stages during spermatogenesis and, when overexpressed, induces apoptosis. Alternatively spliced transcript variants encoding distinct isoforms have been described. [provided by RefSeq, Jul 2008]</p> <p>Transcript Variant: This variant (a) represents the longest transcript and encodes the longest isoform (a).</p> |