

## Product datasheet for **SC317847**

### SLC36A4 (NM\_152313) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	SLC36A4 (NM_152313) Human Untagged Clone
Tag:	Tag Free
Symbol:	SLC36A4
Synonyms:	PAT4
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

**Fully Sequenced ORF:** >SC317847 representing NM\_152313.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGAAGCGGCGCGACGCCGGCGGCTGCCGGGCGCGAGGCGCGAGGAGCTAGATATGGATGTAATG
AGGCCCTTGATAAATGAGCAGAATTTTGTGGACATCAGATGAAGAACATGAGCAAGAGCTTCTGCCT
GTTCAGAAGCATTACCAACTTGATGATCAAGAGGCATTTTCATTTGTACAACTCTTATGCACCTTCTT
AAAGGAAATATTGGAAGTGGCCTTTTAGGACTTCCATTGGCAATAAAAAATGCAGGCATAGTGCTTGGA
CCAATCAGCCTTGTTTATAGGAATTATTTCTGTTCACTGTATGCACATATTGGTACGTTGCAGTCAC
TTTCTATGTCTGAGGTTTAAAAAGTCAACATTAGGTTATAGTGACACTGTGAGCTTTGCTATGGAAGTG
AGTCTTTGGAGTTGTCTTCAGAAGCAAGCAGCATGGGGCGGAGTGTGGTTGACTTTTTTCTGGTGATA
ACACAGCTGGGATTCTGTAGTGTATATTGTCTTCTTAGCTGAAAATGTGAAACAAGTTCATGAAGGA
TTCTGGAGAGTAAAGTGTATTTCAAATAGTACCAATTCATCAAACCCTTGTGAGAGAAGAAGTGT
GACCTAAGGATATATGCTTTGCTTTCTCCATTTATAATTCTTTGGTCTTCATTCGTGAACCTAAAG
AATCTATTTGACTTTTCATTCTTGCCAACGTTTCCATGGCTGTCACTGTGTGATAATTTACCAGTAT
GTTGTCAGGAACATGCCAGATCCCCACAACCTTCCAATAGTGGCTGGTTGGAAGAAATACCCACTCTTT
TTTGGTACTGCTGATTTGCTTTTGAAGGCATAGGAGTGGTCTTCCACTGGAAAACCAATGAAAGAA
TCAAAGCGTTTCCCTCAAGCGTTGAATATTGGCATGGGGATTGTTACAACCTTTGTATGTAACATTAGCT
ACTTTAGGATATATGTGTTTCCATGATGAAATCAAAGGCAGCATAACTTTAAATCTTCCCAAGATGTA
TGGTTATATCAATCAGTGAAAATCTATATTCTTTGGCATTTTTGTGACATATCAATTCAGTTCTAT
GTTCCAGCAGAGATCATTATCCCTGGGATCACATCCAAATTTCACTAAATGGAAGCAAATCTGTGAA
TTTGGGATAAGATCCTTCTTGGTTAGTATTACTGTGCCGAGCAATTCATTCTCCTTTAGACATT
GTGATTTCTTCGTTGGAGCTGTGAGCAGCAGCATTGGCCCTAATCCTGCCACCTTTGGTTGAAATT
CTTACATTTTGAAGAACATTATAATATATGGATGGTCTGAAAAATTTTCTATAGCATTCACTGGA
GTTGTTGGCTTCTATTAGGTACATATATAACTGTTGAAGAAATTTTATCCTACTCCCAAAGTTGTA
GCTGGCACTCCACAGAGTCTTTTCTAAATTTGAATTCACATGCTTAACATCTGGTTTGAATAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

**Restriction Sites:** SgfI-MluI

**Plasmid Map:**



**ACCN:** NM\_152313

<b>Insert Size:</b>	1515 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_152313.3</a>
<b>RefSeq Size:</b>	6024 bp
<b>RefSeq ORF:</b>	1515 bp
<b>Locus ID:</b>	120103
<b>UniProt ID:</b>	<a href="#">Q6YBV0</a>
<b>Cytogenetics:</b>	11q21
<b>Domains:</b>	Aa_trans
<b>Protein Families:</b>	Transmembrane
<b>MW:</b>	56.2 kDa
<b>Gene Summary:</b>	<p>SLC36A4 belongs to the SLC36 family of amino acid transporters based on sequence similarity with other family members (e.g., SLC36A1; MIM 606561). SLC36 proteins contain about 500 amino acids and have 9 to 11 transmembrane domains. Unlike other SLC36 family members, which are proton-coupled amino acid transporters, SLC36A4 is a high-affinity/low-capacity non-proton-coupled amino acid transporter (Pillai and Meredith, 2011 [PubMed 21097500]). [supplied by OMIM, Feb 2011]</p> <p>Transcript Variant: This variant (1) represents the shorter transcript and encodes the longer isoform (1). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>