

## Product datasheet for **SC317868**

### USP30 (NM\_032663) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	USP30 (NM_032663) Human Untagged Clone
Tag:	Tag Free
Symbol:	USP30
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)


**Fully Sequenced ORF:** >SC317868 representing NM\_032663.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGCTGAGCTCCCGGGCCGAGGCGGCGATGACCGCGCCGACAGGGCCATCCAGCGCTTCTCGGGACC
GGGGCGCCGTCAGATATAAAGTCATGAAGAAGTGGGGAGTTATAGGTGGAATTGCTGCTGCTCTTGCA
GCAGGAATATATGTTATTTGGGGTCCCATTACAGAAAGAAAGAAGCGTAGAAAAGGGCTTGTGCCTGGC
CTTGTTAATTTAGGGAACACCTGCTTCATGAACCTCCCTGCTACAAGGCCTGTCTGCCTGTCTGCTTTC
ATCAGGTGGCTGGAAGAGTTCACCTCCCAGTACTCCAGGGATCAGAAGGAGCCCCCTCACACCAGTAT
TTATCCTTAACACTCTTGACCTTCTGAAAGCCTTGTCTGCCAAGAAGTTACTGATGATGAGGTCTTA
GATGCAAGCTGCTTGTGGATGTCTTAAGAATGTACAGATGGCAGATCTCATCTTTGAAGAACAGGAT
GCTCACGAATTATCCATGTCATTACCTCGTCATTGGAAGATGAGCGAGACCGCCAGCCTCGGGTCACA
CATTTGTTTGATGTGCATCCCTGGAGCAGCAGTCAGAAATAACTCCAAACAAATTACCTGCCGCACA
AGAGGGTCACCTCACCTACATCCAATCACTGGAAGTCTCAACATCCTTTTCATGGAAGACTCACTAGT
AATATGGTCTGCAAACACTGTGAACACCAGAGTCTGTTGATTTGATACCTTTGATAGCCTTTCACTA
AGTATTCAGCCGCCACATGGGGTCACCCATTGACCCTGGACCCTGCCTTCACCACTTCATCTCATCA
GAATCAGTGCGGGATGTTGTGTGTGACAACGTGACAAAGATTGAAGCCAAGGGAACGTTGAACGGGGAA
AAGGTGGAACACCAGAGGACCCTTTTGTAAACAGTTAAAACAGGGAAGCTCCCTCAGTGTCTCTGC
ATCCACCTACAGCGGCTGAGCTGGTCCAGCCACGGCAGCCTCTGAAGCGGCATGAGCACGTGCAGTTC
AATGAGTTCCTGATGATGGACATTTACAAGTACCACCTCCTTGGACATAAACCTAGTCAACACAACCT
AACTGAACAAGAACCAGGGCCCTACACTGGAGCTGCAGGATGGGCCGGGAGCCCCACACCAGTCTG
AATCAGCCAGGGGCCCCAAAACACAGATTTTTATGAATGGCGCCTGCTCCCCATCTTTATTGCCAAG
CTGTGAGCGCCGATGCCCTCCCTCTCCAGTTGTTCCGACTACAGCTCCTCCACATACCTCTTCCGG
CTGATGGCAGTTGTCTCCACCATGGAGACATGCACTCTGGACACTTTGTCACTTACCGACGGTCCCCA
CCTTCTGCCAGGAACCTCTCTCAACTAGCAATCAGTGGCTGTGGGTCTCCGATGACACTGTCCGCAAG
GCCAGCTGCAGGAGTCTGTCTCCAGCGCCTACCTGCTGTTCTACGAGCGGTCTTTCCAGGATG
CAGCACAGAGCCAGGAGTGCAAGTCTGAAGAA TGA
ACGCGTACGCGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

**Restriction Sites:** SgfI-MluI

**Plasmid Map:**



**ACCN:** NM\_032663

<b>Insert Size:</b>	1554 bp
<b>OTI Disclaimer:</b>	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at <a href="mailto:custsupport@origene.com">custsupport@origene.com</a> or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a></p>
<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<a href="#">NM_032663.4</a>
<b>RefSeq Size:</b>	3771 bp
<b>RefSeq ORF:</b>	1554 bp
<b>Locus ID:</b>	84749
<b>UniProt ID:</b>	<a href="#">Q70CQ3</a>
<b>Cytogenetics:</b>	12q24.11
<b>Domains:</b>	UCH
<b>Protein Families:</b>	Protease, Transmembrane
<b>MW:</b>	58.5 kDa
<b>Gene Summary:</b>	<p>USP30, a member of the ubiquitin-specific protease family (see USP1, MIM 603478), is a novel mitochondrial deubiquitinating (DUB) enzyme (Nakamura and Hirose, 2008 [PubMed 18287522]).[supplied by OMIM, Dec 2008]</p> <p>Transcript Variant: This variant (1) encodes the longer isoform (1).</p>