

Product datasheet for **SC317886**

NPEPL1 (NM_024663) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	NPEPL1 (NM_024663) Human Untagged Clone
Tag:	Tag Free
Symbol:	NPEPL1
Synonyms:	bA261P9.2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)


Fully Sequenced ORF: >SC317886 representing NM_024663.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGGCGAACGTGGGGCTGCAGTTCAGGGCAGCGCGGGGACTCGGACCCACAGAGCCGGCCCTGCTG
CTGCTCGGGCAGCTGCACCACCTGCACCCGCTGCCCTGGAGCCAGTCCGCGGGAAGCTGCAGCCCGG
GTACCAGGAGACTCTGGCAGGCTGCCCTGAGCAGCTCAACCCAACCCACGGACAGCTGTCCCTC
TACCTGAACTACGCCACCGTGGCTGCCCTGCCCTGCAGGGTGAGCCGGCACAACAGCCCTCGGCCGCC
CACTTCATCACGCGGCTGGTGGGACCTGCCTGCCGCCGGAGCGCATCGCTGCATTGTGATGGTCTGC
GAGCAGCCGGAGGTCTTTGCTTCCGCTGTGCCCTGGCCGGGCCTTCCCGCTGTTACCCACCGCTCA
GGTGCCTCTCGGCGCTTGAGAAGAAGACGGTACCCTGGAGTTTTTCTGGTGGACAAGACAACGGG
CCGGTGGAGGTGCCACATTGCAGTGCTTAGCGAATGCCACAGACGGCTGCGGTAGCAGCCCGCATC
GTGGACACACCCTGCAATGAGATGAACACCGACACCTTCTCGAGGAGATTAACAAGTTGGAAAGGAG
CTGGGGATCATCCAACCATCATCCGGGATGAGGAAGTGAAGACGAGAGGATTTGGAGGAATCTATGGG
GTTGGCAAAGCCGCCCTGCATCCCCAGCCCTGGCCGTCTCAGCCACACCCAGATGGAGCCACGCAG
ACCATCGCTGGGTGGGCAAAGGCATCGTCTATGACACTGGAGGCCTCAGCATCAAAGGGGAGACTACC
ATGCCGGGGATGAAGCGAGACTGCGGGGGTGTGCGGCCGTCTGGGGGCCTTCAGAGCCGAATCAAG
CAGGGTTTCAAAGACAACCTCCACGCTGTGTTCTGCTTGGCTGAGAACTCGTGGGGCCCAATGCGACA
AGGCCAGATGACATCCACTGCTGACTCAGGGAAGACGGTGGAAATCAACAACCGGATGCCGAGGGC
AGGCTGGTGTGGCAGATGGCGTGTCTATGCTTGCAAGACCTGGGGCCGACATCCTGGACATG
GCCACCCTGACCGGGCTCAGGGCATTGCCACAGGGAAGTACCACGCCGGTGTCAACACAGCGCT
GAGTGGGAGCCGCTGTGTGAAGCGGGCAGGAAGTGTGGGACCTGGTGCACCCGCTGGTCTACTGC
CCCGAGCTGCACTTCAGCGAGTTCACCTAGCTGTGGCGGACATGAAGAAGTCACTGGCGGACCGAGAC
AACAGCCCGAGCTCCTGTGCTGGCCTTTCATCGCCTCACACATCGGCTTCGACTGGCCCGAGTCTGG
GTCCACCTGGACATTGCTGCACCGGTGCATGCTGGTGAGCGAGCCACAGGCTTCGGTGTGGCCCTCTG
CTGGCGCTCTTCGGCCGTGCCTCTGAGGACCTCTGCTGAACCTGGTGTCCCACTGGGCTGTGAGGTG
GATGTCGAGGAGGGGGACCTGGGAGGGACTCAAAGAGACGCAGGCTTGTGTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_024663

Insert Size:	1572 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_024663.3
RefSeq Size:	2217 bp
RefSeq ORF:	1572 bp
Locus ID:	79716
UniProt ID:	Q8NDH3
Cytogenetics:	20q13.32
Domains:	Peptidase_M17
Protein Families:	Druggable Genome, Protease
MW:	55.9 kDa
Gene Summary:	Probably catalyzes the removal of unsubstituted N-terminal amino acids from various peptides.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (1) represents the longest transcript and encodes the longest isoform (1).