

Product datasheet for **SC317894**

ZNF35 (NM_003420) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF35 (NM_003420) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF35
Synonyms:	HF.10; HF10; Zfp105
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)


Fully Sequenced ORF: >SC317894 representing NM_003420.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGACTGCAGAATTGAGAGAAGCCATGGCCCTAGCCCCATGGGGCCAGTGAAGGTGAAAAAGGAGGAG
GAAGAAGAAGAAAACCTCCAGGTCCAGGCATCCAGCCAACAAGTGCACCTCCGAGAACATCAAAGTCTGG
GCCCCAGTGCAGGGTCTTCAGACAGGCCTTGATGGATCAGAAGAGGAAGAAAAGGGTCAAGACATATCC
TGGGATATGGCGGTAGTCTGAAAGCAACTCAGGAGGCACCTGCTGCTTCAACCCCTTGGCAGCTACTCA
TTACCAGGGACTCTGGCCAAGAGTGAGATACTGGAGACTCATGGGACCATGAACTTTCTAGGTGCTGAA
ACCAAGAACCTACAGTTACTGGTTCCAAAACTGAGATATGTGAGGAAGCTGAAAAACCCCTCATCATA
TCAGAAAGAAATCCAGAAAGCTGATCCTCAAGGACCTGAGTTAGGAGAAGCTTGTGAAAAGGGAAACATG
TTAAAGAGGCAGAGAATAAAGAGAGAAAAGAAAGATTTTCAGACAAGTGATAGTGAATGACTGTCACCTTA
CCTGAAAGCTTCAAAGAAGAGGAAAACCCAGAAATGTAAGAAATCTGGAGGAAAATATAGCCTTAATTCT
GGCGCTGTAAAAATCCAAAACCCAGCTTGGACAAAAGCCTTTTACGTGTAGCGTGTGTGGGAAAGGA
TTTAGTCAGAGTGCAAACCTCGTTGTGCATCAGCGAATCCACACTGGAGAGAAAACCCCTTTGAATGTCAT
GAGTGTGGGAAGGCCCTTCATTCAGAGTGCAAACCTCGTTGTGCATCAGAGAATCCACACTGGACAGAAA
CCTTATGTTTGCTCAAAATGTGGGAAAGCCCTTCACTCAGAGTTCAAATCTGACTGTACATCAAAAAATC
CACTCCTTAGAAAAAACTTTTAAAGTGAATGAATGTGAGAAAAGCCTTTAGTTACAGCTCACAACTTGCT
CGGCACCAGAAAGTCCACATTACGGAAAAATGCTATGAATGTAATGAATGTGGGAAAACATTTACTAGG
AGCTCAAACCTCATTGTCCACCAGAGGATCCACACTGGGAGAAGCCCTTTGCCTGTAACGACTGTGGC
AAAGCCTTACCAGAGTGCAAATCTTATTGTACATCAGCGAAGCCATACTGGTGAAGCCATATGAG
TGTAAGAGGTGTGGGAAAGCCCTTAGTTGTTTTTACACCTTATTGTGCACCAGAGAATTCACACTGCA
GAGAAACCTTACGACTGCAGCGAATGTGGGAAAGCCCTTCACTCAGCTCTCTTGCCTTATTGTCCACCAG
AGAATTCACAGTGAGATCTTCTTACGTGTGAATGAATGTGGGAAAGCCCTTCACATGTAGCTCATAAC
CTACTTATTCATCAGAGAATTCATAATGGAGAAAAACCTTACACATGTAATGAGTGTGGGAAAGCCCTTC
AGACAGAGGTCGAGCCTCACCGTGCACCAGAGAACCCACACTGGGAGAAGCCCTATGAATGTGAGAAAG
TGTGGTGCAGCTTTCATTTCCAACCTCACACCTCATGCGACACCATAGAACCCATCTTGTGAATAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_003420

Insert Size:	1584 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_003420.3
RefSeq Size:	2664 bp
RefSeq ORF:	1584 bp
Locus ID:	7584
UniProt ID:	P13682
Cytogenetics:	3p21.31
Domains:	zf-C2H2
Protein Families:	Transcription Factors
MW:	59.1 kDa
Gene Summary:	May be involved in transcriptional regulation. Involved in cell differentiation and/or proliferation.[UniProtKB/Swiss-Prot Function]