

## Product datasheet for **SC317897**

### **C14orf45 (BBOF1) (NM\_025057) Human Untagged Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	C14orf45 (BBOF1) (NM_025057) Human Untagged Clone
Tag:	Tag Free
Symbol:	BBOF1
Synonyms:	C14orf45; CCDC176; FBB10
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**Fully Sequenced ORF:** >SC317897 representing NM\_025057.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGCCGTGAAGGAAAGGACAAAAAGAAAGGCAAGAGCAAAGGCAAAGACCGAAGAAGTTAATAAAAA
ACAGATGAATCTGTGGTGGACAGAGCCAAGGCCAATGCCTCCCTTTGGGAGGCCAGTTGGAAGTCACA
GAACTCTCTAGGATTAAGTATCGTGATACTTCACGGATACTGGCAAAAAGTAATGAGGACTTAAAGAAA
AAGCAATGTAATAATGGAGAAAGACATAATGTCAGTATTAAGTTACCTGAAGAAGCAGGATCAGGAGAAA
GATAATATGATTGAAAACTGAAACAGCAATTAATGAAACAAAGGAAAAAGCCCAAGAGGAGAAGGAT
AAATTGGAACAAAAGTATACCAGGCAAATTAATGAACTAGAGGGACAGTTCCATCAAAAAGCCAAAGAA
ATTGGCATGATTCACACAGAGCTGAAAGCAGTAAGACAATCCAGAAGAGAAAAATCCAAGTGGAGAGA
GAGTTAGATGATCTGAAAGAGAATTTAAGAAACACAGAGCGGATACATCAGGAGACTCTTAGAAGACTG
GAAAGCAGATTTTTTGAAGAAAAGCACCGACTAGAACAAGAGGCTGAAAAGAAGATAATAATGCTAGCA
GAGAGAGCCACCATGAGGCTATTGTGCAATTGAACGATGCTGGAAGAAATGTTTTTAAAGAGAATGAT
TATCTTCAGAAAGCTCTGGCATATCACCTGAAGGAAACTGACGCTCTACAAAAAACTCCAGAAGTTG
CAAGAGAGTCATACTTTACTTTTACATCAAAAGGAGATCAATGATCTGTTGGTTAAGGAAAAAGATTATG
CAACTTGTCCAGCAGAGATCACAAATCCAAACCTTCAGAAGAAGGTAGTAACTTGGAGACTGCTCTG
AGTTACATGACCAAAGAGTTTGAGAGTGAAGTTTTAAACTGCAGCAACACGCAATGATAGAGAACCAA
GCAGGTGAGGTAGAAATTGACAAGCTGCAGCACCTTCTCAGATGAAGGACAGGGAAATGAATCGTGTG
AAGAAGCTGGCCAAGAACATACTGGATGAGAGAACAGAAGTGGAAAGATCTTTTTAGATGCTCTGCAC
CAAGTGAAGCAACAGATCCTAATTAGCAGGAAGCATTATAAGCAGATAGCACAAGCTGCTTTCAATTTA
AAAATGAGAGCAGCATGTACAGGAAGAACAATAATCCCAAAATCAGAACATTTGATGGCAGAGAGCAC
AGCACC AATAGTGAATCAGGATCTTCTGGAGGCCGAAAAATGGACACATATTGAAGGAAATGTGGAT
ATTGGAGATTTGACCTGGGAGCAGAAGGAAAAAGTATTGCGATTGCTCTTTGCAAAAATGAATGGCTGT
CCTTCTAGGAAATACAACCAGAGTTCTAGGCCTCCAGTCCAGACTATGTTGTTTCTGACAGTGGGGAA
ACAAAGGAATTTGGGGATGAAAGTAAGCTTCAAGATAAAATCTTCATCACCCAGCAAATGCAATATCA
GACTCTTCTGGTGAAGTGGTGTACCCACTATTCCAAAAGAACCTCAGGAGTCTGACACAGGAACTTC
TGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
  
```

**Restriction Sites:**

Sgfl-MluI

**Plasmid Map:**



<b>ACCN:</b>	NM_025057
<b>Insert Size:</b>	1590 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_025057.2</a></u>
<b>RefSeq Size:</b>	3123 bp
<b>RefSeq ORF:</b>	1590 bp
<b>Locus ID:</b>	80127
<b>UniProt ID:</b>	<u><a href="#">Q8ND07</a></u>
<b>Cytogenetics:</b>	14q24.3
<b>MW:</b>	62 kDa
<b>Gene Summary:</b>	Basal body protein required in multiciliate cells to align and maintain cilia orientation in response to flow. May act by mediating a maturation step that stabilizes and aligns cilia orientation. Not required to respond to planar cell polarity (PCP) or flow-based orientation cues (By similarity).[UniProtKB/Swiss-Prot Function]