

Product datasheet for **SC317903**

PARP3 (NM_005485) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PARP3 (NM_005485) Human Untagged Clone
Tag:	Tag Free
Symbol:	PARP3
Synonyms:	ADPRT3; ADPRTL2; ADPRTL3; ARTD3; IRT1; PADPRT-3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF:

>SC317903 representing NM_005485.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCTCCAAAGCCGAAGCCCTGGGTACAGACTGAGGGCCCTGAGAAGAAGAAGGGCCGGCAGGCAGGA
AGGGAGGAGACCCTTCCGCTCCACCGCTGAGGCCCTCAAGCCATACCCGAGAGAAGCGCATAATC
CGCGTGGATCCAACATGTCCACTCAGCAGCAACCCGGGACCCAGGTGTATGAGGACTACAACCTGCACC
CTGAACCAGACCAACATCGAGAACAACAACAAGTTCTACATCATCCAGCTGTCTCAAGACAGCAAC
CGCTTCTTACCTGCTGGAACCGCTGGGGCCGTGTGGGAGAGGTCGGCCAGTCAAAGATCAACCACTTC
ACAAGGCTAGAAGATGCAAAGAAGGACTTTGAGAAGAAATTTGGGAAAAGACCAAGAACAACCTGGGCA
GAGCGGGACCCTTTGTGTCTACCCGGGCAAGTACACACTTATCGAAGTACAGGCAGAGGATGAGGCC
CAGGAAGCTGTGGTGAAGGTGGACAGAGGCCAGTGAGGACTGTGACTAAGCGGGTGCAGCCCTGTCC
CTGGAACCCAGCCACGAGAAGCTCATCTAACATCTTCAGCAAGGAGATGTTCAAGAACCACCATGGCC
TCATGGACCTGGATGTGAAGAAGATGCCCTGGGAAAGCTGAGCAAGCAACAGATTGCACGGGGTTTC
GAGGCCTTGGAGGCGCTGGAGGAGGCCCTGAAAGGCCCCACGGATGGTGGCCAAAGCCTGGAGGAGCTG
TCCTCACACTTTTACACCGTCATCCCGCACAACTTCGGCCACAGCCAGCCCCGCCCATCAATTC CCT
GAGCTTCTGCAGGCCAAGAAGGACATGCTGCTGGTGTGGCGGACATCGAGCTGGCCAGGCCCTGCAG
GCAGTCTCTGAGCAGGAGAAGACGGTGGAGGAGGTGCCACACCCCTGGACCGAGACTACCAGCTTCTC
AAGTGCCAGCTGCAGCTGCTAGACTCTGGAGCACCTGAGTACAAGGTGATACAGACCTACTTAGAACAG
ACTGGCAGCAACCACAGGTGCCCTACACTTCAACACATCTGAAAGTAAACCAAGAAGGGGAGGAAGAC
AGATCCAGGCCCACTCCAACTGGGTAATCGGAAGCTGCTGTGGCATGGCACCACATGGCCGTGGTG
GCCGCCATCCTCACTAGTGGGCTCCGCATCATGCCACATCTGGTGGGCGTGTGGCAAGGCATCTAC
TTTGCCTCAGAGAACAGCAAGTCAGCTGGATATGTTATTGGCATGAAGTGTGGGGCCACCATGTCGGC
TACATGTTCTGGGTGAGGTGGCCCTGGCAGAGACCATATCAACACGGACAACCCAGCTTGAAG
AGCCACCTCCTGGCTTCGACAGTGTCAATTGCCCGAGGCCACACCGAGCCTGATCCGACCCAGGACACT
GAGTTGGAGCTGGATGGCCAGCAAGTGGTGGTGGCCCGAGGCCAGCCTGTGCCCTGCCAGAGTTGAGC
AGCTCCACATTCTCCAGAGCGAGTACCTCATCTACCAGGAGAGCCAGTGTGCCTGCGCTACCTGCTG
GAGGTCCACCTTGA
ACGGGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites:

Sgfl-MluI

Plasmid Map:



ACCN:	NM_005485
Insert Size:	1602 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_005485.5</u>
RefSeq Size:	2416 bp
RefSeq ORF:	1602 bp
Locus ID:	10039
UniProt ID:	<u>Q9Y6F1</u>
Cytogenetics:	3p21.2
Domains:	PARP, PARP_reg
Protein Families:	Druggable Genome
Protein Pathways:	Base excision repair
MW:	60.1 kDa
Gene Summary:	<p>The protein encoded by this gene belongs to the PARP family. These enzymes modify nuclear proteins by poly-ADP-ribosylation, which is required for DNA repair, regulation of apoptosis, and maintenance of genomic stability. This gene encodes the poly(ADP-ribosyl)transferase 3, which is preferentially localized to the daughter centriole throughout the cell cycle. Alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, Jul 2008]</p> <p>Transcript Variant: This variant (2) differs in the 5' UTR and uses a downstream translational start codon, compared to variant 1. The resulting isoform (b) is shorter at the N-terminus, compared to isoform a.</p>