

Product datasheet for **SC317906**

CCDC9B (NM_207380) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CCDC9B (NM_207380) Human Untagged Clone
Tag:	Tag Free
Symbol:	CCDC9B
Synonyms:	C15orf52
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC317906 representing NM_207380.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGATCTCGTGTGCTGAGCAGCGAAGCCGGCAGGGAGAGGCCGGCAGAGGCCGGCTCCGGTGGCTCCA
GCTTTCTCCCACTCTGGCTCCCAGGGGCTGCTTGAATTCTCTCGGTGCCCGCGGTTGCCATGCAC
TCGGCTGGAACCTCCAGAGCCGAGTCCCCATGAGCAGGCAGGAGAAGGACGCAGAGCTGGATCGGAGG
ATAGTTGCCCTGCGCAAGAAGAACCAGGCCTTGCTCCGCAAGTACCAGGAGATCCAGGAGGACCGTCGG
CAGGCAGAGCAGGGGGGATGGCTGTGACCACACCAGCACTCCTCCAGCCTGATGGCCTACCGTTACC
ATCAGCCAGGTTCCCGGTGAAAAGCGGGTGGTTAGCAGGAAGTGGCAAGGGGTACCTGTGGACCCAGA
GTGACCAACGAGATGCTTGAGGATGAGGATGCTGAGGACCACGGGGTACTTTCTGCTTAGGGGAGCTG
GTGGAGCTGGCTGTGACCATGGAGAACAAGCAGAGGGCAACCGATTGTAAGTAAAAGCCTACCAGA
GCAAGGAACCAAGGCATAGAGGGGTACCTGGAGGGCTGTGACCCGAAGCCCCCACGAGGTGGCC
ATCAGCTCAGATTCTGCACGGAAGGGTTCTTGGAGCCCTGGAGCCGGCCGGTGGGGGAGCCCCGGAG
GCGGGCTGGGACTATGCCAGTGAAGCAGGAGCGGGAGCAGATCGACCTAGCCCGCCTCGCCCGGCAC
AGAGACGCACAGGGTGACTGGCGCCGCCGTGGGACCTGGACAAGGCCAAGTCCACGCTACAGGACTGC
AGCCAGCTGAGGGGAGAAGGCCCGGCCAGGGCAGGCAGCAGAAGGGGTCCAGGAGCCACCAGAACTA
CAGCCCCCACCATTGCTCCCTGATGAAAAGGTGCGGGCGGGCAAGCCAGCAGACCCTCGGTGGCACC
GCCACAGGCAGCAAGGCCGGGGCAAGGAGAGGCTGACTGGCAGGGCCGAAGGTGGGATATGAAGGAA
GACAAGGAGGAGCTGGAAGGTGAGGAGGAAGCCAAAGCACCAGAGAGACTCCAGTGAGGAGGAGCAA
GCCCAGAGCAGAGCGGGATGGAGCAGGGCCGACTGGGGAGCGCCCTGCAGCCAGCCCAGCCCTGGCA
TCCCCAGAGGGGCGGAAGGGGGAGTCAGTGGCTTCCACAGCCAGCTCAGTCCCCTGCTCTCCACAGGAG
CCTGACTTGGCTCCTCTTGACCTCTCCCTAGGAGGGGCTGGCATCCCTGGGCCCAGGGAGAGCGGGTGT
GTGCTCGGTCTGAGGCCTGGGGCCAGGAGAGCCCTGTGTCTTGCCAGAGGGCTCTAAGCAGCAGCCC
CTGGGGTGGAGCAATCACCAGGCTGAGCTGGAAGTACAGACTTGCCCTGAGCCACAGAGAGGAGCAGGG
CTCCCAGAGCCCGGAGAAGACAGGTCTGGCAAGTCTGGGGCCAGCAGGGCCTGGCCCCGAGAAGCCGG
CCCACGAGAGGAGGCAGCCAAAGGTGAGAGGCACAGCAGGTGTGAGGCGCAGGACAGGGCGCCCTGGC
CCGGCAGGAAGATGCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_207380

Insert Size: 1605 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: NM_207380.2

RefSeq Size: 5344 bp

RefSeq ORF: 1605 bp

Locus ID: 388115

UniProt ID: Q6ZUT6

Cytogenetics: 15q15.1

MW: 57.3 kDa