

Product datasheet for **SC317908**

DYNC2I2 (NM_052844) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	DYNC2I2 (NM_052844) Human Untagged Clone
Tag:	Tag Free
Symbol:	DYNC2I2
Synonyms:	bA216B9.3; CFAP133; DIC5; FAP133; SRTD11; WDR34
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC317908 representing NM_052844.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCAACCCGCGCCAGCCGGGGCCACTCAGCCAGGCGGGAAGCGCTGGTGTTCGGCGCTGGCGACA
GTCGGGGTTGCGAGCGGCCCGGGCCGGCCAGGGCCGCTGCAGGACGAGACCCTGGGTGTGGCG
TCCGTGCCCTCGCAGTGGAGGGCCGTCCAGGGCATCCGCTGGGAGACGAAAAGTTGCCAGACGGCCAGC
ATTGCCACTGCCAGTGCATCCGCCAGGCCAGGAATCATGTGGACGCCAGGTGCAGACGGAGGCCCC
GTGCCTGTGAGCGTGCAGCCCCGTCCCAGTATGACATACCCAGGCTCGCAGCCTTTCTTCGGAGAGTG
GAGGCCATGGTATCCGAGAGCTGAACAAGAATTGGCAGAGCCACGCTTTGATGGCTTCGAGGTGAAC
TGGACCGAGCAGCAGCAGATGGTGTCTTGTCTGTATACCCTGGGTACCCGCCAGCCCAAGCGCAGGGT
CTGCATGTGACCAGCATCTCTGGAACCTCACTGGCTCTGTGGTGGCCTGTGCCTACGGCCGGCTGGAC
CATGGGACTGGAGCAGCTTAAGTCCTTCGTGTGTGCCTGGAACCTGGACCGGCGAGACCTGCGTCCC
CAGCAGCCGTGCGCCGTGGTGGAGTCCCAGCGCTGTCTGTGTCTGGCCTTCCACCCACGCAGCCC
TCCCACGTCGCAGGAGGGCTGTACAGTGGTGAGGTGTTGGTGTGGGACCTGAGCCGTCTTGAGGACCCG
CTGCTGTGGCGCACAGGCCTGACGGATGACACCCACAGACCCTGTGTCCCAGGTGGTGTGGCTGCC
GAGCCTGGGCACAGCCACCGCTTCCAGGTGCTGAGTGTGGCCACCGAGCGGAAGGTGCTACTCTGGCAG
GGCATCGGGTAGGCCAGCTGCAGCTCACAGAGGGCTTCGCCCTGGTCATGCAGCAGCTGCCACGGAGC
ACCAAGCTCAAGAAGCATCCCGCGGGGAGACCGAGGTGGGCGCCACGGCAGTGGCCTTCTCCAGCTTT
GACCCTAGGCTGTTTATTCTGGGCACGGAAGCGGCTTCCCGTCAAGTGTTCCTGGCAGCTGGAGAG
GCAGCCCTCAGCGGATGCCAGCTCCGTGCCCTGCGGGCCCGCAGCAGTTTACCTTCTCCCCCAC
GGCGTCCCCTACTCTGTGAGCTGTTCCCTTCCACAGGAATCTTCTCTGAGCGCTGGGACTGAC
GGGCATGTCCACCTGACTCCATGCTGCAGGCCCTCCCTTGACTTCGCTGCAGCTCTCCCTCAAGTAT
CTGTTTGTGTGCGCTGGTCCCCAGTGGCGCCCTTGGTTTTTGCAGCTGCCTCTGGGAAAGTGACGTG
CAGCTGTTTGTCTCCAGAAAAGCTCCCAGAAAACCCACAGTTTTGATCAAGCAAACCCAGGATGAAAGC
CCTGTCTACTGTCTGGAGTTCAACAGCCAGCAGACTCAGCTTTGGTGCGGGCGATGCCAGGGCACA
GTGAAGGTGTGGCAGCTGAGCACAGAGTTCACGGAACAAGGGCCCGGGAAGCTGAGGACCTGGACTGC
CTGGCAGCAGAGGTGGCGGCC TGA
AGCGGACCGACGCTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGAT
ATCCTGGATTACAAGGATGACGACGATAAGGTTTAA
  
```

Restriction Sites:

Sgfl-RsrII

Plasmid Map:



ACCN:	NM_052844
Insert Size:	1611 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_052844.3</u>
RefSeq Size:	1818 bp
RefSeq ORF:	1611 bp
Locus ID:	89891
UniProt ID:	<u>Q96EX3</u>
Cytogenetics:	9q34.11
Domains:	WD40
MW:	57.8 kDa
Gene Summary:	<p>This gene encodes a member of the WD repeat protein family. WD repeats are minimally conserved regions of approximately 40 amino acids typically bracketed by gly-his and trp-asp (GH-WD), which may facilitate formation of heterotrimeric or multiprotein complexes. Members of this family are involved in a variety of cellular processes, including cell cycle progression, signal transduction, apoptosis, and gene regulation. Defects in this gene are a cause of short-rib thoracic dysplasia 11 with or without polydactyly. [provided by RefSeq, Mar 2014]</p>