

Product datasheet for **SC317915**

ANKRD13C (NM_030816) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ANKRD13C (NM_030816) Human Untagged Clone
Tag:	Tag Free
Symbol:	ANKRD13C
Synonyms:	dj677H15.3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC317915 representing NM_030816.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGACCGGGGAGAAGATCCGCTCACTGCGGAGGGACCACAAGCCCAGCAAAGAAGAGGGGACCTGCTG
GAGCCCGGGGATGAGGAAGCGCGGCTGCCCTCGGCGGTACCTTTACCAGAAGCAGGATTGGCAAGGGC
GGCAAAGCTTGTCTATAAGATCTTCAGTAACCATCACACCGGCTACAGCTGAAGGCAGCTCCGGCCTCC
TCCAATCCCCCGGCGCCCGGCTCTGCCGCTGCACAATTCCTCCGTGACTGCCAACTCCAGTCCCGG
GCCCTTCTGGCCGGCACCAACCCGTTGCTGTCTGCGGATGGAGGCAGTTGCCCGCACACTACCCG
GTGCACGAGTGCGTCTTCAAGGGGATGTGAGGAGACTCTCTCTCATCCGCACGCACAATATCGGG
CAGAAAGATAATCACGAAATCTCCTTTACACCTTGCTGTGATGTTAGGAAATAAAGAATGTGCCCAT
TTACTTTTGGCTCACAATGCTCCAGTCAAGGTGAAAAATGCTCAGGGATGGAGCCCTCTGGCGGAAGCC
ATCAGCTATGGAGATAGGCAGATGATTACAGCTCTTTGAGGAAGCTTAAGCAGCAATCCAGGGAAAGT
GTTGAAGAAAAACGACCTCGATTATTAAGGCCCTGAAAGAGCTAGGTGACTTTTATCTAGAACTTAC
TGGGATTTTCAAAGCTGGGTGCCTTTACTTTCCCGAATTCGCTTCCGATGCATGTAAATATACAAA
CAAGGTATCAATATCAGGCTTGACACAACCTCTCATAGACTTTACTGACATGAAGTGCCAACGAGGGGAT
CTAAGCTTCATTTTCAATGGGGATGCGGCGCCCTCTGAATCTTTTGTAGTATTAGACAATGAACAAAA
GTTTATCAGCGAATACATCATGAGGAATCAGAGATGGAACAGAGAAGAGGTGGATATTTAATGAGC
AGTGATATTTACTCTGCAACTTTATCAACAAAATCAATTTCTTTACGCGTGCCAGACAGGATGGCTT
TTTCGGGAAGATAAAACAGAAAGAGTAGGAACTTTTGGCAGACTTTTACCTGGTGAATGGACTTGT
TTAGAATCAAGGAAAAGAAGAGAACATCTCAGTGAAGAGGATATCTTCGAAATAAGGCCATCATGGAG
AGTTTGAGTAAAGGTGGAACATAATGGAACAGAATTTTGAGCCGATTCGAAGACAGTCTCTTACACCT
CCTCCTCAGAACACTATTACATGGGAAGAATATATCTGCTGAAAATGGAAGCTCCTCATCTGGGT
AGAGAATTGGTGTGCAAAGAGAGTAAGAAAACGTTTAAAGCTACGATAGCCATGAGCCAGGAATTTCCC
TTAGGGATAGAGTTATTATTGAATGTTTGAAGTAGTAGCTCCCTTCAAGCACTTTAACAAGCTTAGA
GAATTTGTTTCAAGCTTCTCCAGGCTTTCTGTAAATTAGATATACCTGTGTTTCCACAATC
ACAGCCACTGTGACTTTTCAAGGAGTTTCGATACGATGAATTTGATGGCTCCATCTTTACTATACCTGAT
GACTACAAGGAAGACCAAGCCGTTTCTGATCTTAA
AGCGGACCGACGCTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGAT
ATCCTGGATTACAAGGATGACGACGATAAGGTTTAA
  
```

Restriction Sites: SgfI-RsrII

ACCN: NM_030816

Insert Size: 1626 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_030816.4</u>
RefSeq Size:	5682 bp
RefSeq ORF:	1626 bp
Locus ID:	81573
UniProt ID:	<u>Q8N6S4</u>
Cytogenetics:	1p31.1
MW:	60.8 kDa
Gene Summary:	Acts as a molecular chaperone for G protein-coupled receptors, regulating their biogenesis and exit from the ER.[UniProtKB/Swiss-Prot Function]