

## Product datasheet for **SC317920**

### CSGALNACT2 (NM\_018590) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	CSGALNACT2 (NM_018590) Human Untagged Clone
Tag:	Tag Free
Symbol:	CSGALNACT2
Synonyms:	beta4GalNAcT; ChGn-2; CHGN2; GALNACT-2; GALNACT2; PRO0082
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**Fully Sequenced ORF:** >SC317920 representing NM\_018590.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGCCTAGAAGAGGACTGATTCTTACACCCGGACCCACTGGTTGCTGTTGGGCCTTGCTTTGCTCTGC
AGTTTGGTATTATTTATGTACCTCTGGAATGTGCCCCAGACTGATGAAATGCATCTCTTCCTGGT
GTTGTTGGGAAAATTATGGTAAAGAGTATTATCAAGCCCTCTACAGGAACAAGAAGAATTATCAG
ACCAGGGCAACCAAGTCTGAAACGCCAAATTGCCAACTAAAAACAAGAATTACAAGAAATGAGTGAGAAG
ATGCGGTCACTGCAAGAAAGAAGGAATGTAGGGGCTAATGGCATAGGCTATCAGAGCAACAAGAGCAA
GCACCTAGTGATCTTTAGAGTTTCTTCATTCCCAAATTGACAAAGCTGAAGTTAGCATAGGGGCCAAA
CTACCCAGTGAGTATGGGTCAATCCCTTTGAAAGTTTACCTTAATGAAAGATTTCAATTGAAATG
GGTCTCACTCGCCATCTGAAGAAAAGCCAGTTAGAAAAGACAAACGAGATGAATTGGTGGAAAGTTATT
GAAGCGGGCTTGGAGTCAATAATCCTGATGAAGATGATGAACAAGAAGATGAGGAGGGTCCCCTT
GGAGAGAACTGATATTTAATGAAAATGACTTCGTAGAAGTTATTATCGCACTGAGAGAGATAAGGGC
ACACAGTATGAACTCTTTTTAAGAAAGCAGACCTTACGGAATATAGACATGTGACCCTCTCCGCCCT
TTTGGACCTCTCATGAAAGTGAAGAGTGAGATGATTGACATCACTAGATCAATTATTAATATCATTGTG
CCACTTGCTGAAAGAAGTGAAGCATTGTACAATTTATGCAGAAGTTCAGGGATGTTTGTATTATCAA
GACAAGAAGATTCATCTCACAGTGGTGTATTTGGTAAAGAAGGACTGTCTAAAGTCAAGTCTATCCTA
GAATCTGTCACCAGTGAGTCAATTTTACAATTACACCTTGGTCTCATTGAATGAAGAATTTAATCGT
GGACGAGGACTAAATGTGGTGCCCGAGCTTGGGACAAGGAGAGGTCTTGATGTTTTCTGTGATGTT
GATATCTATTTCTCAGCCGAATTCCTTAACAGCTGCCGTTAAATGCTGAGCCAGGTAAGAAGGTGTTT
TACCCTGTGGTGTTCAGTCTTTACAATCCTGCCATTGTTTATGCCAACAGGAAGTCCACCACCTGTG
GAGCAGCAGCTGGTTACAAAAAGGATTCTGGCTTTTGGCGAGATTTTGGCTTTGGAATGACTTGTGAG
TATCGTTCAGATTTCTGACCATTGGTGGATTTGACATGGAAGTGAAGGTTGGGGTGGAGAAGATGTT
CATCTTTATCGAAAACTTACATGGTGACCTCATTGTGATTCCGACTCCGGTTCCTGGTCTTTCCAC
CTCTGGCATGAAAAGCGCTGTGCTGATGAGCTGACCCCGAGCAGTACCGCATGTGCATCCAGTCTAAA
GCCATGAATGAGGCTCTCACTCCACCTGGGAATGCTGGTCTTCAGGGAGGAAATAGAGACGCATCTT
CATAAACAGGCATACAGGACAAACAGTGAAGCTGTTGGTTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
  
```

**Restriction Sites:**

Sgfl-MluI

**Plasmid Map:**



<b>ACCN:</b>	NM_018590
<b>Insert Size:</b>	1629 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<u><a href="#">NM_018590.4</a></u>
<b>RefSeq Size:</b>	3782 bp
<b>RefSeq ORF:</b>	1629 bp
<b>Locus ID:</b>	55454
<b>UniProt ID:</b>	<u><a href="#">Q8N6G5</a></u>
<b>Cytogenetics:</b>	10q11.21
<b>Protein Families:</b>	Transmembrane
<b>Protein Pathways:</b>	Chondroitin sulfate biosynthesis, Metabolic pathways
<b>MW:</b>	62.6 kDa
<b>Gene Summary:</b>	This gene encodes a member of the chondroitin N-acetylgalactosaminyltransferase family. The encoded protein is involved in elongation during chondroitin sulfate synthesis. Alternative splicing of this gene results in multiple transcript variants. Two related pseudogenes have been identified on chromosome X. [provided by RefSeq, Feb 2016] Transcript Variant: This variant (1) encodes the longest isoform (1).