

Product datasheet for **SC317922**

ZCCHC7 (NM_032226) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZCCHC7 (NM_032226) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZCCHC7
Synonyms:	AIR1; HSPC086
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC317922 representing NM_032226.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGATGTTTGGTGGCTATGAGACTATAGAAGCATACGAAGATGATCTTTATCGAGATGAGTCATCTAGT
GAAGTGAAGTGGTATGAGTGAAGTGAATTTCAACTCTATAGCCAATTCATTATGCCAAGATCTTGAT
GATGTCATCAGGGAGGAAGAGCATGAAGAAAAGAAGCTCTGGGAATTCGGAATCTTCGAGTAGTAAACCA
AATCAGAAGAAGCTAATCGTCCTTTTCAGATAGTGAGGTCATCCAGCTGTGAGATGGGTCAGAGGTCATC
ACTTTGTCTGATGAAGACAGTATTTATAGATGTAAGGAAAGAATGTTAGAGTTCAAGCACAAGAAAAT
GCCCATGGTCTTTCTTCTCTCAATCTAATGAGCTGGTTGATAAGAAATGCAAGAGTGATATTGAG
AAGCTAAATCTGAAGAGAGATCAGGTGTAATCCGAGAGGTCATGATTATAGAGGTCAGTTCAAGTGAA
GAGGAAGAGAGCACCATTCAGAAGGTGATAATGTGGAAGCTGGATGCTACTGGGATGTGAAGTAGAT
GATAAAGATGATGATATCCTTCTCAACCTGTGGGATGTGAAAAGCTGTTACTGAAGGAGAAGATGGT
ATAAACTGGTCCATCAGTGACAAAAGACATTGAGGCCAGATAGCTAATAACCGAACACCTGGAAGATGG
ACCCAGCGGTACTATTACGCCAACAAAAACATTATCTGTAGAAAATGTGACAAACGTGGTCAATTTATCA
AAAAACTGCCCTTACCACGAAAAGTTCGTGCTGCTTCTCTGTGCTCCAGGAGAGGACATCTCCTGTAT
TCTGTCCAGCCCCCTTTGCGAATACTGTCTGTGCCTAAGATGTTGGACCACTCATGTCTTTTCAGA
CATTCTGGGATAAACAGTGTGACCGATGTCATATGCTAGGCCACTATACAGATGCTTGCACAGAAAATC
TGGAGGCAGTATCACCTAACGACCAAACCTGGACCACCCAAAAAGCCGAAGACCCCTTCAAGACCATCA
GCCTTAGCATATTGCTATCACTGCGCGCAAAAAGGCCATTATGGACACGAATGTCAGAAAAGAGAAGTG
TATGACCCGTCTCCAGTATCTCCATTCTGCTACTATGATGACAAATATGAAATTCAGGAGAGAGAA
AAGAGACTAAAAACAAAAATAAAAGTACTCAAGAAAAATGGGGTTATCCCAGAGCCATCCAAGCTACCT
TATATAAAAGCAGCAAAATGAGAACCCCAACCATGATATAAGGAAGGGCCGTGCCTCATGGAAAAGCAC
AGGTGGCCTCAAGAAAATAAAGAAAACAAAAAGAAATGAAGAACAAGAATAGAACTGGGAGAAGCAC
AGGAAGGCTGACAGACATCGTGAAGTGGATGAGGATTTTCCAGGGGCCCAAAACCTACTCTTCTCCT
GGCAGTTTTAAAACCCAGAAGCCTTCTAAGCCCTTTCACCGTTCATCACATTACCACACGTCAAGAGAA
GACAAGTCTCCAAGGAAGGCAAGAGGGGCAAGCAGAAGAAAAAGGAGAGGTGCTGGGAAGATGATGAC
AATGATAACTTATTTCTTATTAAGCAGAGAAAAAAAAGTCTTAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
  
```

Restriction Sites:

Sgfl-MluI

Plasmid Map:



ACCN:	NM_032226
Insert Size:	1632 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_032226.2</u>
RefSeq Size:	2701 bp
RefSeq ORF:	1632 bp
Locus ID:	84186
UniProt ID:	<u>Q8N3Z6</u>
Cytogenetics:	9p13.2
Domains:	zf-CCHC
Protein Pathways:	RNA degradation
MW:	63.1 kDa