

Product datasheet for **SC317929**

GSTCD (NM_024751) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	GSTCD (NM_024751) Human Untagged Clone
Tag:	Tag Free
Symbol:	GSTCD
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC317929 representing NM_024751.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGAAAGCCATAAAGAAAAGTCTTACAGAAGAAGAATACCTGTACCTGGACTTTTCTACCAAACAGAA
GGATGCATCTTTCTCTCATACATCTGTAACCTTTATTTCTGTTATCTTACTGTGACTGTAAAATCTTT
AAAATTTGCTTAGTTGTACCAAAGAGGTTAGTCAGTGGACCAGGCTATGTGAACCTACCATCCCTTTG
GCTATTGAGAATTTTCTCAGAGAATCTTCTGACCAGCCCCCAACTATACCTGTAGAAATACTACAGCTA
GAGAAAAAGCTTAGTGAGCCTGTAGAGTGCATAATGATGATAAACTCCGCAGGCAGAAGCTCAAGCAA
CAGAAGGCTGATGGAGTTGGCCCTCCCCTACTAAGGGAAAGGCAAGAGCAAGTCCACACACAGGAA
ACATCTGAAGGGTTGGATTCTTCATCCAAGAGTCTGGAAGTAAAAGTGGCATTCTCAAAGCTCACAGTA
CAGGAAGAACCAGCTACTACCAACAGAGAGCCTTCTCACATCAGAAAAGCAAAGCCTCCGACCTCCA
CCTCTGGAGCATGTGTTGCAGAAGGGCTTACTTCACTCTGGCAGATATTGTGCTTTGCCCTGTATC
CATCATTTCTTGGTAATTATCAGCAGGAAATTTCTGAGAAGCTGGTAGAATTTCCATTGCTAGCCTCT
TGGTACCAGAGGATTAGGAAAGTCCAGGAGTAAAACAGCAGCTCTAAGTGTGGGATCCAATTTCTC
CATTTACCAAAGTTGTTGACAACCTCAACTGAACAGCATCCAAACTTATGTGAAGTCCCAGGTGTAGAA
GAGCAAAGCGATCCTTTATTTATAGGAGGACCAAGACCAACATGGCCAAGTAAATGAAAAGGGCATT
GAAGTGATGTTTTCTCCCACCTTGGCCCTACTTGGACTCTTGATTGGAATGTTCTCCCTGCAGCAGTC
AGCCCAAAGGAAGGTAATATGTCAGTGCATCGAGCTTTGAGGAAGCAGCAACAGTTGAACAACTTGTG
TATGTGGTAACAAATCAGGCCAAACCTGGTGACAGAATTGTGGATTTCTGCAGCGGTGGGGGCCATGT
GGAATGTCTTGGCTCACATGCTGCCATCATGTCAGGTTACATTAATAGAAAACAAGGAATTATCATTAA
ATTCGTGCTAAGAAGAGAAGTGAATGAAGTGGGTTAAGCAACATTTGGTTCATTCAAGCAAAATATGGAA
TATTTTACTGGGATGTTAATATTGGAGTAGCATTGCATGCTTGTGGAGTGGCAACAGACATGGTGATT
GAGCACTGTATCAAAACACGGGCTTCTTCGTCACATGCCCTTGCTGTTATGGTTTCATTGAGAACC
TCAAAGTCAATTTCCAAAAAGTGAACAATTCAGAAAACCTTTATCATAACAAGGAACACATGATTCTG
TGCAGATTTGCAGACCAGACAGCTGTCCAGCTCCACCCCAACGAAGGCTCATAGAAAACAGTGCATG
TGCTTGGTGGATCTGGATCGAGCAAGAGCTGCAGAAGAATGTGGATACTCCGTTCAAGTGATATCCATG
GAGCCAGAGAGTGTCTCCAAAAATAACATGATTGTGGGAGTCCCATTTAA
ACGGCTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCTGGAT
TACAAGGATGACGACGATAAGGTTAAACGGCCGCGC
  
```

- Restriction Sites:** Sgfl-Mlul
- ACCN:** NM_024751
- Insert Size:** 1641 bp
- OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
- OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
- Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_024751.3](#)

RefSeq Size: 4043 bp

RefSeq ORF: 1641 bp

Locus ID: 79807

UniProt ID: [Q8NEC7](#)

Cytogenetics: 4q24

Protein Families: Druggable Genome

MW: 61.3 kDa