

## Product datasheet for **SC317932**

### **SLC22A15 (NM\_018420) Human Untagged Clone**

#### **Product data:**

|                           |                                           |
|---------------------------|-------------------------------------------|
| Product Type:             | Expression Plasmids                       |
| Product Name:             | SLC22A15 (NM_018420) Human Untagged Clone |
| Tag:                      | Tag Free                                  |
| Symbol:                   | SLC22A15                                  |
| Synonyms:                 | FLIPT1; PRO34686                          |
| Mammalian Cell Selection: | Neomycin                                  |
| Vector:                   | pCMV6-Entry (PS100001)                    |
| E. coli Selection:        | Kanamycin (25 ug/mL)                      |



[View online »](#)

**Fully Sequenced ORF:** >SC317932 representing NM\_018420.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGAGGTGGAGGAGCGTTCCAGGCGGTGGGGGAGATGGGCATCTACCAGATGTACTTGTCTCCTG
CTGGCCGTGCTGCTGCAGCTCTACGTGGCCACGGAGGCCATCCTCATTGCACTGGTTGGGGCCACGCCA
TCTACCACTGGGACCTGGCAGAGCTCCTGCCAAATCAGAGCCACGGTAACCACTCAGCTGGTGAAGAC
CAGGCCTTTGGGGACTGGCTCCTGACAGCCAACGGCAGTGAGATCCATAAGCACGTGCATTTTCAGCAGC
AGCTTCACCTCCATCGCCTCGGAGTGGTTTTTAATTGCCAACAGATCCTACAAAGTCAGTGCAGCAAGC
TCTTTTTCTTCAGTGGTGTATTTGTTGGAGTTATCTCTTTTGGTCAGCTTTTCAGATCGCTTCGGAAGG
AAAAAAGTCTATCTCACAGGTTTTGCTCTTGACATCTTATTTGCAATTGCAAATGGATTTTCCCCCTCA
TATGAGTTCTTTGCAGTAACCTCGCTTCTGGTGGGCATGATGAATGGAGGGATGTCGCTGGTGGCCTTT
GTCTTGCTTAATGAATGTGTGGGCACCGCTACTGGGCACTTGCAGGATCGATTGGCGGCCTGTTCTTT
GCAGTTGGCATTGCCAATATGCCCTGTAGGATACTTCATCCGCTCCTGGAGGACCCTAGCCATTCTG
GTTAACCTGCAGGAACGGTGGTCTTTCTCTTATCTTTATTCATTCTGAATCACCTCGTTGGTTATAC
TCCCAGGGTCGACTGAGTGAGGCTGAAGAGGCGCTGTACCTCATTGCCAAGAGGAACCGCAAACCTCAAG
TGCACGTTCTCACTAACACACCCAGCCAACAGGAGCTGCAGGGAGACTGGAAGTTTCTGGATCTCTTT
CGTTACCGGGTCTGTTAGGACACACTTGTATCCTGATGTTTCTGTTTGTGTGCAGCTTGGTGTAT
TATGGCCTAACTCTGAGTGCGGGTGATCTAGGTGGAAGTATTTATGCCAACCTGGCCCTGTCTGGCCTC
ATAGAGATTCCATCTTACCCTCTCTGTATCTACTTGATTAACCAAAATGGTTTGGTGGGAAGCGAACA
TTATCAGCATTCTGTGCCTAGGAGGACTGGCTTGCTTATTGTAATGTTTCTCCAGAAAAGAAAGAC
ACAGGTGTGTTTGCAGTGGTGAACAGCCATTCTTGCTTCTGCTGGGAAGCTGACCATCAGTGTGCC
TTTAACATTGTTTATATCTACACCTCTGAGCTTTACCCTACAGTCATCAGGAATGTTGGGCTTGGAACT
TGTTCCATGTTCTCCCGAGTTGGTGGGATTATTGCTCCCTTCATCCCTCACTGAAATATGTGCAATGG
TCTTTACCATTATTGTCTTCGGAGCCACGGGTCTGACCTCCGGCCTCCTGAGTTTGTATTGCCGGAG
ACCCTTAACAGTCCGCTGCTAGAAACATTCTCCGACCTTCAGGTGATTTCGTATCGCAGGCTGGGAGAA
GAAGCATTATCTTTACAGGCTTTGGACCCCAACAGTGTGTGGACAAGGAGAGCTCTTTAGGGAGTGAG
AGTGAGGAAGAGGAAGAATTTTATGATGCAGATGAAGAGACTCAGATGATCAAGTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
  
```

**Restriction Sites:** SgfI-MluI

**ACCN:** NM\_018420

**Insert Size:** 1644 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

|                               |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
|-------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Reconstitution Method:</b> | <ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol> |
| <b>RefSeq:</b>                | <u>NM_018420.2</u>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>RefSeq Size:</b>           | 4729 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>RefSeq ORF:</b>            | 1644 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Locus ID:</b>              | 55356                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>UniProt ID:</b>            | <u>Q8IZD6</u>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Cytogenetics:</b>          | 1p13.1                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Domains:</b>               | sugar_tr                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Protein Families:</b>      | Transmembrane                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>MW:</b>                    | 60.5 kDa                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Gene Summary:</b>          | Organic ion transporters, such as SLC22A15, transport various medically and physiologically important compounds, including pharmaceuticals, toxins, hormones, neurotransmitters, and cellular metabolites. These transporters are also referred to as amphiphilic solute facilitators (ASFs).[supplied by OMIM, Apr 2004]                                                                                                                                                                                               |