

Product datasheet for **SC317934**

SLC22A7 (NM_153320) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SLC22A7 (NM_153320) Human Untagged Clone
Tag:	Tag Free
Symbol:	SLC22A7
Synonyms:	hOAT11; NLT; OAT2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC317934 representing NM_153320.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGGGCTTTGAGGAGCTGCTGGAGCAGGTGGCGGGCTTTGGGCCCTTCCAACCTGCGGAATGTGGCACTG
CTGGCCCTGCCCGAGTGTCTACCACTGCACCTCTCTCTGCCATCTTCTGGCTGCCGTGCCCTGCC
CACCGATGTGCCCTGCCGGTGCCCTGCCAACTTCAGCCATCAGGATGTGTGGCTGGAGGCCCATCTT
CCCCGGGAGCCTGATGGCACGCTCAGCTCTGCCTCCGCTTTGCCTATCCCCAGGCTCTCCCAACACC
ACGTTGGGGGAAGAAAGGCAGAGCCGTGGGGAGCTGGAGGATGAACCTGCCACAGTGCCTGCTCTCAG
GGCTGGGAGTACGACCACTCAGAATTCCTCTACCATTGCAACTGAGTCCCAGTGGGATCTGGTGTGT
GAGCAGAAAGGTCTGAACAGAGCTGCGTCCACTTTCTTCTCGCCGGTGTCTGGTGGGGCTGTGGCC
TTTGGATATCTGTCGACAGGTTTGGGCGCGCGCTGTCTGCTGGTAGCCTACGTGAGTACCCTGGTG
CTGGCCCTGGCATCTGCAGCCTCCGTAGCTATGTAATGTTTCCATCACCCGCACCCTTACTGGCTCA
GCCCTGGCTGGTTTTACCATCATCGTGATGCCACTGGAGCTGGAGTGGCTGGATGTGGAGCACCGCACC
GTGGCTGGAGTCTGAGCAGCACCTTCTGGACAGGGGGCGTGATGCTGCTGGCACTGGTTGGGTACCTG
ATACGGGACTGGCGATGGCTTCTGCTAGCTGTCAACCCTGCCTTGTGCCCCAGGCATCCTCAGCCTCTGG
TGGGTGCCTGAGTCTGCACGCTGGCTTCTGACCAAGGCCATGTGAAAGAGGCCACAGGTAATTGCTC
CACTGTGCCAGGCTCAATGGGCGGCCAGTGTGTGAGGACAGCTTCAGCCAGGAGGCTGTGAGCAAAGTG
GCCGCCGGGAACGGGTGGTCCGAAGACCTTCATACCTAGACCTGTTCCGCACACCACGGCTCCGACAC
ATCTACTGTGCTGCGTGGTGGTGGTTCGGAGTGAACCTTCTCTATTACGGCTGAGTCTGGATGTG
TCGGGGCTGGGGCTGAACGTGTACCAGACACAGCTGTTGTTTCGGGGCTGTGGAAGTGCCTCCAAGCTG
CTGGTCTACTTGTCCGTGCGCTACGAGGACCGCCCTCACGAAGCCGGGACACTGCTGGGACCGGCC
CTGGCGTTTCGGCACTAGACTGCTAGTGTCTCCGATATGAAGTCCCTGGAGCACTGCTTGGCAGTGATG
GGAAAGCTTTTTCTGAAGCTGCCTTACCACCTGCCTACCTGTTCACTTACAGAGTTGTACCTACGGTG
CTCAGACAGACAGGGATGGGGCTGACTGCACTGGTGGGCCGGCTGGGGGGCTCTTTGGCCCCACTGGCG
GCCTTGTGGATGGAGTGTGGCTGTCACTGCCAAGCTTACTTATGGGGGATCGCCCTGCTGGCTGCC
GGCACCGCCCTCTGCTGCCAGAGACGAGGCAGGCACAGCTGCCAGAGACCATCCAGGACGTGGAGAGA
AAGAGTGCCTCAACCAGTCTTCCAGGAGGAAGAGATGCCCATGAAGCAGGTCCAGAATAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites:

Sgfl-MluI

Plasmid Map:



ACCN:	NM_153320
Insert Size:	1647 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_153320.2</u>
RefSeq Size:	2559 bp
RefSeq ORF:	1647 bp
Locus ID:	10864
UniProt ID:	<u>Q9Y694</u>
Cytogenetics:	6p21.1
Domains:	sugar_tr
Protein Families:	Druggable Genome, Transmembrane
MW:	60 kDa
Gene Summary:	<p>The protein encoded by this gene is involved in the sodium-independent transport and excretion of organic anions, some of which are potentially toxic. The encoded protein is an integral membrane protein and appears to be localized to the basolateral membrane of the kidney. Alternatively spliced transcript variants encoding different isoforms have been described. [provided by RefSeq, Jul 2008]</p> <p>Transcript Variant: This variant (2) contains an alternate segment in the 5' coding region compared to variant 1. The encoded isoform (b) is longer compared to isoform a.</p>