

Product datasheet for **SC317943**

SPTLC3 (NM_018327) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SPTLC3 (NM_018327) Human Untagged Clone
Tag:	Tag Free
Symbol:	SPTLC3
Synonyms:	C20orf38; dj718P11; dj718P11.1; hLCB2b; LCB2B; LCB 3; LCB3; SPT 3; SPT3; SPTLC2L
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC317943 representing NM_018327.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCACGATCGCC
ATGGCTAACCCGTGAGGTGGTGTCTGTTTCAACGGGAACTTCACAATCACAAGAAACAGAGCAATGGC
TCACAAAGCAGAAACTGCACAAAGAATGGAATAGTGAAGGAAGCCAGCAAAATGGGAAGCCACATTTT
TATGATAAGCTCATTGTTGAATCGTTTGAGGAAGCACCCCTTCATGTTATGGTTTTCACTTACATGGGA
TATGGAATTGGAACCTGTTTGGCTATCTCAGAGACTTTTTAAGAACTGGGGAATAGAAAAATGCAAC
GCAGCTGTGGAACGAAAAGAACAAAAGATTTTGTCCACTGTATCAAGACTTTGAAAATTTTATACA
AGAAACCTTTACATGCGAATCAGAGACAACTGGAACCGGCCATCTGCAGTGCCCGAGGGCTCTGTTT
GATTTGATGGAGAGGTATCAGACGACTATAACTGGACGTTTAGGTTTACTGGAAGATCATCAAAGAT
GTCATCAACATGGGCTCCTATAACTTCTTGGTCTTGACGCCAAGTATGATGAGTCTATGAGGACAATA
AAGGATGTTTTAGAGGTGATGGCACAGGCGTGGCCAGCACCAGGCATGAAATGGGCACCTTGGATAAG
CACAAGGAGTTGGAGGACCTTGTGGCTAAGTTCCTGAATGTGGAAGCAGCTATGGCTTTGGGATGGGA
TTCCGCAACTAACTCAATGAATATCCAGCATTAGTTGGAAGGGATGCCTCATTTAAGTATGAGTTA
AACCACACATCGCTTGTGCTTGGGGCCGACTCTCAGGTGCAACCATAAGAACTTCAAACACAACAAC
ACACAAAGCCTAGAGAAGCTCCTGAGAGATGCTGTATGATGATGATGATGATGATGATGATGATGATGAT
AAAAAGATTCTCATCCTGGTGGAGGGTGTCTACAGCATGGAAGGTTCCATCGTGCATCTGCCCCAGATC
ATAGCTCTAAAGAAGAAATACAAGGCTTACCTCTACATAGATGAAGCTCACAGTATTGGGGCCGTGGGC
CCAACCGCCGGGGTGTACGGAGTTCTTTGGACTAGACCCTCATGAAGTTGATGTGCTCATGGGCACA
TTCACAAAAGTTTTGGAGCTTCCAGGAGTTACATAGCTGGAAGGAAGGACCTCGTGGATTATTTACGG
GTTCACTCGCATAGTGTGTTTATGCTTCCATCCATGAGCCACCGATAGCAGAGCAAATCATCAGATCA
CTAAAACCTTATCATGGGACTGGATGGGACCACTCAAGGGCTGCAGAGAGTACAGCAACTTGCGAAAAAC
ACAAGATACTTCAGACAAAGACTGCAGGAAATGGGATTATTATCTATGGCAATGAGAATGCTTCTGTT
GTTCTCTGCTTCTTTATATGCCTGGTAAAGTAGCGGCTTTTGAAGGCATATGCTAGAGAAAAAATT
GGAGTGGTGGTGGGATTTCCAGCCACTCCCCTCGCAGAAGCTCGGGCTCGGTTTTGTGTTTCAGCG
GCACATACCCGGGAGATGTTAGACACGGTTTTAGAAGCTCTTGATGAAATGGGTGATCTCTTGAACCTG
AAATATCCCGGCACAAGAAGTCAGCACGTCCTGAGCTCTATGATGAGACGAGCTTTGAACTCGAAGAT
TAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```

Restriction Sites: Sgfl-Mlul

Plasmid Map:


ACCN: NM_018327

Insert Size: 1659 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_018327.2](#)

RefSeq Size: 3855 bp

RefSeq ORF: 1659 bp

Locus ID: 55304

UniProt ID: [Q9NUV7](#)

Cytogenetics: 20p12.1

MW: 62 kDa

Gene Summary: This gene encodes a subunit of the serine palmitoyltransferase complex which catalyzes the rate-limiting step in sphingolipid biosynthesis. This subunit metabolizes lauroyl- and myristoyl-CoA and generates C14 and C16-sphingoid bases. [provided by RefSeq, Mar 2017]