

## Product datasheet for **SC317953**

### UTP18 (NM\_016001) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	UTP18 (NM_016001) Human Untagged Clone
Tag:	Tag Free
Symbol:	UTP18
Synonyms:	CGI-48; WDR50
Vector:	<u>pCMV6 series</u>
Fully Sequenced ORF:	>NCBI ORF sequence for NM_016001, the custom clone sequence may differ by one or more nucleotides

```

ATGCCGCCGAGCGGAGGAGACGAATGAACTGGACCGGAGAACCGGAGCGAAGCCGAAG
CGGAAGCCCGAATGAGGCCGACTGGAAGCCGGAGCGGGCCAGGCGGGCCTCCCAA
AAGCCTGCCCTTCATCCAGCGGAAACCGCCGGCCCGGAGCGGCGGCGCTGCG
ATTGCAGTCGCGGCGGAGGAAGAGACGGCTCCGGCAGCGGAACCGCTGAGGCTG
GAGGAGGACAAACCGCCGTGGAGCGGTGCTTGGAGGAGCTGGTCTTCGGCGACGTCGAG
AACGACGAGGACGCTTGTGCGGCGTCTGCGAGCCCGAGGGTCAAGAACATGAAGAC
TCGGGTGACTCAGAAGTGGAGAATGAAGCAAAAGGTAATTTCCACCTCAAAGAAGCCA
GTTTGGGTGGATGAAGAAGATGAAGATGAGGAAATGGTTGACATGATGAACAATCGGTTT
CGGAAGGATATGATGAAAATGCTAGTGAAGTAACTTTCGAAAAGACAACCTTAAAAAG
AGACTTAAAGAAGAAATCCAACATGCCATGGGAGGAGTACCTGCCTGGGCAGAGACTACT
AAGCGGAAAACATCTTCAGATGATGAAAGTGAAGAGGATGAAGATGATTTGTTGCAAAGG
ACTGGGAATTTTCATATCCACATCAACTTCTTCCAGAGGAATCTTGAAGATGAAGAAC
TGCCAGCATGCGAATGCTGAACGTCCTACTGTTGCTCGGATCTCATCTGTGCAGTCCAT
CCCGGTGCACAGATTGTGATGGTTGCTGGATTAGATAATGCTGTATCACTATTTTCAGTT
GATGGGAAAACAAATCCTAAAATTCAGAGCATCTATTTGGAAAGGTTTCCAATCTTTAAG
GCTTGTTTTAGTGCTAATGGGGAAGAAGTTTTAGCCACGAGTACCCACAGCAAGGTTCTT
TATGTCTATGACATGCTGGCTGGAAAGTTAATCCTGTGCATCAAGTGAAGGTTTGAAA
GAGAAGATAGTGAGGAGCTTTGAAGTCTCCCCAGATGGGTCTTCTTGCTCATAAATGGC
ATTGCTGGATATTTGCATTTGCTAGCAATGAAGACCAAAGAAGTATTGGAAGCATGAAA
ATTAATGGAAGGTTGCAGCATCCACATCTCTTCAGATAGTAAGAAAGTATACGCCTCT
TCGGGGATGGAGAAGTTTATGTTGGGATGTGAACCAAGGAAGTGCCTAACAGATTT
GTTGATGAAGGCAGTTTATGATTAAGCATTGCCACATCTAGGAATGGACAGTATGTT
GCTTGTGGTTCTAATTGTGGAGTGGTAAATATATAACAATCAAGATTCTTGTCTCCAAGAA
ACAAACCCAAAGCCAATAAAAGCTATAATGAACTTGGTTACAGGTGTTACTTCTCTGACC
TTCAATCTACTACAGAAATCTTGGCAATTGCTTCAGAAAAAATGAAAGAAGCAGTCAGA
TTGGTTTACTTCTTCTGTACAGTATTTTCAAACCTCCAGTCATTAATAAAGAAT
ATTTCTCATGTTACATACCATGGATTTTTCTCCGAGAAGTGGAATACTTGCCTTGGGGAAT
GAAAAGGGCAAGGCCCTGATGTATAGGTTGCACCATTACTCAGACTTC

```



[View online »](#)

<b>Restriction Sites:</b>	Please inquire
<b>ACCN:</b>	NM_016001
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<u><a href="#">NM_016001.2</a></u> , <u><a href="#">NP_057085.2</a></u>
<b>RefSeq Size:</b>	1893 bp
<b>RefSeq ORF:</b>	1671 bp
<b>Locus ID:</b>	51096
<b>UniProt ID:</b>	<u><a href="#">Q9Y5J1</a></u>
<b>Cytogenetics:</b>	17q21.33
<b>Domains:</b>	WD40
<b>Gene Summary:</b>	Involved in nucleolar processing of pre-18S ribosomal RNA.[UniProtKB/Swiss-Prot Function]