

Product datasheet for **SC317958**

LRRC71 (NM_144702) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: LRRC71 (NM_144702) Human Untagged Clone
Tag: Tag Free
Symbol: LRRC71
Synonyms: C1orf92
Vector: pCMV6 series
Fully Sequenced ORF: >NCBI ORF sequence for NM_144702, the custom clone sequence may differ by one or more nucleotides

```
ATGTCGAGCGAGCAGAGCGCGCCGGGGCCTACCCAGGGCCCCGCTCCGGGGACCCAG
AAGTCTTCTGGCGCGGTGACAAAAAGGGAGAGCGCGCGCCAAAGAGAAGCCAGCGACC
GTTCTGCCTCCCGTGGGGGAGGAGGCCAAAAAGCCCTGAGGAGTACCAGTGCTCCGGG
GTCCTCGAGACCGACTTCGCCGAGCTCTGCACGCGGTGGGGCTACACGGACTTCCCCAAA
GTTGTCAACCGGCCCGCCCCACCCGCCCTTCGTCCTCCGCCTTTTGTGGAAAAG
GCCACCTTAGACGATCCGCGGTGTCGGGTCTGCAGCCTCAATAGCCTGGAGAGCAAA
TACGTGTTCTTCCGGCCACCATCCAGGTGGAGCTGGAGCAGGAGGACAGCAAGTCAGTG
AAGGAAATCTACATCCGCGGTTGGAAGTTGAGGAACGGATTCTGGGTGCTTCTCTAAA
TGTCGTGCCCCGCTTACCCAGCTACAGGCCATCAACTTGTGGAAGTGGGGCTGACCGAT
AAGACCCTGACCACCTTCATCGAGCTCCTGCCTCTGTTCATCCACGCTCAGGAAGGTG
TCTCTGGAGGGGAACCCACTGCCGGAGCAGTCCTATCACAAGCTCATGGCCTTGGACAGC
ACGATTGCGCACTTGTCTCTGCGGAACAATAACATCGACGACCGCGGGCGCAACTCCTG
GGCCAGGCGCTGTCCACGCTGCACAGCTGCAACCGGACCTCGTCTCGCTCAACCTGGGT
TTCAACCACATCGGTGACGAGGGCGCAGGCTACATCGCGGACGGCCTCCGGCTGAACCGT
TCCCTGCTCTGGCTGTCCCTGGCCACAACCGCATCCAGGACAAGGGCGCCCTGAAGCTG
GCTGAGGTCTGCGCGCCTTCGAGCTGACACACACCGAAGTGGTGGAGCGCCGACGCTC
CTGCTGGAAAAAGGGACACAGGAGCGCTCGCGATCGCCCTCCTCCTCTCGACACGGGGAC
TCCAAAACGGACCGTGAGAAGAGTCAGATGGTAGGGATCAGCAATAGTGCATTGGTGGAC
AAGACAGACAAGACGCAGACAATGAAAACCCCTAAGGGCTGGGCAAGAAAAAGGAGAAA
TCATGGGAATTGGCCAAGAAAGAGGAGAAGTTGGGGTCTGGGCAGTACCCACACAAGGA
ACCCCTAAGAAGGAAGATGCCACAAAGGCAGGCAAGGGGAAGGTAACCATCCCTGAACAG
AAGCCAAGCAGGGCAAAAGGGATCAAGATCGGGAGCAGAGAGAAGCGCAGCATCCTCCTG
GAGTCCGAGCTGGTTGTTGAGGCTACTGAGGTGGTCAACCCTCCTCTGGAGCCTGTGGAG
CACCGAGATGGGAAAGTTTTTCATGCCTGGGAACAAGGTCTTTTGCACCTCAACCTCATC
CGGAACCGCATCACAGAGGTGGGGCTGGAGGGTCTCCTCGCCACGGTGCAGTATCAGATG
CAGTTCTCCAAGGCCAAGAGTGCATCCAAGGGTCCAGTGGGGCTGCTGTGGCTGTCCCTG
GCTAAAAATTGCTTCGCCCCACAATGTCCTGCGTACGCCATAATCCAGGAGCTGATGTTG
CCAAGGGATCCCATCAAGGCCAACTCAGGGAGGATGAGGCCATGGCATTCTTCCC
```



[View online >](#)

Restriction Sites:	Please inquire
ACCN:	NM_144702
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_144702.2</u> , <u>NP_653303.2</u>
RefSeq Size:	1989 bp
RefSeq ORF:	1680 bp
Locus ID:	149499
UniProt ID:	<u>Q8N4P6</u>
Cytogenetics:	1q23.1