

Product datasheet for **SC317961**

ZNF34 (NM_030580) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF34 (NM_030580) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF34
Synonyms:	KOX32
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC317961 representing NM_030580.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGCTTCTGCTCCTCAGTGACCAGTCTCTGCTCACTGCCCTCAGAAAGCCCAACCCTCAGGCAATGGCG
GCCTTGTCTGTCTGCCCCACCCAGGCCGAGGTGACCTTCGAGGACGTGGCTGTACCTCTCCCGG
GAGGAATGGGGCCGCTGGGCCCTGCTCAGAGGGGCTCTACAGGACGTGATGCTGGAGACCTACGGG
AACCTAGTCTCACTGGGAGTAGGACCTGCAGGCCCAAGCCTGGAGTGATCTCGCAGTTGGAGCGAGGG
GATGAGCCCTGGGTCTGGATGTTGAGGACCTCTGGGAAAGAGCACCTGAGAGTCAACAGCCCAGCT
CTTGGGACCAGAAGTGTGACAAGGAGTTGACTTCACAGGAGACATTTGGTGAGGAAGATCCCCAGGGA
TCTGAGCCAGTAGAAGCCTGTGACCACATCAGTAAAGTCAAGGGGAGCCTGGAAAAGCTAGTGGAGCAG
AGAGGCCCCAGGGCAGTCACTGACCAACGGGGAGAGCAGCAGGGAGTCTGGGGAAACCTCAGTTG
CTGTCAAGACCTGTTCTGATCAGAGACCTCACAATGTGATATATGTGAGCAAAGTTTTGAACAGAGA
TCATATCTCAACAACCATAAGCGTGTACACAGGTCAAAAAAACAAATACAGTTCGTAACCTGGGGAA
ATCTTCAGTGAAACTTAGTTGTTAAAGAAGATCAGAAAATTCCTACTGGGAAAAAATTGCATTATTGC
AGTTACTGTGGGAAAACATTCAGGTACAGTGCCAACCTTGTCAAGCATCAGCGGCTTCACACTGAAGAG
AAGCCCTACAAATGTGATGAGTGTGGGAAAGCCTTCAGCCAGAGCTGCGAGTTTCAATCACCGAAGG
ATGCACTCAGGAGAGATTCCTACCGGTGTGACGAGTGTGGGAAGACATTCACCCGGAGGCCAACCTC
ATGAAGCACCAGAGGATTCACACTGGGGAGAAACCCTACAAGTGTGGGGAGTGTGGGAAGCACTTTAGC
GCCTACTCTTCCCTGATTTATCACCAGAGAATCCACACCGGAGAGAAACCCTATAAATGTAATGACTGC
GGGAAAGCCTTCAGTGATGGCTCAATCCTTATCCGACATCGTCGGACTCACACCGGAGAGAAGCCATTT
GAGTGCAAGGAATGTGGCAAAGGCTTTACACAAAGTTCTAACCTTATCCAACATCAGAGAATTCACACT
GGAGAGAAACCCTATAAATGTAATGAATGTGAGAAAGCTTTCATTCAAAAAACCAAACTCGTGGAACAT
CAGAGAAGCCACACTGGAGAGAAGCCCTATGAATGCAATGACTGTGGCAAAGTTTTTCAGCCAAAGCACA
CACCTCATCCAGCACCAGAGAATCCACACAGGAGAGAAGCCCTACAAGTGCAGCGAGTGTGGGAAGGCC
TTCCACAACAGTTCAGACTCATCCACCACCAGAGGCTGCACCACGGAGAGAAACCCTACAGATGCAGC
GATTGCAAGAAAGCCTTCAGCCAGAGCACGTACTTGATTGAGCACCAGGAGATCCACACCGGGGAGAG
CCCTACAAGTGCAGCGAGTGTGGGAAGGCCTTCCGGCACAGTTCCAACATGTGTCAGCATCAGCGGATT
CACCTCCGGGAGGACTTCTCCATGTAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```

Restriction Sites: Sgfl-Mlul

Plasmid Map:


ACCN: NM_030580

Insert Size: 1683 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_030580.4](#)

RefSeq Size: 2841 bp

RefSeq ORF: 1683 bp

Locus ID: 80778

UniProt ID: [Q8IZ26](#)

Cytogenetics: 8q24.3

Domains:	KRAB, zf-C2H2
Protein Families:	Transcription Factors
MW:	64 kDa
Gene Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (1) encodes the longest isoform (1). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.