

Product datasheet for **SC317962**

KLHL30 (NM_198582) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	KLHL30 (NM_198582) Human Untagged Clone
Tag:	Tag Free
Symbol:	KLHL30
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC317962 representing NM_198582.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGTGCAGAACGTGGATGACCTGGATTTCACCTGCCCTCGCATGCCAGGACATGCTGGATGGCCTG
CAGCGCCTGCGCTCTCAGCCCAAGCTGGCCGACGTCACTGCTGGTGGCGGCCGGGAGCTGCCATGC
CACCGCGGCTCCTGGCGCTCAGCAGCCCTACTTCCATGCCATGTTTGGGGTGACTTCGCCGAGAGC
TTCTCTGCGCGCGTGGAGCTGCGGGACGTGGAGCCCGCCGTGGTGGGACAACCTGGTGGACTTCGTGTAC
ACAGGCCGGCTGACCATCACGACGGGAACGTGGAGGCGCTGACACGCACGGCTGCGCGCTGCACTTC
CCCTCGGTGCAGAAGGTCTGCGGCCGTACCTGCAGCAGCAACTGGATGCCGCCAACTGCCTGGGCATC
TGTGAGTTCGGGAGCAGCAAGGCTGCTGGCGTGGCTGCCAAGGCTGGGCTTCCTGCGAGAGAAC
TTTGAGGTGTGGCAGTGGAGCAGTTCCTGCAGCTTCCCGAGAGCGGCTGGTCACTGTCTGGCC
GGCGACCTGCTGCAGGTACAGCCGGAGCAAAGCCGACTCGAGGCCCTGATGCGCTGGGTGCGCCATGAC
CCGAGGCCCGGCCGCCACCTGCCCGAGCTGCTCAGCTAGTGCACCTGGACGCCGTGCCAGGCC
TGGGTGACGAACCTGCTGGCCTCAGAGCCCTGATCCAGGAGTCAGAGGCATGCCGGGACGCCCTGTCC
CAGGGCCATGATGGGGCACCACTCGCCCTCCAGCAGAAGCTGGAGGAGTCTGTTGGTGGTGGCGGG
CAGGCGCTGGAGGAGGAGGAGGAGGAGGTGAGGAGCCACCCCGGCTTGGGAACCTTGCCTTCTACAAC
AGCAAGGCCAAGAGGTGGATGGCACTTCCAGACTTCCCGACTATCACAAGTGGGTTTCTCCCTGGCG
GCCCTGAACAACAACATCTATGTACAGGTGGCTCTCGGGGCACAAAGACAGACACCTGGTCAACCACC
CAGGCCTGGTGTCTCCCTGAAGGAGGCTCTGGAAGCCGTGGCGCCATGCTGAAGCCCGCACC
AACCACGCCAGCGCGGCCCTCAATGGGGAGATCTACGTTATCGGCGGCACCACTGGACGTGGTGGAG
GTGGAGAGCTATGACCCCTACACGACAGCTGGACGCCGTGAGCCCGCCCTCAAATACGTGAGCAAC
TTCTCGGCTGCCGCTGCGGGGCCGCTCTACCTGGTGGGCTCCAGCGCTGCAAGTACAACGCCCTG
GCCCTGCAGTGTACAACCTGTACAGATGCGTGGAGTGTGATCGCTCGCCCTTCTGCCAAGTAC
CTGTCTCGCCTCGCTGTGCTGCACTGCACGGGAGCTTACCTATTGGGGACAACACCAAGAAGTTC
TACGTGTACGACCCCGGGGCCAAGCTGTGGCAGAAGGTGAGTACAGCACAGCTGCATGAGAATGGC
GCGCTGGTGCCACTGGGTGATGCGCTGTACGTGACGGCGGCCGCTGGCAGGGCATGGAAGGTGACTAC
CACGTGGAGATGGAGGCTACGACACGGTTCGGGACACCTGGACCCGCCACGGCGCCCTGCCCGGCTC
TGGCTCTACCACGGGCTCCACCGTCTTCTGGATGTCTCAAGTGGACCCAGCCCTCCGGCCCCACC
CAGGAGCACTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_198582

Insert Size: 1737 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_198582.3](#)

RefSeq Size: 3685 bp

RefSeq ORF: 1737 bp

Locus ID: 377007

UniProt ID: [Q0D2K2](#)

Cytogenetics: 2q37.3

MW: 64 kDa