

Product datasheet for **SC317967**

ZNF257 (NM_033468) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF257 (NM_033468) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF257
Synonyms:	BMZF-4; BMZF4
Vector:	<u>pCMV6 series</u>
Fully Sequenced ORF:	>NCBI ORF sequence for NM_033468, the custom clone sequence may differ by one or more nucleotides <pre> ATGGGACCACTGACAATTAGGGATGTGACTGTAGAATTCTCTCTGGAGGAGTGGCATTGC CTGGACACTGCACAGCAGAATTTATATAGGGATGTGATGTTAGAGAACTACAGAAACCTG GTCTTCTGGGTATTGCTGTCTCTAAGCCAGACCTGATCACCTGTCTGGAGCAAGGAAAA GAGCCCTGTAATATGAAGAGACATGAGATGGTAGCCAAACCCCAAGTTATGTGTTCTCAT ATTGCTGAAGACCTTTGCCAGAGCGAGACATAAAATATTTTTTCCAAAAAGTCATACTG AGGAGATATGATAATGTGAACATGAGAATTTACAATTAAGAAAGGGCTGTAAAAGTGTG GATGAGTGAAGGTGTGCAAAGGAGGTTATAATGGACTTAACCAATGTCTGATAACTACC CAGAGCAAAATGTATCAATGTGATAAATATGTAAAAGTCTTCTATAAGTTTTCAAATTCA GATAGACATAAGATAAGACATACTGAAAAGAAAAGCTTGCAATGTAAAGAATGTGGCAAA TCATTTTGCATGCTTTCACAACTAACTCGACATAAGAGAATTCATATTAGAGAGAATTCC CACAAATGTGAAGAATGTGGCAAAGCCTTTAACCAGTCCTCAGCTCTTACTCGACATAAG ATGACTCATACTGGAGAGAAACCTACAAATGTGAAGAGTGTGGAAAAGCTTTTAACCGG TCTTCACACCTTACTCAACATAAGGTAATTCATACTAGAGAGAAACCTACAAATGTGAA GAGTGTGGCAAAGCCTTTAACCAGTCTTCACACATTACTCAACATAAGAGAATTCATAAT AGAGAGAAGCCCTTCAAATATGATGAATGTTGCAAAGCCTTTAAGTGGTCCTCAGCTCTT ACTACCCTTACTCAACATAAGAGAATTCATACTGGAGAGAAACCTACAAATGTGAAGAG TGTGGCAAAGCCTTTAACCAGTCCTCAGCCCTTACTCGACATAAGATGATTCTACTGGA GAGAAACCTTCCAATGTGAAGAGTGTGGCAAAGCTTTTAACCGGTCTTCACACCTTACT CAACATAAGATAATTCATACTAAAGAGAAACCTACAAATGTGAAGAGTGTGGAAAAGCC TTTAACCGGTCTTCACACCTTACTAAACATAAGAGAATTCATACTAGAGAGAAGGCTAC AAATGTGATGAATATTGCAAAGCTTTTAAGTGGTCCTCAGCTCTTACTACCCTTACTCAG CATAAGATAATTCATACTGGAGAGAAACCTACAAATGTGAAGAGTGTGGCAAAGCCTTT AACCGGTCTTCATACCTATTTCGACATAAGATAATTCATACTGGAGAGAAACCTACAAA TGTGAAGAGTGTGGCAAAGCCTTTAACCAGTCTTCACACCTTACTCAACATAAAATAATT CATACTGGGGAGAAGCCCTACAAATGTGAAGAATGTGGCAAAGCCTTTAACCGGTCTTCA CACCTTTCTCAACATAAGATAATTCATACTGGAGAGAAACCTACAAATGTGAAGAATGT GGCAAAGCCTTTTAATCGTTTCTCATACCTTACCGTACATAAGAGAATTCATGCTGGAGAG AACCCCAACAAATATGAAGAATGTGGCAAAGCTTGTAAACCATTCTCAAACCTTACTAAA CATAATTCA </pre>



[View online »](#)

Restriction Sites:	Please inquire
ACCN:	NM_033468
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_033468.2</u> , <u>NP_258429.2</u>
RefSeq Size:	3534 bp
RefSeq ORF:	1692 bp
Locus ID:	113835
UniProt ID:	<u>Q9Y2Q1</u>
Cytogenetics:	19p12
Domains:	KRAB, zf-C2H2
Protein Families:	Transcription Factors
Gene Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (1) encodes the longest isoform (a).