

Product datasheet for **SC317981**

MROH9 (NM_025063) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	MROH9 (NM_025063) Human Untagged Clone
Tag:	Tag Free
Symbol:	MROH9
Synonyms:	ARMC11; C1orf129
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC317981 representing NM_025063.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGTTGACAAGGAATCCAAAAACAAAGAGTAGTCTCCAGATTCTGCAAGACAGTGTGAAATGGCACCAC
ATGGCACATAAAGTTAACAGCCTATTGGATGCATACTCAGGCCTGTTAAGTAATGAATCCATGATCCTG
GCTGTGAAGTCCAGCTTTGTGGATCCCTTACTGCAGTTTGAATCTCAGTTGAAGATAATAGAGTCATCC
TTTGAATGCTAGTTGTCTGATGCCAAGTCTTGACAAAGTAAAAGAAATGGGAGCAGTTATGAGTACATT
GAGGACATGGAGAATTTATACCACAATATTTTAAATATATATGAGAACATTCTTACGAGCTTGGTGTCT
AAAGACCTCTACAACTACAGATCTTAAAGGAAATGCTCGTGTGGATGAGTAAAGATAGCTCATATCTG
CAAGAGAGAATAATGGTGATCATCAACAAGGTGTTAAGATTTACAGTCACAAAAGTCAGAAAATACATA
AGTGTGATGCTCCATGTTTGGGTCTCCTGGCAGCAGAGCTGTCTCTTTGTGTTCCCATGAAGATCCC
TCGATTGTAACAAGCATCATTGGGAATGTGTCACCTCCTCTACATTGCACGGTGTGAGAACGATATT
GGAACAAATAAGCCTACAAACGGTAAAAGTCATAGCCTCCAGTTTCTTCTCTGATGTAGAATTTCTA
CCCAAGGAGTTTCAACAAGACGAAAGTAAAATAGCTCAGAGAGTAGGGCAAACCTTACTGCCTCCCTTG
CTGACTGACTTTGTGCAGAGTCTCCTGATGAAACTCTCTTACCTGATGATAAAATCGCATCTGATGCA
GCATCCATACTGATATTTACTCTGGAATTTTATGCCGAGAAGGTCACCATGGTGTCTAAGATCGTGGAT
GCTATTTACAGGCAACTGTGTGATAACAATTGTATGAAGGATGTTATGTTGCAGGTTATCACTTTGTTG
ACATGCACTTCACCCAAGAAGGTCATCTTTCACTTATGGACTACCCAGTTCCAGCAGACGACACTCTG
ATACAGATGTGGAAGGCGGCATGTTCTCAGGCGAGCGTGGCCCTCAGTGCTGAAGACAATCTTATTG
ATACTGAAAGGAAAGCCTGGGGAAATGGAGGACACCGTAACGGAAGGGAACGTTTCTCTCTGATATT
ACCAACTTGATGCCTTTGGCGGCATGCCAGGCCCTGTGCACCTTCTGCCTCTTGGTTCCTACAGGAAA
GCGGTGGCCAGTATTTCCCCAGCTCTTGACGACTCTTATGTTCCAAGTCTTCTACAACAGTGAGCTG
AAACCGATACTCAAGGACAGGGCTTTGTATGCCAGGATGCCCTGAGAGTTCTGTGAATTGTTCTGGA
CTGCAACAGGTGGATATTACTCTAATGAAGGAGAATTTCTGGACAGTTATCTGAAGATCTGTGTAC
TATCATGGAGTCTGCTTTATTGTAAACTCTCAGTGAATATAACTTTCCACAGTTTCCGGAGACCTG
AGTTATCTCTATAAGCTCTCAGTAGAAGGTCTAGAAGGTGAGAAGACTGTATCGTATTAATATTC
CTCACTGAAGTGAGTTTGTAGACTGTGAACAGCTGTGTTCTCATTTCCTTTCTACCCAAATTTAAA
TCAAAATTCAGTTTCTCGTATCTCTCCCTCTGAATGTTGGTTCTTACCAAGACTTAAGATCATGA
ACGGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_025063

Insert Size: 1722 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_025063.2](#)

RefSeq Size: 2258 bp

RefSeq ORF: 1722 bp

Locus ID: 80133

UniProt ID: [Q5TGP6](#)

Cytogenetics: 1q24.3

MW: 65.1 kDa