

Product datasheet for **SC317992**

INSC (NM_001031853) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	INSC (NM_001031853) Human Untagged Clone
Tag:	Tag Free
Symbol:	INSC
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC317992 representing NM_001031853.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGAGACGGCCCCCTGGCAATGGAGAGGGCGCCAGCGAAGGTCCAGGTGGCTGGGTCTATGGGAGTC
CAGGAGTCCAGGAGGCTGTGCTGTGCCGGTCACGACCGCTGCAAGCAGGCTTTGCTGCAGATTGGGATC
AACATGATGGCACTGCCTGGAGGTCGCCACCTGGACTCCGTACCCTGCCGGGTACAGCGGTACACCTG
ATGCAGGTGGACTCAGTCCAGCGCTGGATGGAAGATCTGAAGCTCATGACCGAGTGCAGTGCATGTGT
GTCCTGCAGGCCAAGCCCATCAGCCTGGAAGAGGATGCACAGGGTGACCTCATCCTGGCAGGTGGCCCT
GGCCCTGGAGACCCCTGCAGCTGCTGCTCAAACGGGGTTGGGTGATTAGCACAGAGCTGCGCAGGATC
GGGCAGAAGCTGGCCAGGACCGCTGGGCACGGGTGCACAGCATGAGCGTGCCTGACCTGCCATGCC
CGCTCCATGGTCAGCGAGTACAGTGTGTGACGAGGAACCTCTGAAGGAAATGGCGGAGATTGAGAAG
CTGCTAATGGAGAAATGCTCGGAGCTCTCGGCAGTACAGAGAGGTGCCTTCAGGTTGAGAAATGAGCAT
GTCCTGAAGTCAATGAAGCCTGCGTGAGTGCAGCCCTGAGCATGCTGGCCAGCACTTTGGTCAGCTG
CTGGAGCTGGCCCTGACACGGGAGGTTACGGCACTGGTGAGAAAAATTGATGCCTCAGACAAATATCTAC
ACCACAGAGTCCACCACAGGGAACCTGTTGACCTGACCCAGGAGGGGGTCCCTTGTGCCGCATCATA
GCCAAGGAGGGTGGGGTTCGTAGCACTCTTCAAGGTTTGCCGGCAGGACAGTTTCCGGTGCTGTACCCC
CAGGCGCTCCGCACGCTGGCCCTCATCTGCTGCGTGGAAGAGGGTGTCCACCAGCTGGAGAAGGTGGAT
GGCGTTCTGTGCTTGGCCGACATCCTGACCGACAACAGCCACTCAGAGGCCACACGGGCTGAGGCTGCG
GCTGTGGTGGCCAGGTACCTCCCCACACCTGCCCGTACCCAGCACCTCAGTAGCTTCTGGAGAGC
ATGGAGGAGATCGTGACAGCCCTCGTCAAACGTGCAAGAGGCCTCATCAGGGGAAGTCTTCTACTG
GCCTCTGCGGCCCTTGCCAACATCACGTTCTTTGACACAATGGCCTGCGAGATGCTCCTGCAGTTGAAT
GCCATCCGTGTTCTCCTGGAAGCCTGCAGTGACAAGCAGAGAGTGGACACGCCTTACACTCGGGACCAG
ATTGTGACCATCTTGGCAAACATGTCTGTCTAGAACAGTGTGCCTCTGACATCATTAGGAAAAATGGG
GTCCAGCTTATCATGGGCATGCTGTCTGAAAAACCAAGGTCTGGGACTCCTGCTGAAGTGGCAGCCTGT
GAGCGAGTCCAGCAGAAAGCTGCAGTGACCCTGGCTCGTCTCAGCCGAGACCCAGATGTGGCACGGGAG
GCCGTGCGGCTCAGCTGCATGTCCCGTCTCATCGAGCTCTGCAGATCCCCATCAGAGAGGAACAGCAGT
GACGCCGTGCTTGTGGCCTGCCTGGCTGCTCTGCGTAGATTGGCTGGGTCTGCCCTGAAGGCCTCCAG
GACTCTGACTTTCAGCAGTTGGTCCAGCCTCGGCTGGTGGACTCCTTCTACTCTGCAGCAACATGGAG
GAGAGTTTGTGTAG
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_001031853

Insert Size: 1740 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001031853.4](#)

RefSeq Size: 2970 bp

RefSeq ORF: 1740 bp

Locus ID: 387755

UniProt ID: [Q1MX18](#)

Cytogenetics: 11p15.2

MW: 63.5 kDa

Gene Summary: In *Drosophila*, neuroblasts divide asymmetrically into another neuroblast at the apical side and a smaller ganglion mother cell on the basal side. Cell polarization is precisely regulated by 2 apically localized multiprotein signaling complexes that are tethered by Inscuteable, which regulates their apical localization (Izaki et al., 2006 [PubMed 16458856]).[supplied by OMIM, Mar 2008]

Transcript Variant: This variant (1) encodes the longest isoform (a).