

Product datasheet for **SC318005**

HEXD (NM_173620) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	HEXD (NM_173620) Human Untagged Clone
Tag:	Tag Free
Symbol:	HEXD
Synonyms:	HEXDC
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC318005 representing NM_173620.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGTCAGGTTCCACTCCATTTAGATGAGATTAGTTTACCTTAAAGGAGCTCCACCAAAGGTC
TCGTACCTCTCAGAGATTTTCCCTCTGTTCCGTGCGCTAGGTGCAAACGGCCTCCTCATTGAGTATGAA
GACATGTTTTCCCTACGAGGGCCCTCTGAGGCTGCTGAGGGCCAAGTACGCCTACAGCCCCCTCTGAAATC
AAAGAGATCTTGATCTGGCTGGACTCAATGAGCTGGAGGTGATTCCCTTGGTGCAGACATTTGGACAC
ATGGAGTTTGTGCTGAAGCACACGGCCTTCGCCACCTGCGGGAGGTGGGCTCCTTCCCCTGCACCTG
AACCCCCACGAGGCAGAGTCCCTGGCGCTGGTGGCGCCATGATTGACCAGGTCTGGAGCTACACCCA
GGCGCCACGCGGTGCACATCGGTGTGATGAGGTCTATTACCTCGGAGAGGGGGAGGCTCGCGCCGG
TGGCTACAGCAAGAGCAGAACAGCACGGGGAAGTTGTGCCTGTACACATGCGGGCGGTGCCAGCGGC
GTGAAGGCCCGGCCCGCCAGCGTGACACCCCTGGTGTGGGACGACATGCTCCGAGACCTGCCTGAGGAC
CAGCTCGCAGCGTCCGGGTGCCGAGCTGGTGGAGCCGTGCTCTGGGACTACACGGCCGACCTGGAT
GTGCACGCAAGTCTCCTCATGCAGAAGTACCGGCGGTGCGGCTTTCCGACGCTGTGGGCAGCCAGT
GCCTTCAAGGGTGCCACGGGGCCAGCCAGGCCGTGCCCTGTTGAGCACCACTCAGGAACACAGTG
CAGTGGCTGCAGGTGGCGGGCAGCGGGCCACGGACTCACTGCAGGGCATCATCTGACCGGCTGGCAG
AGGTACGACCACTACTCTGTGCTGTGCGAGCTGCTGCCCGCAGGAGTCCCGTCCCTGGCCGCTGCCTG
CAGTTGCTTCTACGCGGAGGATTTGATGAAGATGTTAAAGCGAAAGTGGAGAACCTTCTCGGGATTTCC
AGCCTGGAAAAACGGACCCTGTTAGGCAAGCACCTGCAGCCCTCCCTGTCCCCTTCTTCCCCTCCCC
TTCCCCCGCCCGTGGAGACAGCTGTTCTCAGCAGGGCTCTCCGAGGGAGGGGGCCGGCTCCTTCCCTG
GCAGCAACATCCTTGCCCTTGTACACAAGTCAGCCTCCATCTGCGCAGCTCTGTGGATGCGCTGCTGG
AGGGCAACAGGTATGTCACTGGCTGGTTACGCCCTACCACCGCCAGCGGAAGCTCATCCACCCGGTCA
TGGTTCAGCACATCCAGCCCGCAGCGCTCAGCCTCCTGGCACAGTGGAGCACCTCGTGCAGGAGCTGG
AGGCTGCCCTGCAGCTGGCTTTTACCCGGATGCCGTGGAGGAGTGGCTGGAGGAAAACGTGCACCCCA
GCCTGCAGCGGTGCAAGCTCTGCTGCAGGACCTCAGCGAGGTGTCTGCCCCCGCTGCCACCCACCA
GCCCTGGCAGGGACGTTGCTCAGGACCCCTGAGGGGAGAGCTCATGCCAGGGGGCTCCTGCTGGAGGCT
GGGGGGGCTCTGCACTGCCAAATGGCTGGGCAATACGGGCCACGTGGGCGTCTGCCCCTCTGGCCCA
GCAGTGTCTTGCCCACTCAGTTCTGAGGGCCCTGGGCAGCCCTGGGGGAGAGACTAGAAAACACA
GAAGGAAGCAGCACAGGGAGACCCGCTTTGTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: Sgfl-MluI
ACCN: NM_173620
Insert Size: 1758 bp

OTI Disclaimer:	Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery. The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_173620.2
RefSeq Size:	2272 bp
RefSeq ORF:	1758 bp
Locus ID:	284004
UniProt ID:	Q8WVB3
Cytogenetics:	17q25.3
MW:	63.5 kDa
Gene Summary:	Has hexosaminidase activity. Responsible for the cleavage of the monosaccharides N-acetylglucosamine (GlcNAc) and N-acetylgalactosamine (GalNAc) from cellular substrates. Has a preference for galactosaminide over glucosaminide substrates (PubMed:27149221). [UniProtKB/Swiss-Prot Function]