

Product datasheet for **SC318013**

ATG16L1 (NM_017974) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ATG16L1 (NM_017974) Human Untagged Clone
Tag:	Tag Free
Symbol:	ATG16L1
Synonyms:	APG16L; ATG16A; ATG16L; IBD10; WDR30
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC318013 representing NM_017974.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGTCGTCGGGCCTCCGCGCCGCTGACTTCCCCCGCTGGAAGCGCCACATCTCGGAGCAACTGAGGCGC
CGGGACCGGCTGCAGAGACAGGCGTTCGAGGAGATCATCCTGCAGTATAACAAATTGCTGGAAAAGTCA
GATCTTCATTAGTGTTGGCCAGAAACTACAGGCTGAAAAGCATTGACGTACCAAACAGGCACGAGATA
AGTCCCGGACATGATGGCACATGGAATGACAATCAGCTACAAGAAATGGCCCAACTGAGGATTAAGCAC
CAAGAGGAACTGACTGAATTACACAAGAAACGTGGGGAGTTAGCTCAACTGGTGATTGACCTGAATAAC
CAAATGCAGCGGAAGGACAGGGAGATGCAGATGAATGAAGCAAAAATTGCAGAATGTTTGCAGACTATC
TCTGACCTGGAGACGGAGTGCCTAGACCTGCGCACTAAGCTTTGTGACCTTGAAAGAGCCAAACCAGACC
CTGAAGGATGAATATGATGCCCTGCAGATCACTTTACTGCCTTGAGGGGAAAACAGGAAAACACG
GAAGAGAACCAGGAGCTGGTCACCAGATGGATGGCTGAGAAAGCCAGGAAGCCAATCGGCTTAATGCA
GAGAATGAAAAGACTCCAGGAGCGGCAAGCCCGGCTGCAGAAAGAGCTTGCAGAAGCAGCAAAGGAA
CCTCTACCAGTCGAACAGGATGATGACATTGAGGTCATTGTGGATGAAACTTCTGATCACACAGAAGAG
ACCTCTCCTGTGCGAGCCATCAGCAGAGCAGCCACGAGACGCTCTGTCTTCTTCCCAGTCCCCCAG
GACAATGTGGATACTCATCCTGGTTCTGGTAAAGAAGTGAGGGTACCAGCTACTGCCTTGTGTGTCTTC
GATGCACATGATGGGGAAGTCAACGCTGTGCAGTTCAGTCCAGGTTCCCGGTTACTGGCCACTGGAGGC
ATGGACCCGAGGGTTAAGCTTTGGGAAGTATTTGGAGAAAAATGTGAGTTCAAGGGTCCCTATCTGGC
AGTAATGCAGGAATTACAAGCATTGAATTTGATAGTGTGGATCTTACCTCTTAGCAGCTTCAAATGAT
TTTGCAAGCCGAATCTGGACTGTGGATGATTATCGATTACGGCACACACTCACGGGACACAGTGGGAAA
GTGCTGTCTGCTAAGTTCCTGCTGGACAATGCGCGGATTGTCTCAGGAAGTCACGACCGGACTCTCAA
CTCTGGGATCTACGCAGCAAAGTCTGCATAAAGACAGTGTTCAGGATCCAGTTGCAATGATATTGTC
TGCACAGAGCAATGTGTAATGAGTGGACATTTTGACAAGAAAATTCGTTTCTGGGACATTCGATCAGAG
AGCATAGTTTCGAGAGATGGAGCTGTTGGGAAAGATTACTGCCCTGGACTTAAACCCAGAAAGGACTGAG
CTCCTGAGCTGCTCCCGTATGACTTGTAAAAGTTATTGATCTCCGAACAAATGCTATCAAGCAGACA
TTCAGTGCACCTGGGTTCAAGTGCAGCTCTGACTGGACCAGAGTTGTCTTCCAGCCCTGATGGCAGTTAC
GTGGCGGACAGGCTCTGCTGAGGGCTCTGTATATCTGGAGTGTGCTCACAGGAAAGTGGAAAAGGTT
CTTTCAAAGCAGCACAGCTCATCCATCAATGCGGTGGCGTGGTCGCCCTCTGGCTCGCACGTTGTCAGT
GTGGACAAGGATGCAAAGCTGTGCTGTGGGCACAGTACTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_017974

Insert Size: 1767 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_017974.3](#)

RefSeq Size: 3354 bp

RefSeq ORF: 1767 bp

Locus ID: 55054

UniProt ID: [Q676U5](#)

Cytogenetics: 2q37.1

Domains: WD40

MW: 66.3 kDa

Gene Summary: The protein encoded by this gene is part of a large protein complex that is necessary for autophagy, the major process by which intracellular components are targeted to lysosomes for degradation. Defects in this gene are a cause of susceptibility to inflammatory bowel disease type 10 (IBD10). Several transcript variants encoding different isoforms have been found for this gene.[provided by RefSeq, Jun 2010]
Transcript Variant: This variant (2) lacks an alternate in-frame exon compared to variant 1. The resulting isoform (2) has the same N- and C-termini but is shorter compared to isoform 1.