

Product datasheet for **SC318017**

NCKAP5 (NM_207481) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	NCKAP5 (NM_207481) Human Untagged Clone
Tag:	Tag Free
Symbol:	NCKAP5
Synonyms:	ERIH1; ERIH2; NAP5
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC318017 representing NM_207481.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGAGGGAAAGAGACAGCTTGAGAAAAGGACTTTGAAAAAGGCTGTCTCTAGACAGCAGTCTTG
GAATACATGGACTCCAATAAATACATTGAGCATCTGCTGACTCAGCTTGAGGAGCAACACAGGAGTCTC
TGGAGGGAGAAGTTGGCAGTGGCCCGACTACAGCGAGAAGTTGCCCAAAGAACAAAGTGAGGGGGCCATG
CATGAGAAGCTGATACATGAAGTGAAGAGGAGAGACACTTACGTCTTCAAAGCGAGAAGCGGTTGCAG
GAGGTGACCTAGAGTCTGAACGCAACAGAATTCAGATGCGTAGCTTGACGACGAGTCTCCAGGATG
GAAGAAACAGTACGAAATCTATTGCAGAGTCAAGGATCTCCAGAGCAGAAAAAGAAGAACTGTTAAT
ATAATGGTCTATCAGGAAAAGCTGTGAGAGGAAGAGAGAAAACATAAGGAAGCTTTGGAAGATCTTAC
ATGGTGGTTGATGAGGATTCGAGGAGTGAAGCAGCAGTACAGATGAGGGAAAAGAAAAGACCAATTG
CTATTAGAGAGATTGAAAGCTCTAGAGGCAGAGAATTCAGCGTTGGCTTTGGAGAATGAAATCAAAGG
GAACAATATGAGCGATGTCTTGATGAGGTAGCCAACCAAGTTGTTCAAGCTCTGCTTACTCAAAGGAT
CTAAGAGAGGAATGTGTGAAGTTGAAAACAAGAGTGTTTGATTTGGAACAGCAGAATCGAACACTAAGC
ATCCTATTCCAACAGCGAGTCAACCCACTTCTGATCTGCTCCTCCAGAACTTCACTCAGTCTCTTG
GATCTTTTCATCTGGAGATTTGCTTTAGAGGTGGAAGAAACCGAAGTCTGACACAGTCAAGAACTGAT
GCTGAGGTGCACGAACACCAACTAAATACAAAATCGGCCCTTGAATGCCCTGGCTTGGGGGCTGTCATC
CCTGGTCACTCTGTCTCGAAACAGCTACAGTAGCAGCAGTGAAGTGTCTTTCAAGCACCTGCAGC
GAGTACTCCAGTGGCTCCTCCTACAGTGGCAGCATGGGAAGAACCTCAGGAAAAGGGTAAAAGAGGCA
AATGAAAACCTTGAAGAGGATGAAGACGATGCAGTTGCAGATTCTGTATTTAGAGCCACATCATAGAA
TCCAACCTGCCAGATGAGAACATTGGACAGTGGGATCGGAACCTTCCACTCCAGACTCGGGAAATCGC
TCGACAGGACGCTACCTATGCCAGCCAGACTCCCGAGAGGACGCTGAGCCTCTCCTGCTCTCCAGTCA
GCCCTTTCTGAGTTTCTTCCATGAGAGCCCAACCTTGAACGTGAAGTGCCTTCTCCACAGACGGC
CAGCGCCCTGCAGATAGCGCCATTGTTCAATCCACATCCGACCCCATCATGACCGCCAGAGGGATGAGG
CCTCTTACAGCCGCTCCCAAACAGCTTCTCAGGAAAAGTCAAGTCCCAAAGCAGAATGAAGCA
GAGCCAAGGCCTCAGACATGCTCATCTTCGGATATGCTGAAGACCCAATGGCAAGCCAGCCGCTTCCA
GACTGGGGGAGTGAAGTTGCTGCCACGGGACCCAGGACAAGGCACCCAGAATGTGTACGTACTCTGCC
AGCGGTGGCAGTAATAGTGACAGTGACCTGGACTATGGAGATAATGGTTTGGAGCTGGAAGGGGACAG
TTAGTGAAGCACTGAAGAGCGCTGCCCAAGAAATTGAGACAATTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_207481

Insert Size: 1773 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_207481.3](#)

RefSeq Size: 3654 bp

RefSeq ORF: 1773 bp

Locus ID: 344148

UniProt ID: [O14513](#)

Cytogenetics: 2q21.2

MW: 66.4 kDa