

Product datasheet for **SC318028**

DCAF13 (NM_015420) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	DCAF13 (NM_015420) Human Untagged Clone
Tag:	Tag Free
Symbol:	DCAF13
Synonyms:	GM83; HSPC064; Sof1; WDSOF1
Vector:	<u>pCMV6 series</u>



[View online »](#)

Fully Sequenced ORF:	>NCBI ORF sequence for NM_015420, the custom clone sequence may differ by one or more nucleotides <pre> ATGCAGTGGCCGCCACCGTGGCGACGGTCGCCCCCTTGTGAGCGCAACCCACCTCCGGGA CCAACGAGAGGACTCTTATGCCCAAGGCGCGTGAGCGAAGCCGGAGACTCTAGTACTGAG GGGGCAAGAACGGGGCACAGCCACTACGTACGCCGGCGGAAAGCGCGATGGAGGGGAGG TCCCAGCCTCTCCCGAGTCTCCGCCCTGCCTCGCCACAAATGCTCCAGCCCACGGACGC GGAAGCGACTGCACAGCTAAAGAACC GCCCTTGCTAACGTTTGCCGGCAAAACTACGGA GCGGCGCAGGATCCAGTTGAGCTAGCAGGCCGCCGCTCCGCGAGTCACGTGACTGGAAGC AGTCTGGGAAAAGCGGAAGTCGCCTGTGGGAGGAGGTGGCGGTGGGCGGAACCTCAGCG GACACCTCGTGGAGTCCGGCCGGAAGAGCAACCGAGATGAAGGTGAAGATGCTGAGCCGG AATCCGGACAATTATGTCCGCAAAACCAAGTTGGACTTACAGAGAGTTCCAAGAACTAT GATCCTGCTTTACATCCTTTTGGAGTCCCACGAGAATATGTAAGAGCTTTAAATGCTACC AAACTGGAACGAGTATTTGCAAAACCATTCTTGCTTCGCTGGATGGTCACCGTGATGGA GTCAATTGCTTGGCAAAGCATCCAGAGAAGCTGGCTACTGTCTTTCTGGGGCGTGTGAT GGAGAGGTTAGAATTTGGAATCTAACTCAGCGGAATTGTATCCGTACAATACAAGCACAT GAAGGCTTTGTACGAGGAATATGTACTCGCTTTTGTGGGACTTCTTTTTCTACTGTTGGT GATGACAAAACCTGTGAAGCAGTGGAATGGATGGGCCAGGCTATGGAGACGAGGAAGAG CCATTACATACAATATTAGGAAAGACAGTGTATACTGGGATTGATCATCTGGAAGAA GCTGTTTTTGGCAGATGTGGACAGCAAGTAGACATTTGGGATGAACAAAGAACTAATCCT ATATGTTCAATGACCTGGGGATTTGACAGTATAAGTAGTGTTAAATTTAACCCAATTGAG ACATTTCTCTTGGGAAGTTGTGCATCTGACAGGAATATAGTACTGTACGATATGAGGCAA GCTACTCCTTTGAAAAAGTTATCTTAGATATGAGAACAATAACAATCTGTTGGAACCTT ATGGAAGCTTTTCATTTTACAGCAGCAAATGAAGATTATAACTTATATACTTTTGATATG CGTGCACTGGACACTCCTGTAATGGTCCATATGGATCATGTATCTGCAGTGCTTGATGTG GATTACTCTCCCACTGGGAAGGAGTTTGTGTCTGCTAGTTTCGATAAATCTATTCAATC TTTCTGTAGACAAAAGTCGAAGCAGGGAGGTATATCATACAAAGAGAATGCAACATGTT ATCTGTGTAAAATGGACTTCTGACAGCAAGTATATTATGTGTGGATCTGATGAAATGAAC ATTGCGCTGTGGAAAGCTAATGCTTCTGAAAAATTGGGTGTGCTTACATCACGAGAAAAA GCAGCCAAGGATTATAACCAGAAATTGAAGGAGAAATTCAGCATTATCCTCATATAAAA CGTATAGCTCGTCATCGACATCTACCAAAATCTATCTATAGCCAGATTGAGAACAGCGC ATCATGAAAGAAGCTCGTCGACGAAAGGAGTGAATCGTATTAACACAGCAAGCCTGGA TCTGTGCCACTTGTGTCAGAGAAGAAGAACACGTAAGTGGCAGTTGTAAAA </pre>
Restriction Sites:	Please inquire
ACCN:	NM_015420
Insert Size:	2657 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_015420.5, NP_056235.3</u>
RefSeq Size:	2657 bp
RefSeq ORF:	1794 bp
Locus ID:	25879
UniProt ID:	<u>Q9NV06</u>
Cytogenetics:	8q22.3
Domains:	WD40, Sof1
Gene Summary:	Possible role in ribosomal RNA processing (By similarity). May function as a substrate receptor for CUL4-DDB1 E3 ubiquitin-protein ligase complex.[UniProtKB/Swiss-Prot Function]