

Product datasheet for **SC318041**

SLC26A11 (NM_173626) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SLC26A11 (NM_173626) Human Untagged Clone
Tag:	Tag Free
Symbol:	SLC26A11
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC318041 representing NM_173626.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGCCTTCTTCGGTGACGGCGCTGGGTACAGCCAGGTCTCTGGCCCCGGGATGGCCCCAGCGCCTGC
TGCTGCTCCCTGCGGCCCTGCAGAGGAGGTGCCATCCTGGCGTGGCTGCCAGTACTCCCTGCAG
TGGCTGAAGATGGATTTGTCGCCGGCCTCTAGTTGGCCTCACTGCCATTCCCCAGGCGCTGGCCTAT
GCTGAAGTGGCTGGACTCCCGCCCCAGTATGGCCTCTACTCTGCCTTCATGGGCTGCTTCGTGATTTTC
TTCCTGGGCACCTCCCGGGATGTGACTCTGGGCCACCACCATTATGTCCCTCCTGGTCTCCTTCTAC
ACCTTCCATGAGCCCGCTACGCTGTGCTGCTGGCCTTCTGTCCGGCTGCATCCAGCTGGCCATGGGG
GTCCTGCGTTTGGGGTTCCTGCTGGACTTCATTTCTACCCCGTCATTAAGGCTTACCTCTGCTGCT
GCCGTACCATCGGCTTGGACAGATCAAGAACCTGCTGGGACTACAGAACATCCCCAGGCCGTTCTTC
CTGCAGGTGTACCACACCTTCTCAGGATTGCAGAGACCAGGGTAGGTGACGCCGTCTGGGGCTGGTC
TGCATGCTGCTGCTGCTGGTGTGAAGCTGATGCGGGACCACGTGCCTCCCGTCCACCCGAGATGCC
CCTGGTGTGCGGCTCAGCCGTGGGCTGGTCTGGGCTGCCACGACAGCTCGCAACGCCCTGGTGGTCTCC
TTCGACGCCCTGGTTGCGTACTCCTTCGAGGTGACTGGATACCAGCCTTTCATCCTAACAGGGGAGACA
GCTGAGGGGCTCCCTCCAGTCCGGATCCCGCCTTCTCAGTGACCACAGCCAACGGGACGATCTCCTTC
ACCGAGATGGTGCAGGACATGGGAGCCGGGCTGGCCGTGGTGGCCCTGATGGGCCCTCCTGGAGAGCATT
GCGGTGGCCAAAGCCTTCGCATCTCAGAATAATTACCGCATCGATGCCAACCAGGAGCTGCTGGCCATC
GGTCTACCAACATGTTGGGCTCCCTCGTCTCCTACCCGGTACAGGCAGCTTGGACGGACAGCC
GTGAACGCTCAGTCGGGGGTGTGCACCCCGCGGGGGCCTGGTGACGGGAGTGTGGTGTCTGTCT
CTGGACTACCTGACCTCACTGTTCTACTACATCCCCAAGTCTGCCCTGGCTGCCGTCATCATCATGGCC
GTGGCCCCGCTGTTGACACCAAGATCTTCAGGACGCTCTGGCGTGTAAAGAGGCTGGACCTGCTGCC
CTGTGCGTGACCTTCTGCTGTGCTTCTGGGAGGTGCAGTACGGCATCCTGGCCGGGGCCCTGGTGTCT
CTGCTCATGCTCCTGCACTCTGCAGCCAGGCCTGAGACCAAGGTGTGAGAGGGCCGGTTCTGGTCTG
CAGCCGGCCAGCGGCTGTCTTCCCTGCCATGGAGGCTCTGCGGGAGGAGATCCTAAGCCGGGCCCTG
GAAGTGTCCCGCCACGCTGCCTGGTCTGGAGTGCACCCATGTCTGCAGCATCGACTACACTGTGGTG
CTGGGACTCGGCGAGCTCCTCCAGGACTCCAGAAGCAGGGCGTGCCTGGCCTTGTGGGCTGCAG
GTCCCCGTTCTCCGTGCTCTGTCCGCTGACCTGAAGGGGTTCCAGTACTTCTCTACCCTGGAAGAA
GCAGAGAAGCACCTGAGGCAGGAGCCAGGGACCCAGCCCTACAACATCAGAGAAGACTCCATTCTGGAC
CAAAGGTTGCCCTGCTCAAGGCATAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_173626

Insert Size: 1821 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_173626.3](#)

RefSeq Size: 2927 bp

RefSeq ORF: 1821 bp

Locus ID: 284129

UniProt ID: [Q86WA9](#)

Cytogenetics: 17q25.3

Protein Families: Transmembrane

MW: 65.3 kDa

Gene Summary: This gene encodes a member of the solute linked carrier 26 family of anion exchangers. Members of this family of proteins are essential for numerous cellular functions including homeostasis and intracellular electrolyte balance. The encoded protein is a sodium independent sulfate transporter that is sensitive to the anion exchanger inhibitor 4,4'-diisothiocyanostilbene-2,2'-disulfonic acid. Alternate splicing results in multiple transcript variants.[provided by RefSeq, Oct 2009]
Transcript Variant: This variant (2) differs in the 5' UTR compared to variant 1. Variants 1, 2, 3 and 4 encode the same protein.