

Product datasheet for **SC318046**

Rhotekin 2 (RTKN2) (NM_145307) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Rhotekin 2 (RTKN2) (NM_145307) Human Untagged Clone
Tag:	Tag Free
Symbol:	Rhotekin 2
Synonyms:	bA531F24.1; PLEKHK1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC318046 representing NM_145307.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCACTG
GATCCGGTACCGAGGAGATCTGCCGCCACGATCGCC
ATGGAGGGGCGGAGCCTGAGGGGTCTGCGCTCCGCTGGCGGGGCTTCCACCCAGCAGGACTGCAAC
ATTCAAGAAAAATAGACTTAGAAATTCGAATGCGAGAAGGAATATGAAACTCCTTTCTCTGAGCACT
CAGAAAGATCAAGTTTTACATGCAGTTAAGAATCTCATGGTGTGCAATGCTCGACTAATGGCCTATACA
TCGGAGCTACAGAAATTAGAAGAACAGATTGCAATCAGACTGGAAGATGTGATGTGAAATTTGAAAGT
AAAGAACGAACAGCATGTAAAGGAAAGATTGCCATATCAGATATTCGAATACCACTAATGTGAAAGAC
TCTGATCACTTCAGCAATAAAGAACGATCAGCAGCTATGCCATTTTTTGTATTCAAAATGGGAGCT
AATGTGTTTGATACTGATGTGGTGAATGTGATAAAACAATCACAGATATATGTTTTGAAATGTAACC
ATATTTAATGAAGCAGGGCCAGACTTTCAGATAAAGGTGGAAGTGTATAGTTGCTGTACAGAAGAATCT
TCTATAACCAACACACAAAAAAGCTAGCTAAGAACTCAAGACATCTATAAGTAAAGCTACAGGAAAG
AAAATAAGTTCAGTGCTCAAGAAGAGGATGATGAAATGTGCTTCTCCTCAGCTCTGCTGTTTTGGT
GTAAAGTATAATTTGCTAGCTCACACTACCCTAACCTTGAAAGTGTGAGGATAGTTTCAAGACCCAT
AATCTGTCTATTAATGGAATGAGGAGTCTTCATTTTGGCTTCCCTGTATGGCAACATGTGCTGCCGA
CTAGTTGCCAGCCAGCTTGTATGGCTGAGGATGCATTTGCAGGATTTCTTAATCAGCAGCAATGGTA
GAAGGTCTGATTAGTTGGAGAAGGTTGATTGTGTTTTGCGAGGAGGTAACTCTATTGTTTTACAGT
CCAGAGGAAATGAAGCTAAAGTGAACCCAGCTTTGGTAGTACCCATTAAACAGGAAACCAGAATCCGG
GCAATGGATAAGGATGCCAAGAAAAGATCCATAATTTCTCTGTCATCAATCCTGTTCTGGACAAGCT
ATAACTCAGATTTTGCAGTTGACAATAGAGAAGATCTCAGAAGTGGATGGAAGCCTTCTGGCAGCAT
TTCTTTGATCTTAGCCAATGGAAGCACTGTTGTGAAGAATTATGAAATTGAGATTATGTCACCACGG
AAACCACCTTTGTTCTTGACAAAAGAGGCAACCTCAGTCTACCATGATATGAGCATTGATTCACCTATG
AAACTTGAAGTTTAACGATATAATACAAAAAAATGAAGAGACAAATGGGAGTTCCTTATTGGT
CAGCATGAAGAATCCTTACCACCTCCTTGGGCCACACTCTTGTATGGTAACCATCAATGGTCATCCAG
AAAAAGGTACTGTATCCTGCAAGTGAGCCATTACATGATGAAAAAGGAAAAAGAGACAAGCTCCCTT
CCTCCTTCTGATAAACTCCATTAGCTTAAATCACAGAGTAACACAGATCAATTGGTTAAGGACAAC
TGGGAAAAACAAGTGTATCTCAGACATCGTCTTGGATACCAAACTATCAACTCTAATGCATCACTTA
CAGAAACCAATGGCTGCTCCTCGAAAATCTGCTGCCAGGAGAAATAGATTAAAGTATGGTGAGCAC
ACAGACACTAAACCAATTTTGAAGCAAGCCAGTGCCAGCTCCAAGGCAGAAATCCATCAAAGACATT
CTGGACCTTAGATCATGGCTGCAGGCACAAGTATAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: Sgfl-MluI

ACCN: NM_145307

Insert Size: 1830 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_145307.3</u>
RefSeq Size:	6915 bp
RefSeq ORF:	1830 bp
Locus ID:	219790
UniProt ID:	<u>Q8IZC4</u>
Cytogenetics:	10q21.2
Domains:	PH
MW:	69.3 kDa
Gene Summary:	May play an important role in lymphopoiesis.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (1) represents the longer transcript and encodes the longer isoform (1).