

Product datasheet for **SC318051**

PKD2L2 (NM_014386) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PKD2L2 (NM_014386) Human Untagged Clone
Tag:	Tag Free
Symbol:	PKD2L2
Synonyms:	TRPP5
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC318051 representing NM_014386.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCTGAGGCGTCACGGTGGCACCGAGGCGGGCTTCGAAACATAAGTTGCATTACAGAAAGGAAGTA
GAAATTACAACCACACTTCAGGAATTGTTACTCTACTTTATTTTTTAATAAACCTATGTATATTGACT
TTTGGGATGGTAAACCCACATATGTATTACTTAAACAAGTTATGTCATCTCTATTTTTGGACACTTCT
GTGCCTGGTGAAGAAAGAACCAACTTTAAGTCCATTTCGACGATAACTGATTTTTGGAAGTTTATGGAA
GGACCCCTTTTGAAGTCTGTACTGGGATTCATGGTACAATAACCAGCAGCTGTATAATTTAAAGAAC
AGCAGTCGCATCTACTATGAAAATATACTTCTAGGAGTCCCAGAGTTCGCAACTAAAAGTCCGCAAC
AACACATGCAAAGTCTATTCATCTTTTCAGTCTTTGATGAGTGAATGTTATGGCAAATATACTTCTGCA
AATGAAGACCTCTAATTTTGGCCTTCAAATTAATACTGAATGGAGATTTCTACTTCTAATACCAAC
TCCCTTGGCACTGGGGATTTCTTGGTGTTCACGAAATGGGGATACATTTTCACTTTATCAAAATCG
AAATCTGAAACAAAAACAAGTTCATTGACCTTCGACTGAACAGCTGGATCACAAGAGGGACTAGAGTT
ATTTTTATTGATTTTTCTTATATAATGCTAATGTAAATCTATTTTGTATTATCAGATTGGTGGCAGAA
TTCCCTGCAACTGGAGGAATACTTACTTCATGGCAGTTTTACTCTGTGAAGCTCCTCAGATATGTTAGC
TACTATGACTATTTTATTGCTTCTGTGAAATCACATTCTGTATTTTTCTTTTTGTCTTCAACACAAA
GAAGTCAAAAAATAAAGAATTTAAGTCTGCCTATTTCAAAGTATTTGGAAGTGGCTAGAAATTGCTA
CTTTTGTGTGTGTGTGTGTGGCTGTTTCTTCAACACATACTATAATGTACAAATTTTTCTTACTT
GGACAGCTGTTGAAAAGTACTGAAAAATTTTGCATGGATAAAGATTTCAAATTCATAAGCTTTAACAAG
ACAATGTCTCAGCTGTCATCAACCTTGTCCCCTTGTGTTAAAGACATAGTAGGATTTGCCATCATGTTT
TTTATAATATTCTTTGCTTATGCCAGTTAGGATTTCTTGTTTTTGGATCACAAGTTGATGACTTTTCC
ACTTTTCAGAATTCATATTTGCACAATTTGCAATTGTTCTTGGAGATTTTAATTTTGTGTTATTCAG
CAAGCCAATCCTATCTTGGGACCCATTTACTTCATCACTTTCATCTTTTTTGTGTTCTTTGTCTGCTG
AATATGTTCTTGGCAATTATTAATGATACCTATTCTGAAGTAAAAGCTGACTATTCAATAGGCAGAAGG
CTAGATTTTGAAGTTGGCAAAATGATTAACAGAGTTACAAAAATGTTCTCGAGAAATTCAGACTGAAG
AAAGCTCAAAAAGATGAAGACAAGAAAACCAAAGGCAGCGGAGATTTGGCTGAACAAGCCAGAAGAGAA
GGCTTTGACGAAAATGAGATTCAAACGCAGAGCAGATGAAAAATGAAAAGAGAGGCTTGAGAAAAAG
TATTATTCTATGAAAATCAAGATGACTACCAGCTGTCACTCAAGAAGAATTTGAGATGGCACCACA
ACCAAGTACAAAATGAGATTCTCTGAGTGCCTGACAAAACGAATTTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_014386

Insert Size: 1842 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_014386.3](#)

RefSeq Size: 2218 bp

RefSeq ORF: 1842 bp

Locus ID: 27039

UniProt ID: [Q9NZM6](#)

Cytogenetics: 5q31.2

Protein Families:	Druggable Genome, Ion Channels: Transient receptor potential, Transmembrane
MW:	72.5 kDa
Gene Summary:	<p>May function as a subunit of a cation channel and play a role in fertilization.[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (1) encodes isoform 1 of 602aa. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>