

Product datasheet for **SC318054**

ANKRD55 (NM_024669) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ANKRD55 (NM_024669) Human Untagged Clone
Tag:	Tag Free
Symbol:	ANKRD55
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC318054 representing NM_024669.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGATGAGACAGGCTACCATGGATTTCAGCACCCCTTCTGTGTTGATCAGCAAAGAGGTGACTCATCT
GAGGAAGTTGACCTGACCATGGTTTATCAAGCAGCCTCTAATGGAGATGTCAATGCTCTGACTGCAGTG
ATTCGGGAAGACCCTTCTATCCTAGAATGCTGTGACAGTGAAGGATGCACGCCCTTGATGCATGCGGTT
TCTGGACGTCAAGCGGACACAGTGAAGCTGCTGTTGAAGATGGGAGCCAATATTAACATGCAGGATGCT
TATGGCCGCACAAGTTTATGCCTGGCCACCTACCTGGGCTGGCTTGAAGGCTGTGTGAGTCTACTCAGA
AACGGTGCCAAGCACAATATCCCAGATAAAAAATGGCCGCCTGCCACTGCATGCTGCCACTGCTGAGCCC
GATATGAGGCTCCTCACGGTCTGTTGCAACAGTCAACATCAGCGAGATTAATCACCAGGACAATGAG
GGATGACACCACTCCACTGGGCGGCTTCCACAACCAGCCTCAACACACAAAAATGCTGCTGAAGAAG
GGGCAGACCCACCCTTGTGGATAAAGACTTTAAAACCGCTCTCCACTGGCAGTCCAGAGTGGAAAT
AGGATTCTGTGCTCCATCATTCTGAGCCATCACCAGGGCCGTCATAATCAACTATGATGATGAGAGT
GGGAAGACATGTGTACATATCGCAGCGGACGGGCTTACGCGATATTATTCATGAGCTGGCAAGAGTC
CCTGAGTGTAACTGCAGGCTCTGGATGTGGATGACAGGACACCTCTGCACTGGGCTGCAGCTGCAGGG
AAGGCCGAATGTGTCCAGTCACTGCTGGAGTTGGGAATGGACAGCAACCTGCGGGACATCAATGAGAGC
ACGCCCTTGGCCTATGCCCTGACTGCGGTACACGGCGTGTGTCAAACCTCTCCCAAGAGAGCAGA
ACAGAGCCTACTCGACCCCTCCCTCCCAGAGCAGTGGGCCCCAGAAGAAGGAGAGACGGTTCAACGTG
CTCAACCAAAATTTCTGCAAAAACAAGAAAGAAGAGCAGAGAGCCCATCAGAAGGATCCCAGCAGGGAC
CGATACAGAGAGGAGGACACCTCAGAAGTCAATGACATCATCACCACCTTGGATAGCATCGTGGGTACC
AACTGCCAAGAACAAGCCCTGGTGATCAGGTGGCTATGGTTGAATTTAAGAAGAAAACCTCAGACAATTCA
AAATATCTTTACCAGAAAAGAAACCGCTGGCCGTAAGGGGCTTCCACCAATCAGAACGCAGAGTCTC
CCACCCATCACCTGGGCAATAACTTCCCTAACAGCCTCCCATAGGGCCACTTCCCATGCAGGCTGAGC
TCTGCTCCTCATCATATGGCCCAGCGATCTCAGAAAAGTCAAGTGAAGTGAAGGATTTATTAATAACAGA
ACTGGCTGCCAGATGTTACTAGATAACCCCTGGAAGAGTGATTCTAATCAGGATTTTCTACAAAGTT
TGGACTGTGCTTCTTCTGATAAGCTGCTGGACAGATTGCTCAGTGTCCGGCTGGTCCCAAGAGGTC
TCCGTGCCACCACACCTTCGCCATCTACATAATCCATCATCAGGACAAAAATTTTCCAGCATCTTCCCA
AACAGACAAAAATCAGGGATCTTCTTTCACTCGGAACAACCTAGCTCCCCTACCAGATCAAAAAATTT
CTATCTGGAGAACCTCTGCGGACAAACCGAGTGTCTTCTGCAATTCGAAGTCAACGAAGACACAGACA
GCAGCAGAAGAGAGTGAACATTCTGCCAACCCACCAGTGATGAAAATTAAA
ACGCGTACGCGGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTAAACGGCCGGC
  
```

- Restriction Sites:** SgfI-MluI
- ACCN:** NM_024669
- Insert Size:** 1845 bp
- OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
- OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
- Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_024669.2](#)

RefSeq Size: 2502 bp

RefSeq ORF: 1845 bp

Locus ID: 79722

UniProt ID: [Q3KP44](#)

Cytogenetics: 5q11.2

MW: 68.4 kDa