

Product datasheet for **SC318068**

KLHL6 (NM_130446) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	KLHL6 (NM_130446) Human Untagged Clone
Tag:	Tag Free
Symbol:	KLHL6
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC318068 representing NM_130446.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGTTGATGGCAGGACAAAGGGGCGCCTGGACCATGGGTGATGTGGTTCGAGAAGAGTTTGAAGGGCCC
GTGGCACCTTCTACAGATGAGCCCTCCAGAAAACAGGAGACTTGGTCGAGATCTTAAATGGGAAAAG
GTCAAATTTGACGACGCGGGACTCTCCTTAATTCTTCAGAATGGCCTGGAAACCCCTGCGAATGGAAAAC
GCTCTGACAGATGTCATCTTGTGTGGACATTAGGAATTCTCCTGCCACCGCGTGGTGCCTGGCCGA
GCCAGCAACTATTTAGGGCCATGTTCTGCAACGATTTAAAGGAAAAGTATGAGAAAAGGATCATTATT
AAAGGGGTTGATGCTGAGACCATGCACACTCTGTTGGACTACACGTACACCAGCAAGGCGCTGATCACC
AAGCAGAAATGTCCAGCGGGTCTGGAGGCTGCTAACCTCTCCAGTTCCTGCGGATGGTGGATGCCTGT
GCCAGTTCCTCACTGAAGCCTGAACCCAGAAAACGCGTTGGAATACTGAGGCTGGCTGACACACAC
TCGCTGGACAGTCTAAAGAAGCAGGTTTCAGAGTTACATCATTCAAACCTTTGTGCAGATTCTGAACCT
GAGGAGTTTCTTGACCTGCCCGTGGACACTCTGCACCACATCTTGAAGAGTGATGACCTTACGTGACC
GAGGAGGCTCAGGTGTTGAGACCGTGATGAGCTGGGTCCGGCACAAGCCATCAGAACGACTCTGCTTA
CTCCCCTATGTCCTCGAGAACGTGCGCTTACCGCTTCTGGACCCGTTGGTACTTTGTGGAGACGGTGAA
GCAGATCCTCTCATCAGGCAGTGCCAGAGGTCTTCCCGCTGCTCCAGGAAGCCAGGATGTACCACCTT
TCTGGCAATGAGATCATTTCGGAACGCACCAAGCCAGGATGCATGAGTTCAGTCTGAGGTGTTTCATG
ATCATTGGCGGCTGCACGAAGGATGAACGGTTTGTGGCAGAGGTGACCTGCCTGGACCCCTGAGGCGC
AGCCGCTGGAGGTGGCAAGCTCCCGCTAACAGAGCATGAGCTGGAGAGTGAGAATAAGAAGTGGGTG
GAGTTTGCATGCGTGACATTGAAAAATGAGGTTACATCTCAGGTGGCAAAGAAACACAGCATGATGTT
TGGAAATATAATTCTTCGATCAACAAGTGGATTTCAGATTGAGTATTTAAACATAGGCCGCTGGAGGCAT
AAGATGGTGGTGTGGGTGGCAAGGTCTATGTGATCGGAGGCTTTGACGGCTTACAGAGAATCAACAAT
GTGGAGACCTACGACCCCTTTCACAACCTGCTGGTTCAGAGGCCGACCCCTCCTTGTCCATGTCAGTTCC
TTTGCAGCCACCAGCCATAAGAAGAAGTGTATGTGATCGGGGAGGGCCCAATGGGAACTGGCCACA
GACAAGACTCAGTGTTATGACCCCTCCACCAACAAGTGGAGTTTGAAGGCGGCCATGCCCGTGGAGGCT
AAATGCATCAATGCAGTGAGTTTCCGGGACCGCATCTATGTCGTTGGTGGGGCCATGAGAGCGCTGTAC
GCCTACAGCCCGCTGGAAGACAGCTGGTGCCTGGTACCAGCTCAGCCACGAGCGGGCCAGCTGCGGT
ATCGGCCCTGCAACAACCGCTCTACATCACCAGCGGGGACGAGAAGAACGAGGTTATCGCCACG
GTGCTGTGCTGGGACCCGAGGCCAGAACTGACAGAGGAGTGCGTCTGCCCCGGGCGTGTGCGCAC
CAGGGCAGCGTCACCATCAGGAAGTCGTACCCACATCCGACAGGATCGTCCCCGGAGCAGTGTCTGTC
TGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_130446

Insert Size: 1866 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_130446.2](#)

RefSeq Size: 6320 bp

RefSeq ORF: 1866 bp

Locus ID: 89857

UniProt ID: [Q8WZ60](#)

Cytogenetics: 3q27.1

Domains: BTB, Kelch

MW: 70.4 kDa

Gene Summary: This gene encodes a member of the kelch-like (KLHL) family of proteins, which is involved in B-lymphocyte antigen receptor signaling and germinal-center B-cell maturation. The encoded protein contains an N-terminal broad-complex, tramtrack and bric a brac (BTB) domain that facilitates protein binding and dimerization, a BTB and C-terminal kelch (BACK) domain, and six C-terminal kelch repeat domains. Naturally occurring mutations in this gene are associated with chronic lymphocytic leukemia. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Feb 2017]